- A. Pala
B. Display
C. Botón On/Off Tara
D. Botón Prop
E. Botón Hold

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
1. Technical specifications

- Functions:** Tare adjustment; on/off key; display blocking function; "proportion calculation" function.
- Batteries:** 3V 1*CR2032 button type lithium battery (Included in the kit).
- Working temperature:** 0°C - 40 °C

The kit includes two interchangeable spoons: one of them for small doses, and a pan type spoon for bigger doses.

2. Description

This spoon with weighting function has been specially designed for the dental sector. Compact and attractive, this reliable and multifunction device weighs accurately plasters, composites and other dental materials.

It allows storing a given proportion in the memory, and afterwards applying the proportion to the element we want to weight. For instance, if we store in the memory a 24/100 type proportion (24 g of element "A" for 100 g of element "B", or in other words, a 24% of "A" over "B"), whenever we press the "prop" key, the display will show the amount of the other element required to maintain the proportion.

3. Operation

- Press the "ON/OFF" key. The spoon will start up and will begin its internal test. Once the test is over, the display will show "0.0". The spoon is ready to weight.
- Whenever you want to put the display in a "0.0" position, just press the "TARE" key (the "ON/OFF" key and the "TARE" key are the same key). Whatever the value the display is showing at the moment, the controls will carry out a reset, and the display will show a "0.0" weight.
- Weight holding function: sometimes it will be advisable to maintain a value in the display in order to prevent annoying oscillations. If you want to hold the value shown in the display, –once it remains static– just press the "HOLD" key.
- If you want to store in the memory the amount of the second element to be mixed with the element that is now being weighted, proceed as follows:
 - Press the "PROP"+"HOLD" keys for a few seconds. You will notice that the display shows the value of the last proportion introduced. The first digit will start flickering. Then, press the "PROP" key and adjust the value of the first digit.
 - Then press the "HOLD" key. The second digit of the proportion will start flickering. Adjust the value of the second digit. Then, press the "HOLD" key and you will see that the controls return automatically to the "Stand-by" situation. The proportion has been stored in the memory.
- If you want to know what amount of the second element should be mixed with the value shown in the display, first of all press the "HOLD" key to hold the value in the display. Afterwards, whenever you press the "PROP" key, the display will show the amount in grams of the second element to be mixed in order to maintain the proportion stored in the memory.

4. Preventive measures


Important:

The controls need some time to process the sequence of instructions to be programmed. The controls have to wait for each step to finish, before starting with the next one. If the keys are pressed too quickly, the controls might get clogged up.

- Do not force the spoon of the device, and protect it from bumps and falls. The accuracy of the scales depends to a great extent on how the weight sensor is cared for.
- Do not use the scales in case of high humidity, heat, dust or vibrations.
- Always use the scales in a horizontal position (handle horizontal). Otherwise, there might be errors in the weighting.
- Whenever the display shows the "low battery" sign, proceed to replace the battery.
- Keeping the device clean is essential for a proper maintenance. Clean it with a cloth. Do not use solvents or other chemical products.
- Do not get the weighting spoon wet with water or other liquids. Liquids could cause internal short-circuits.
- If you are not going to use the spoon for a long time, please, remove the battery.

Important:

The spoon is unsuitable for weighting in economic transactions.

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