R-010098

E.....	ARTICULADOR MODELO "ARCON"	4
GB	ARTICULATOR "ARCON"	8
F.....	ARTICULATEUR MODELE "ARCON"	12
P.....	ARTICULADOR "ARCON"	16

MESTRA®

TALLERES MESTRAITUA S.L.


Txori-Erri Etorbidea, 60

48150 SONDIKA - BILBAO - ESPAÑA

Tfno. + 34 944530388 - Fax + 34 944711725

E-mail: mestra@mestra.es - www.mestra.es

ARTICULADOR / ARTICULATOR / ARTICULATEUR

Nº	Descripción / Description / <i>Description</i> / Descrição
1	Tornillo de fijación de la mesa incisal / Incisal table screw / <i>Vis de fixation de la table incisale</i> / Parafuso de fixação da mesa incisal
2	Rama inferior / Lower frame / <i>Branche inférieure</i> / Ramo inferior
3	Mesa incisal padrón de acrílico / Incisal table (acrylic) / <i>Table incisale en acrylique</i> / Mesa incisal padrão acrílica
4	Tapa de acabamiento de la columna derecha / Finishing cover of right column / <i>Enjoliveur du bras droit</i> / Tampa de acabamento da coluna directa
5	Aguja indicadora del plano oclusal / Occlusal plan indicator / <i>Repère indicateur du plan occlusal</i> /
6	Pino incisal con agujero / Incisal pin with hole / <i>Jauge incisale avec trou</i> / Pino incisal
7	Clavija traba de la posición central del articulador / Lock pin of centric position (articulator) / <i>Goupille de blocage de la position centrale de l'articulateur</i> / Pino trava de posição cêntrica do articulador
8	Clavija guía del arco facial / Face bow's guide pin / <i>Goupille guide de l'arc facial</i> / Pino guía do arco
9	Rama superior / Upper frame / <i>Bras supérieur</i> / Ramo superior
10	Clavija de apoyo de la rama superior / Upper frame support pin / <i>Tige de soutien du bras supérieur</i> / Pino de apoio do ramo superior
11	Tapa de acabamiento del imán derecho / Magnet finishing cover (Right) / <i>Cache de finition de l'aimant droit</i> / Tampa de acabamento do imã direito
12	Imán estabilizador / Stabilizer magnet / <i>Aimant de stabilisation</i> / Imã estabilizador
13	Etiqueta de identificación del modelo / Model identification label / <i>Etiquette d'identification du modèle</i> / Etiqueta de identificação do modelo
14	Tapa de acabamiento del imán izquierdo / Magnet finishing cover (Left) / <i>Cache de finition de l'aimant gauche</i> / Tampa de acabamento do imã esquerdo
15	Imán de fijación de la placa de montaje / Mounting plate magnet / <i>Aimant de fixation de la plaque de montage</i> / Imã de fixação da placa de montagem
16	Placa de montaje riel / Mounting plate (track) / <i>Plaque de montage</i> / Placa de montagem trilho
17	Tornillo de fijación del pino incisal / Incisal pin screw / <i>Vis de fixation de la jauge incisale.</i> / Parafuso de fixação do pino incisal
18	Tornillo de fijación del elemento condilar / Condylar element screw / <i>Vis de fixation de l'élément du condyle</i> / Parafuso trava do elemento condilar
19	Anillo O-ring de la clavija traba / O-ring of lock pin / <i>O-ring de la goupille de blocage</i> / Anel o´ring do pino trava
20	Elemento condilar / Condylar element / <i>Élément condylien</i> / Elemento condilar
21	Tapa de acabamiento de la columna izquierda / Finishing cover of left column / <i>Enjoliveur du bras gauche</i> / Tampa de acabamento da coluna esquerda
22	Encaje del imán / Magnetic fitting / <i>Capsule de l'aimant</i> / Caneca do imã
23	Clavija guía de la placa de montaje / Mounting plate guide pin / <i>Butée guide de la plaque de montage</i> / Pino guía da placa de montagem
24	Tornillo de fijación de la placa de montaje / Mounting plate screw / <i>Vis de fixation de la plaque de montage</i> / Parafuso de fixação da placa de montagem

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INTRODUCTION

This manual is a general presentation of your product and it will give you important details to help you to solve possible problems. Please, read it and keep this with you.

Attention: Before using your MESTRA articulator, please read carefully all the Instructions Manual.

Attention: The use of the equipment by professionals not qualified can cause damage or harm to the product and/or patient.

PRODUCT DESCRIPTION

The articulator is an instrument used to simulate the maxillo-mandibular relation and movements of a patient in a laboratory, with the purpose of studying the occlusion and production of dental devices that will be used by the patient. These devices include complete dentures, partial dentures, bridges, crowns and bite plate, among others.

This Manual contains only basic and simplified information regarding the use of the equipment, which do not replace the need of a specific course of Occlusion and/or Prosthesis, offered in many Teaching Institutions and Dental Faculties, as a compulsory discipline.

The articulator use is restricted to qualified professionals.

It is important to emphasize, however, that the R-010098 articulator provides a simple, fast and highly precise way of reproducing the human mandible movements, enabling the dental professional to carry out corrective and restorative dentistry tasks more easily, rapidly and at lower costs than those of traditional, time-consuming techniques involving expensive and highly complex equipment. Furthermore, the semi-adjustable articulator produces far more accurate results than those produced by “simple hinge” articulators, which involve arbitrary mounting of the mandibular cast on the device and whose movements are also limited.

The use of the semi-adjustable articulator is therefore recommended for most prosthetic, occlusal and rehabilitation work. The technique is simple, fast and easy, offering highly satisfactory results for both the patient and professional.

MOUNTING THE CASTS ON THE ARTICULATOR

Preparing the articulator

Aiming as main characteristic simplicity and easiness to work, the R-010098 articulator was developed with the Bennet and Condyle angles fixed at the average of 15° and 30°, respectively. The intercondylar distance was also fixed at the average of 110 mm. Therefore, this model does not require preliminary adjustments for these components, such as the other models previously mentioned. In order to lock the R-010098 articulator in the centric position, there is an innovative system with two lateral pins, denominated: "centric lock pin". This pin has 3 main positions: first position - totally closed (locks the articulator in the centric position – fig. 01); second position - intermediate (releases the articulation movements – fig. 02) and third position (allows to remove the upper from the lower frame – fig. 03).

Fig. 01

Fig. 02

Fig. 03

Interocclusal registration

To set up the mandibular (lower) cast in the articulator, you should have a record correlating the upper and lower dental arches, according to the purpose of the mounting, in one of the following ways:

- Maximum Intercuspatation (MI);
- Centric Relation (CR);

These records can be made using material of the operator's preference: wax, resins, addition or condensation silicon, etc. (fig. 04 and 05).

To obtain the M.I registration, place the selected material to register the patient's arcade and ask him to occlude.

There are several techniques to obtain the R.C registration, i.e. Peter Dawson's technique (bilateral manipulation of the mandible), Lucia's JIG technique and James Long (use of plastic spacer) among others.

Wax registration

Fig. 04

Silicone registration

Fig. 05

Note: For totally edentulous patients, this registration is made when the upper and lower wax rollers are joined based on the correct dimensions of the patient.

Fig. 06

Fig. 07

Fig. 08

Fig. 09

Fig. 10

Mounting the mandibular cast

- a) Place the Incisal Guide Pin in the Upper Frame of the articulator with its rounded tip pointing downward so that the upper and lower frames are parallel, i.e. on the zero marking of the Incisal Guide Pin (fig. 06).
- b) Now turn the articulator upside down on the laboratory bench and affirm the lower cast, with retentions and previously hydrated, upon the interocclusal register that should be placed in the mounted upper cast (fig. 07). In totally edentulous patients, this was obtained by joining the wax rollers together.
- c) We recommend fixing the casts with rubber bands or staples (fixed with "godiva" or wax) so that they remain in position until the plaster hardens. Place a small amount of plaster on the lower part of the mandibular cast and a small amount on the mounting plate of the lower frame of the articulator to fill in the gap between them (fig. 08).
- d) Make sure the condylar elements are duly positioned and close the lower frame until the Incisal guide pin touches the Incisal table. Afterwards, fasten the articulator frames using rubber bands to prevent possible distortion occurred by the plaster expansion (fig. 09).

Note: The sides of the mounting plates, which are in contact with the articulator frames, should be devoid of plaster.
- e) After the plaster hardens, turn the articulator back to its upright position (lower frame resting on the laboratory bench) and complete the work, filling the cast fixation towers with plaster for the finishing touches (fig. 10).

Features of the railmounting plate

The possibility of reusing the Mounting Plate without breaking the model requires care with the side finishing of the plaster on the plate, i.e., the plate should not have plaster at the lateral border (fig. 11).

In order to take off the model, the plaster should be removed of from the retention hole (fig. 12) and then the model should be hold while tapping the plate towards the indicated arrows located on the lower face of it (fig. 13).

PS.: Once the model has been taken off the Mounting Plate, we can not guarantee the returning of it in the initial position with precision, since it depends on the technical quality of the work and the materials used (especially the plaster).


Fig. 11

Fig. 12

Fig. 13

MAINTENANCE AND CLEANING

- This manual only contains basic and simplified information regarding the use of the instrument and under no circumstances it is a substitute for a proper training course. The articulator and face-bow are products aimed at the exclusive use of qualified dental Professionals.
- The R-010098 articulator is not interchangeable, i.e. the casts mounted on one articulator model should not be transferred (mounted) to any other articulator. Therefore, Mestra does not guarantee precision when casts mounted on one articulator are transferred to another.
- Before using the articulator and face-bow, the professional should check the instrument for possible damage, distortion of the incisal pin, face-bow's parallelism, proper centricity of the incisal pin in relation to the lower frame, etc. Should any abnormality be found, contact your distributor.
- Mestra recommends only warm water for the general cleaning of the articulator and facebow.
- When using "godiva" or any other registration material, special care should be taken to avoid excessive heating of the material, which might burn the patient's mouth. Mestra urges the user to carefully follow the instructions of the registration material manufacturer.
- Because the articulator is a precision instrument, it should be handled, transported and stored carefully.

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