R-080130

E.....	HORNO DE PRECALENTAMIENTO DIRECTO 50D.....	3
GB.....	50D DIRECT PREHEATING FURNACE.....	6
F.....	FOUR DE CHAUFFE SANS PALIER 50D.....	9

MESTRA®




TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

48150 SONDIKA - BILBAO - ESPAÑA

Tfno. + 34 944530388 - Fax + 34 944711725

E-mail: mestra@mestra.es - www.mestra.es

- | | |
|---|---|
| A | Conexión extractor / Extractor connection / <i>Connexion extracteur</i> |
| B | Control / Controls / <i>Contrôle</i> |
| C | Interruptor general / Main switch / <i>Interrupteur général</i> |
| D | Piloto resistencia / Resistor pilot light / <i>Témoin résistance</i> |
| E | Piloto fin extractor / Extractor stop pilot light / <i>Témoin fin extraction</i> |
| F | Tecla SET / SET key / <i>Touche SET</i> |
| G | Display consigna / Set point temperature display / <i>Display consigne</i> |
| H | Display mufla / Muffle display / <i>Display moufle</i> |
| I | Tecla selección / Sélection key / <i>Touche sélection</i> |
| J | Tecla + / + key / <i>Touche +</i> |
| K | Tecla - / - key / <i>Touche -</i> |
| L | Magnetotérmicos / Magnetothermic switches / <i>Disjoncteurs thermiques</i> |

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The R-080130 is a practical furnace specially designed for heating cylinders and crucibles within the field of dental prosthesis. It is perfect for those works that do not require a classical programming cycle of ramps and plateaus. Just select the set point temperature, and the furnace will start heating directly until this temperature is reached. Then, this temperature will be kept for a time programmed by the user.

Its digital controls ensure an amazing temperature accuracy and prevent the effects of thermal inertia. Besides, it features a connection for an external catalytic converter, and this element will operate from the beginning of the cycle up to the stopping temperature programmed by the user. This is an especially quick furnace, capable of reaching 1000 °C in 50 minutes approximately.

Caution: before starting the unit, read these instructions.
Store these instructions for further consultation.

INSTALLATION

- Carefully unpack the unit. Inside the box there should be an SiO₂ hearth plate. Check also that the muffle lining the inside of the furnace has suffered no damage during transport.
- Set the furnace up on a smooth, firm, horizontal base away from any source of vibration. Allow enough space for ventilation. There should be no other objects less than 15 cm from the unit. Ensure especially that the ventilation tube is free and at least 20 cm from the closest object or wall.
- During normal furnace operation fumes may filter through the ventilation tube. This tube should therefore be placed below an extractor hood.
- Connect the unit to a 230 V, 50/60 Hz, earthed power outlet which protected by a thermal magnetic cut-out with a suitable amp rating (consult your installer).
- Before using the furnace heat it up to 300 °C and hold it there for one hour to evaporate off any moisture in the muffle.

PROGRAMMING

Set point temperature

It is the temperature that the muffle will reach.

- Keep the door open to prevent the furnace starting to heat. Then, turn the main switch (C) on. The upper display (G) will show the temperature inside the muffle
- Press the SET key (F). The lower display (H) will show the message “SP”, the set point temperature.
- Use the ▲▼ keys to select the set point temperature. You can use the ◀ key to shift the digit.

Waiting time

It is the time that the furnace will keep the set point temperature once this is reached. The unit is factory pre-programmed for 1 hour waiting time.

- Once the set point temperature has been inserted, press the SET key again. The lower display will show the message “St”, waiting time.
- Use the ▲▼ keys to select the heating time in hours and minutes (*hh:mm*). You can use the ◀ key to shift the digit.
- Press the SET key to store the programming.

Extractor stop

The furnace features a power outlet on its side (A) to connect directly a fume extractor (R-080125) or a catalytic converter (R-080124). The extractor will operate until the temperature stored in the controls is reached. The device is factory supplied with a stopping temperature of 850 °C. This temperature can be easily modified as follows:

- Keep the SET key pressed for 3 seconds. The lower display will show the message "Al1", stop of the extractor.
- Use the ▲▼ keys to select the extractor stopping temperature. You can use the ◀ key to shift the digit.
- Keep the SET key pressed for 3 seconds to store the programming.

Once the extractor stopping temperature has been reached, it will stop and the pilot light (E) will light up.

OPERATION

- Place the SiO₂ hearth-plate on the hearth, taking care not to damage the muffle.
- Put in the cylinders and crucibles to be heated. When the furnace is operating, always use gauntlets and tongs. Talleres Mestraitúa, S.L. recommends MESTRA tongs.
- To program the furnace, follow the steps explained in the previous section.

The pilot (D) will light up when the controls are powering the heating resistor.

- Shut the door. The furnace will start heating.
- Once the set point temperature has been reached within the programmed time, this temperature will be kept for the programmed waiting time.
- At the end of the process, remove the cylinders and crucibles carefully. While the unit is running the outer housing may heat up. The cylinders and crucibles also heat up considerably, and so gauntlets and tongs should be used when handling them.
- Use the main (C) switch to disconnect the furnace from the mains.

MAINTENANCE AND CLEANING

To ensure long, reliable operation of the unit, proceed as follows:

- Clean the outside of the furnace regularly with a cloth.
- Take care not to damage the muffle lining the inside walls of the furnace hearth when you put in or take out cylinders and crucibles.
- Always use the SiO₂ hearth plate to prevent filtration, which may damage the muffle.
- Prevent liquids or solids from contaminating or dirtying the muffle of the furnace.
- To make best use of the unit, TALLERES MESTRAITUA, S.L. recommends the use of MESTRA cylinders and crucibles.
- The furnace has two circuit breakers located near the cable (J and K) to protect the furnace against electric overloads. To get the furnace working again, reset the circuit breaker.

PRECAUTIONS

- ◆ Before connecting the unit ensure that the power outlet is for 230 V, 50/60 Hz, and that it is earthed. It must also be properly protected by a thermal magnetic cut-out (consult your installer).
- ◆ Do not allow children or unskilled personnel to handle the unit.
- ◆ Always keep the door closed while the furnace is operating.

- ◆ Take care when opening the door: the hot fumes pouring out can cause burns.
- ◆ Do not overfill the furnace or place objects in it which might fall when the door is opened.
- ◆ Do not put inflammable materials or liquids in the furnace, or substances which may give off toxic fumes.
- ◆ Place the furnace on a flat, firm, vibration-free surface close to an extractor hood. Ensure that the furnace is well away from any inflammable materials.
- ◆ Place the unit so that its sides are at least 15 cm from any wall or other object. Take particular care to place the ventilation tube in a free area at least 20 cm from the closest object.
- ◆ Always use gauntlets and tongs to remove and put in cylinders and crucibles. Do not touch the metal housing of the furnace when it may be hot.
- ◆ Handle heated materials removed from the furnace with great care: their temperature may be as high as 1100 °C and may therefore cause accidents.
- ◆ Before using the furnace for the first time, hold it at 300 °C for 1 hour to eliminate any moisture in the muffle.
- ◆ Always use the furnace under an extractor hood.

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