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The R-080410 is a practical pantograph type support especially designed for micro-milling work and to facilitate the handling of its built-in parallelometer. The unit is supplied with the following accessories:

- 4 sets of bushings of different sizes.
- 1 tool holder support.
- 1 cutter holder with wrench.
- 2 charcoal scribes and 1 analyser tip.
- 3 wing tips (0.25, 0.5 & 0.75 mm).
- 1 spatula.
- 2 Allen wrenches (3 & 5 mm).

### **WORKPIECE ATTACHMENT**

1. Release the fixing screw (A)
2. Locate the object to be machined on the work bed (B). This bed has two grooved stops to prevent the workpiece from moving: ensure that the workpiece is properly fitted between them.
3. Tighten the fixing screw (A).
4. The work bed has a ball joint so that it can be set at the optimum inclination for each job. To adjust this joint, slacken the lock nut (C), move the bed as required and re-tighten it.
5. The bed can be dismantled by releasing the two screws fixing it to the support with the 5 mm Allen wrench. This enables the workpiece to be oriented by hand for greater precision.

### **USING THE PARALLELOMETER**

The lower con-rod (D) enables either the pantograph-type micro-milling support or the parallelometer to be brought closer to the work bed, thus preventing any problems due to these devices interfering with each other.

To use the parallelometer, turn the lower con rod (D) until it is perpendicular to you and the minimum gap between the con rod end and the bed is reached. To prevent unwanted movements during handling, the lower con rod can be locked by turning the locking screw beside the turning shaft with a 3 mm Allen wrench.

The parallelometer has two more articulated con rods which provide access for the analyser tip of the instrument to any point on the horizontal plane. Vertical movement is obtained by pressing the opposite end to the analyser tip.

### **USING THE MICRO-MILLING SUPPORT**

- Fit the required micro-motor hand-piece on the micro-milling head (E). To ensure a correct fit, select the most suitable of the four sets of bushings supplied with the unit. To attach the micro-motor hand-piece tighten the head fixing screw.
- Bring the micro-milling support assembly close to the work bed by rotating it on the lower con rod (D). To prevent unwanted movements during micro-milling, the lower con rod can be locked by turning the locking screw with a 3 mm Allen wrench.

The micro-milling support has two more con rods which allow the head to move freely in the horizontal plane. One of the two directions of horizontal movement can be disabled by turning the locking screw beside the con rod turning shafts with a 3 mm Allen wrench.

- The head is suspended from a non-rigid parallelogram which enables it to move vertically. To adjust the level to which it drops, slacken the lock nut on the descent stop (F) and turn the stop screw to the desired level. Then re-tighten the lock nut.
- The head height can also be held steady at the desired level. To do this, lock the articulations of the non-rigid parallelogram by tightening the two turning screws with a 5 mm Allen wrench.

For certain tasks, the cutter holder (H) supplied with the unit can be mounted on the micro-milling support head. This enables a wide range of tips and accessories to be used (consult your MESTRA dealer). The tips not being used can be stored in the tool holder support (G) at the base of the unit.

## MAINTENANCE & CLEANING

To ensure long life and service from the unit, we recommend the following:

- After micro-milling, use a fine brush to clean away dirt and dust deposited on the unit. Take particular care to clean the articulations of the con rods.
- Con rods are articulated on bearings. The articulations of the non-rigid parallelogram on the head can be lubricated regularly if required, using fine mineral oil.
- For general cleaning of the unit, use a dry cloth. Never submerge it in liquid.
- Do not force the con rods to move when they are locked: damage to their shafts could ensue.

## PRECAUTIONS

- ◆ When micro-milling use suitable protective gauntlets, goggles and a face-mask.
- ◆ Do not allow children or unskilled personnel to handle the unit.
- ◆ Before starting work, ensure that the workpiece is properly fixed to the bed.

## TECHNICAL CHARACTERISTICS

Weight:	9 kg
Height:	350 mm
Width:	200 mm
Depth:	230 mm

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