
Ref. 080532

ES.....	5	Manual original
EN.....	10	
FR.....	15	
DE.....	20	



Nº	Referencia	Descripción	Description	Description
1	080530-04	Tubo de aspiración	Aspiration tube	Tube d'aspiration
2	080532-14	Microfiltro de motor	Motor micro filter	Micro-filtre moteur
3	080531-30	Filtro de bolsa	Bag filter	Sac d'aspiration
4	000232	Fusible 5A	Fuse 5A	Fusible 5A
5	080532-22	Conjunto motor de aspiración	Aspiration motor kit	Ensemble moteur d'aspiration
6	080530-27	Escobilla de motor (1 U.)	Motor brush (1 U.)	Charbons (1 U.)
7	000250	Cable de conexión	Connection cable	Câble d'aspiration
8*	080530-26	Bifurcación dos puestos	Two places Y-connector	Bifurcation d'aspiration

* Accesorio opcional / **Optional fitting** / *Accessoire optionnel*

- | | |
|---|---|
| A | Boquilla de aspiración / Suction hole / <i>Entrée d'aspiration</i> / Saugmundstück |
| B | Asa / Handle / <i>Poignée</i> / Griff |
| C | Cierre / Lock / Riegel zur / <i>Fermeture</i> / Deckelöffnung |
| D | Tapa del filtro / Filter lid / <i>Couvercle du filtre</i> / Filterdeckel |
| E | Ruedas / Wheels / <i>Roulettes</i> / Rollen |

1. Piloto de aspiración / **Suction pilot** / *Pilote d'aspiration* / *Leuchtanzeige*
2. Display
3. +/-
4. Indicador de estado del filtro / **Filter state indicator** / *Témoin de remplissage filtre* / *Filterfüllstandsanzeige*
5. Automático-manual / **Automatic-manual** / *Automatique-Manuel* / *Automatik-Manuell*
6. Marcha-paro / **Start-stop** / *Marche-Arrêt* / *Start-Stopp*
7. Interruptor general / **Main switch** / *Interrupteur général* / *Hauptschalter*
8. Tomas de corriente / **Power sockets** / *Prises de courant* / *Versorgungsstecker Schuko*
9. Interruptor térmico / **Thermal switch** / *Interrupteur thermique* / *Thermischer Sicherungsautomat*

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Ignacio Mestraitua
Director Gerente
Sondika, 22.04.2020

The Ref. 080532 EOLO PLUS is a practical aspiration unit specially designed to be used within the field of dental prosthesis. It allows the aspiration of aluminium oxide, silicon oxide, resin dust, plasters, and investments and, in general, any material used in the dental lab. The unit's most important features are the following:

- Three-stage filtering circuit, including a filter-bag, a cartridge-filter and a plate-filter in the clean chamber.
- Indicator of the filter's filling level. "Obstruction/Filter Full" alarm. Long operation alarm, which disconnects the unit after operating continuously for more than three hours.
- Electronic regulation of aspiration power.
- Possibility of "automatic" mode operation, synchronizing the machine's start up with 2 micromotors and a sand blast.
- Possibility of simultaneous operation in two working places.

In order to achieve the unit's optimal performance and a long working life, we strongly recommend to read carefully the following operating instructions.

INSTALLATION

1. Carefully unpack your EOLO PLUS aspiration unit. You will find a Ø 35 mm flexible connecting hose inside the packing.
2. The machine has been designed to operate under the workbench, but it can be placed anywhere else. Choose a firm surface to place it, far from heat sources, damp or vibrations.
3. Open the filter lid (D) and make sure the filter-bag is correctly placed. Close the lid again.
4. Connect the unit's aspiration hose to the aspiration nozzle (A). To do this, simply press it firmly. The other end of the hose can be connected to a micromotor, sand blast, dry trimmer or any other similar device. In case you want to use the aspiration unit in two working places simultaneously, you need a "T junction for two working places (ref. 080530-26)", supplied as an accessory.
5. Connect the power cord to a grounded 230 V, 50/60 Hz power source.
6. The unit is ready to be operated.

OPERATION

The aspiration unit features two operational modes:

- **Manual mode:** The aspiration unit operates continuously monitored only by the [START/STOP] key.
- **Automatic mode:** The aspiration unit operates synchronized with a micromotor. When pressing the pedal the aspiration starts up automatically and continues in operation until about five seconds after releasing the pedal.

The operational mode is selected with the [AUTO/MAN.] key placed at the unit's controls desk. To start up the aspiration unit, proceed as follows:

1. Turn the main switch (7) on.
2. The display (2) will show the aspiration power (from 45 to 99%).
3. Select the operational mode you wish by pressing the [AUTO/MAN.] key.
4. If you select the "manual" mode (pilot light off), the aspiration motor will start up automatically.
5. If you select the "automatic" mode (pilot light on) the machine remains in "Stand-by" until you press the pedal of the micromotor connected to the corresponding base of the aspiration unit. When pressing the pedal of these machines the unit will start up and will continue in operation until about five seconds after releasing the pedal.
6. The motor power can be controlled during the operation of the aspiration unit within a range from 45 to 99% by pressing the [+/-] keys.

Every time there is an obstruction in the aspiration hose or the filter is full, the unit will stop and the red LED pilot (4) will light intermittently. If the unit was operating in automatic mode, the aspiration will not start up until it detects that the filter is in good condition, even if the micromotor is in operation. If the unit was working in manual mode, the aspiration will stop and the LED pilot will cease to light. To start the unit up again, press the [START/STOP] key. If after the "Filter-full" time the unit fails to detect that the filter is in good condition, the aspiration will stop again.

IMPORTANT

The "Obstruction/Filter Full" detection system has been designed to protect the aspiration unit from any eventuality. Under no circumstances should you wait to replace the filter until the detection system operates. Check the filling level of the filter-bag periodically and proceed to replace it (see "Filter and replacement" section in this manual).

The aspiration unit also features a motor protection system, which disconnects the motor after operating continuously for more than 3 hours. In this case, the unit will stop for 1 hour. The aim of this system is to prevent motor damage when the unit, by an oversight, is in operation for a too long period. When this occurs, the aspiration LED pilot (1) will light up intermittently. The unit will not aspirate until the main switch (7) is turned off and on again.

The automatic start-up system allows connecting a sand blast up to 400 W. If this power is exceeded, the thermal switch (9) will disconnect the micromotor. When this happens, wait for about 10 minutes before resetting the thermal switch again by pressing the key.

SELF-CALIBRATION FOR MICRO MOTOR DETECTION

1. To adjust the detection of the micro motor, turn the unit on and press the keys [AUTO/MAN.] and [START/STOP] simultaneously for 5 seconds until the message '1-' is displayed, indicating which outlet is to be adjusted. To select another outlet (2- or 3-) press the key [AUTO/MAN.].
2. Once the outlet with the micro motor to calibrate has been selected, keep the micro motor pedal pressed (spin min.). Press the key [AUTO/MAN.] for 2 seconds to start the self-calibration.
3. If the self-calibration succeeds, the message 'P1', 'P2' or 'P3' will be displayed.
4. Press the key [AUTO/MAN.] for the following outlet (or end of calibration).
5. Press the key [START/STOP] to escape in any moment.

FILTER REPLACEMENT


The EOLO PLUS aspiration unit features a three-stage filtering circuit: the first stage is a recyclable filter-bag; the second stage is a cartridge-type absolute filter; and the third stage is a plate-filter placed in the clean chamber. This latter filter has an unlimited duration, so it should only be replaced occasionally and by the Technical Service. To clean and replace the filters proceed as follows:

Replacing the filter-bag:

1. Open the aspiration unit's upper lid (D) by opening the lock (C). Remove the filter-bag together with the lid of the machine's aspiration chamber. Be extremely careful while performing this operation to prevent dust particles falling into the unit's aspiration chamber.
2. Place the new filter inside the aspiration chamber and make sure that it fits over the ring of the aspiration nozzle (A), and that the filter base rests on the lower tray of the aspiration chamber.
3. Close the unit's upper lid again.

Replacing the cartridge-filter:

1. Open the aspiration unit's upper lid (D) by opening the lock (C) and remove it together with the filter-bag.
2. Release the bolt holding the filter frame to access the cartridge-filter. Remove it carefully avoiding dust into the motor aspiration hole.

3. Place the fresh cartridge-filter and proceed to set again the rest of the components in reverse order.



The cartridge-filter cannot be recycled, so it always needs to be replaced. **It is recommended to replace it every 1500 or 2000 operating hours of the unit.**

MOTOR BRUSHES MAINTENANCE



With your Eolo Plus aspiration unit, you will receive two spare brushes for the motor. The unit has a 16000 rpm motor, so brushes must be replaced periodically. After working 1000 h, the electronic controls will show the message 'br', meaning that it is the time to replace the brushes. Then, proceed as following:

1. Unplug the unit.



2. Release the 6 bolts of the inner lid and remove it. You will get access to the motor.
3. Remove the 4 bolts holding the brushes (10).

4. Remove the brushes and disconnect the cables (11).
5. Disassemble metal parts from plastic parts (12).

6. Replace the new brushes into the plastic parts (13).
7. Connect the cables to the new brushes (14).


8. Use a screwdriver to pull the grey part of the brushes and make them go out (15).
9. Replace the new brushes on the motor and fix them with the 4 bolts (10).
10. Close the inner lid and check the correct working of the unit.
11. In order to reset the brush timer, press simultaneously keys [AUTO/MAN] and [-] for 5 seconds. Notice how the display blinks and 2 seconds later the timer resets and the 'br' message will disappear.

NOTE: When the message 'br' appears, you will have to put the unit off and on in order to run again. The message will keep on to remind that brushes must be replaced as soon as possible.

ASPIRATION MOTOR MAINTENANCE

After working 1600 h, the electronic controls will show the message 'En', meaning that it is the time to replace the aspiration motor. Then, proceed as following:

1. Unplug the unit.
2. Pull the 4 pins at the lower lid (16) and remove it. You will gain access to the motor.

3. Unlock the 2 clamps holding the motor (17).
4. Remove the old motor and place the new one. Do not forget the cable connection.
5. Hold the motor with the 2 clamps.
6. Close the inner lid and check the correct working of the unit.
7. In order to reset the motor timer, press simultaneously keys [AUTO/MAN] and [+] for 5 seconds. Notice how the display blinks and 2 seconds later the timer resets and the 'En' message will disappear.

NOTE: When the message 'En' appears, you will have to put the unit off and on in order to run again. The message will keep on to remind that motor must be replaced as soon as possible. If not, after 100 hours the unit will stop until you replace the motor.

MAINTENANCE AND CLEANING

In order to achieve the machine's optimal performance and a long working life, we strongly recommend you:

- To clean the outer casing of the unit by always using a cloth. Do not use solvents or any other inflammable product.
- To check the filling level of the filter-bag regularly and to replace it according to the instructions in this manual. **Do not wait until the "Obstruction/Filter full" detection system operates.**
- To replace the filters periodically. The cartridge-filter has an estimated working life of 1500/2000 hours. Always use MESTRA original spare parts.
- To perform all the filter replacement operations very carefully in order to avoid dirtying the inside of the aspiration chamber. Make sure that the filter bag is correctly fitted into the aspiration nozzle (A) and that it rests on the lower tray of the aspiration chamber.
- In case of overload in a Schuko power outlet (8), the thermal switch (9) will go off. When this occurs, wait for several seconds and set it again by pressing the key.
- Make sure the aspiration hoses installed in your lab do not have losses due to cracks, breaks or faulty connections.
- For the aspiration motor's sake, the unit should not operate continuously for long periods (more than 2-3 hours). In any case, the unit features a protection system for this eventuality.

PRECAUTIONS

- Before connecting the machine, make sure that the power source is a grounded 230 V, 50/60 Hz.
- Do not connect micromotors or other devices with a power consumption of more than 400 W in the Schuko power outlets.
- Do not allow children or unskilled personnel to operate the machine.
- Choose a rigid and well-balanced surface to place the aspiration unit.
- Do not open the upper lid when the unit is in operation.
- Disconnect the machine from the power source during the cleaning and filter replacement operations.

TECHNICAL SPECIFICATIONS

Voltage:	230 V, 50/60 Hz
Power:	900 W
Height:	585 mm
Width:	200 mm
Depth:	485 mm
Weight:	17 kg

DECLARATION OF CONFORMITY UE

The firm **Talleres Mestraitua S.L.**

with social address: Txori-erri Etorbidea, 60
48150 Sondika - Bilbao - Spain

Hereby declares that the product

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Meets the requirements set out in the following directives:

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GUARANTEE

All **MESTRA®** units are guaranteed against every defect in manufacturing for a period of one year from date of purchase. Prior to honouring a guarantee, **MESTRA®** will require that the attached coupon be filled in with date of sale, duly stamped by the distributor.

The guarantee does not cover breakdowns caused by incorrect installation or use.

MESTRA® will not effect any repair outside its own workshops and will require transport costs.

Repairs carried out on the premises of the customer, be they with or without guarantee, are entrusted to the technical assistance of the distributor who has sold the unit or alternatively to the technicians he has contracted

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