R-080246

E.....	CHORREADORA DE ARENA “PLANETARIUM PLUS IV”4
GB.....	“PLANETARIUM PLUS IV” SANDBLASTER7
F.....	MACHINE À SABLER “PLANETARIUM PLUS IV” 10
P.....	JACTO DE AREIA “PLANETARIUM PLUS IV” 13
I.....	SABBIATRICI “PLANETARIUM PLUS IV” 16
D.....	“PLANETARIUM PLUS IV”-SANDSTRAHLER 19

MESTRA®

TALLERES MESTRAITUA S.L.



Txori-Erri Etorbidea, 60

Tfno. + 34 944530388 - Fax + 34 944711725

E-mail: mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA


Rev. 23/02/2010

A	Vaso depósito de esmeril / Emery deposit vessel / Vase réservoir d'émeri / Vaso depósito de esmeril / Vaso deposito smeriglio / Schmirgelbehältergefäß
B	Manorregulador del vaso / Vessel pressure control / Régulateur de pression du vase / Regulador da pressão do vaso / Regolatore pressione vaso / Gefäßdruckregler
C	Interruptor luz cabina / Booth light switch / Interrupteur éclairage cabine / Interruptor luz cabina / Interruttore luce cabina / Kabinenlichtschalter
D	Tubo de entrada aspiración exterior / External inlet tube / Tuyau d'entrée aspiration extérieur / Tubo de entrada da aspiração exterior / Tubo d'entrata aspirazione esterna / Einlassrohr äußere Absaugung
E	Sufridera de goma / Rubber guard / Plaque de choc en caoutchouc / Placa de choque de borracha / Tassello di gomma / Gummidruckplatte
F	Interruptor sistema de calefacción de vasos / Vessel heat system switch / Interrupteur système de chauffe des vases / Interruptor do sistema de aquecimento dos vasos / Interruttore sistema di riscaldamento vasi / Schalter Gefäßaufwärmesystem
G	Conmutador de selección de vaso / Vessel selection switch gear / Commutateur de sélection de vase / Comutador de selecção do vaso / Commutatore di selezione vaso / Gefäßwahlumschalter
H	Manguitos de tela / Fabric sleeve / Manchons de toile / Juntas de tecido / Manicotti di tela / Gewebesschläuche

1.	Interruptor / Main switch / Interrupteur / Interruptor / Interruttore / Schalter
2.	Aire comprimido / Air / Air comprimé / Ar comprimido / Aria compressa / Druckluft
3.	Pedal / Pedal / Pédale / Pedal / Pedale / Pedal
4.	Enchufe / Plug / Prise / Tomada / Presa / Stecker

IMPORTANTE
IMPORTANT
IMPORTANT

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
The Planetarium Plus IV R080246 is a sand blaster specially designed for micro-sandblasting small articles of jewellery and dental prostheses. To get the most out of these apparatuses, carefully read the instructions below.

INSTALLATION

- Carefully unpack the apparatus. Inside you will find the following items:
 1. Quick-fit mains compressed air connection with 6 mm diameter tube.
 2. 90° PVC elbow-bend for fitting on exhauster outlet tube if necessary.
 3. Pedal for operating the apparatus.
 4. Plastic protection screens (2 units).
- Set up the apparatus on a hard, flat surface away from sources of vibration and heat, in an atmosphere that is not too humid. Excess humidity can cause particles of grinding medium in the booth to compact, causing them to clog the micro-sandblasting nozzles. Ensure the unit is at a height where objects inside the booth can be handled comfortably.
- Position the plastic protector inside the glass.
- Connect one end of the quick-fit air connector provided with the unit to the inlet connector on the rear panel of the booth (4) and the other end to a compressed air outlet working at a minimum of 7 bar and a flow rate of at least 100 l/min.
- Plug the electrical lead (2) into an earthed 230 V 50/60 Hz power outlet.
- Press the main switch (1)
- Connect the drive pedal tube to the nozzle on the left rear of the unit (3). Then place the pedal within reach of your foot, in a position where it can be worked easily. This pedal can also be set up to be worked with the knee if this is more comfortable, or if considerations of space make it advisable.
- Remove the threaded stoppers from the micro-sandblasting modules (A) and fill the vessels with grinding medium. Never fill beyond the maximum level marks on each vessel. Screw the stoppers back on and ensure that the seal is leak tight.
- The unit is supplied with colour coded nozzles in the following diameters:

Yellow nozzle	Ø 0.75 mm	
Blue nozzle	Ø 1 mm	
Green nozzle	Ø 1 mm +	
Red nozzle	Ø 1.25 mm	

PRESSURE REGULATION

The blowing pressure of the sandblasting guns can be set to suit the nature and characteristics of the work to be done by adjusting the pressure regulator (B) for each vessel. Pull the pressure regulator knob up and turn it to select the pressure desired. Each vessel has its own pressure gauge, so working pressure can be selected easily. Once the pressure is set, push down the knob to lock it in place.

CONNECTING THE EXHAUSTER

The micro-sandblasting unit must always be used with a forced-draught exhauster. Outlet tube (J) on the left hand side of the unit is fitted for this purpose. The exhauster must be capable of creating a slight vacuum pressure inside the booth so that dust particles move smoothly and continuously towards the outlet tube.

The PVC elbow-bend supplied with the machine can be used if the shape of the external exhauster makes it necessary.

SAND HEATING SYSTEM

The vessels of the Planetarium sandblaster are electrically heated to maintain sand at around 40-50 °C. This eliminates damp and prevents the sand from compacting and clogging the unit. This device is especially recommended for work in humid atmospheres.

To start up the heater just connect the switch on the top of the booth (F). An automatic regulator will then maintain the temperature. The lamp on the switch burns steadily when the heater is on. When connecting the automatic sand heater remember that the pneumatic parts of the unit will expand during the early stages of heating until they reach thermal equilibrium. **During this period there may be whistling or other noises in the electro-valves as they operate. These noises do not affect the operation or service life of the unit, but may be a nuisance to users.** We therefore recommend that the vessel heating system be connected at least half an hour before work is due to begin.

USING THE APPARATUS

These recommendations will help you obtain the best results when micro-sandblasting:

1. If the time of year or the conditions at your laboratory make it advisable to use the sand

- heater, connect it at least half an hour before starting work.
2. Open the booth cover and put the workpiece to be handled inside. Workpieces can also be put in through the openings in the bushings. Check that there is enough grinding medium in the vessels.
 3. Close the booth cover and ensure it is leak-tight (with the glass properly fitted into the frame).
 4. Put on protective goggles, gauntlets and a face-mask as necessary.
 5. Start up the forced draught exhauster. If the workpiece is not clearly visible, connect switch (C) to turn on the light in the booth.
 6. Use the main switch (G) to select the vessel to be used. The colour coding system used makes it easy to see which is the active pan at all times.
 7. Pick up the moveable gun selected from inside the booth (the colours of the gun nozzle holders match the colours on the vessel stoppers, i.e. red, yellow, green or blue). Then press the drive pedal to start the sand flow.
 8. Direct the sand jet carefully: if it hits the booth walls directly, it could damage them. If you need to check the pressure or force of the jet, direct it against the rubber guard (E) on the left of the booth interior.
 9. When you finish sandblasting, release the pedal to stop sand flowing through the nozzles. The pneumatic discharge mechanisms will cut in immediately to relieve the pressure in the vessels. Stop the exhauster and open the booth cover to remove the workpiece.

MAINTENANCE & CLEANING

To ensure the apparatus has a long working life, we recommend the following:

EXTREMELY IMPORTANT
The silencer filter at the bottom of each vessel must be changed every 6 months to prevent perforations.

- Always use top quality grinding medium. Talleres Mestraitua, S.L. recommends MESTRA grinding media. 1.25 or 1 mm diameter nozzles for micro-sandblasting guns work with grinding media up to 250 and 150 micron in diameter, respectively. 0.75 mm guns use media with diameters up to 75 microns (consult the table).
- Dry piston air compressors are preferable. Oil particles in the compressed air will create

lumps in the sand which may clog the gun nozzles.

- Never fill the tanks with grinding material residues from the booth floor: they may contain large-diameter particles which could clog the nozzles, and are contaminated with particles blasted off workpieces.
- Never direct the jet of grinding media from the guns directly at the booth walls, as the resulting abrasion causes wear.
- Clear the grinding media out of the booth regularly via the drainage stopper inside the booth.
- Regularly bleed off any condensation which builds up in the air inlet filter. The bleed control is easily accessible on the rear of the unit.
- Clean the apparatus with a dry cloth or a cloth slightly moistened in soapy water.

If the guns clog up, remove the nozzle holder and remove the nozzle from inside. Clear it by blowing against the normal air flow direction or by inserting a pin or a thin wire, then fit it back into the gun.

The table below shows the maximum diameters of the grinding media for each nozzle diameter.

MAXIMUM DIAMETERS OF THE GRINDING MEDIA FOR EACH NOZZLE DIAMETER	
Nozzle	Maximum diameter
1.25 mm	250 microns
1 mm	150 microns
0.75 mm	75 microns

PRECAUTIONS

- ◆ Never use pressures of more than 5 bar. Optimum working pressure is between 2 and 4 bar.
- ◆ Always work with the booth cover closed, and wear protective clothing, including a face-mask and goggles.
- ◆ Connect the booth outlet tube to an exhauster.
- ◆ Immediately replace any broken protective bushings.
- ◆ Do not allow children or unskilled personnel to use the apparatus.
- ◆ Only use top quality grinding media: we recommend MESTRA media. Low quality grinding media may be toxic if inhaled and may clog the compressed air nozzles.
- ◆ Connect the power cable to a 230 V, 50/60 Hz earthed power outlet.
- ◆ Ensure that the compressed air system used to power the apparatus has at least one protective pneumatic discharge valve.

- ◆ Consult your installer immediately if leaks appear in the compressed air circuit.
- ◆ Before starting work, ensure that the stoppers of the micro-sandblasting modules are properly closed.
- ◆ Never work with the booth door open.

TECHNICAL CHARACTERISTICS

Height:	385 mm
Length:	515 mm
Width:	440 mm
Weight:	19.5 kg
Working pressure:	2-6 Bares
Voltage:	230 V, 50/60 Hz
Power consumption:	115 W

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