**Ref. 080510**

ES.....	MEZCLADORA DE GELATINA .....	6
EN .....	GEL MIXER .....	12
FR.....	MELANGEUR DE GELATINE .....	18
DE .....	GELATINEMISCHGERÄT .....	24
PO .....	MISTURADORA DE GELATINA.....	30
IT .....	MESCOLATORE DI GELATINA .....	36

**MESTRA®**

**TALLERES MESTRAITUA S.L.**

Txori-Erri Etorbidea, 60

Tel. +34 944530388 - Fax +34 944711725

mestra@mestra.es - www.mestra.es


**48150 SONDIKA - BILBAO - ESPAÑA**

Rev. 26/05/2014


## Ref. 080510

- |   |  |
|---|--|
| A | Tapa superior / <b>Upper lid</b> / <i>Couvercle</i> / Oberer Deckel / <b>Tampa superior</b> / <i>Coperchio superiore</i>                         |
| B | Panel de control / <b>Control panel</b> / <i>Panneau de contrôle</i> / Bedientafel / <b>Painel de controlo</b> / <i>Pannello di controllo</i>    |
| C | Patas de asiento / <b>Legs</b> / <i>Pieds</i> / Stützfüße / <b>Patas de encaixe</b> / <i>Piedini d'assetto</i>                                   |
| D | Interruptor general / <b>Main switch</b> / <i>Interrupteur général</i> / Hauptschalter / <b>Interruptor geral</b> / <i>Interruttore generale</i> |
| E | Caño / <b>Spout</b> / <i>Bec verseur</i> / Auslaufstutzen / <b>Bico de descarga</b> / <i>Tube</i>  |
| F | Grifo de vertido / <b>Tap</b> / <i>Robinet</i> / Ablasshahn / <b>Torneira de despejo</b> / <i>Rubinetto di scarico</i>                           |

1. Display / **Display** / *Cadran digital* / Display / **Display** / *Display*
2. Piloto de resistencia encendida / **Resistance light on** / *Temoin de fonctionnement de la resistance* / Anzeigelämpchen Widerstand leuchtet / **Lâmpada de sinalização de resistência acesa** / *Spia di resistenza accesa*
3. Piloto de temperatura °C1 encendida / **°C1 temperature light on** / *Temoin de température °C1* / Anzeigelämpchen Temperatur °C1 leuchtet / **Lâmpada de sinalização de temperatura °C1 acesa** / *Spia di temperatura °C1 accesa*
4. Piloto de batidora en marcha / **Mixer on light** / *Temoin de fonctionnement du malaxeur* / Anzeigelämpchen Mischgerät in Betrieb / **Lâmpada de sinalização de batedor em funcionamento** / *Spia di mescolatura avviata*
5. Tecla de incrementar / **Increase light** / *Touche +* / Taste Erhöhen / **Tecla de incrementar** / *Tasto incrementare*
6. Tecla de decrementar / **Decrease light** / *Touche -* / Taste Verringern / **Tecla de diminuir** / *Tasto diminuire*
7. Tecla de reloj / **Watch key** / *Touche de l'horloge* / Uhrentaste / **Tecla de relógio** / *Tasto orologio*
8. Tecla de selección de temperatura °C1 / **°C1 temperature selection key** / *Touche de selection de la temperature °C1* / Wähltaste Temperatur °C1 / **Tecla de selecção de temperatura °C1** / *Tasto di selezione temperatura °C1*
9. Tecla de selección de temperatura °C2 / **°C2 temperature selection key** / *Touche de selection de la temperature °C2* / Wähltaste Temperatur °C2 / **Tecla de selecção de temperatura °C2** / *Tasto di selezione temperatura °C2*
10. Piloto de arranque retardado programado / **Programmed delayed startup light** / *Temoin de programmation differee* / Anzeigelämpchen programmierter verzögerter Anlauf / **Lâmpada de sinalização de arranque retardado programado** / *Spia di avvio ritardato programmato*
11. Piloto de fase de calentamiento iniciada / **Heating phase initiated light** / *Temoin de programme de chauffe °C1* / Anzeigelämpchen Aufwärmphase gestartet / **Lâmpada de sinalização de fase de aquecimento iniciada** / *Spia fase di riscaldamento iniziata*
12. Piloto de fase de enfriamiento iniciada / **Cooling phase initiated light** / *Temoin de programme de refroidissement jusqu'a °C2* / Anzeigelämpchen Abkühlphase gestartet / **Lâmpada de sinalização de fase de esfriamento iniciado** / *Spia fase di raffreddamento iniziata*
13. Interruptor general / **Main switch** / *Interrupteur general* / Hauptschalter / **Interruptor geral** / *Interruttore generale*

Nº	Ref.	DESCRIPCIÓN	DESCRIPTION
0	080510-01	Carcasa	Housing
1	080510-34	Aspas	Lever
2	080510-35	Retén de eje	Shaft retainer
3	080510-36	Casquillo de eje	Shaft bushing
4	080510-37	Porta aspas	Lever holder
5	080510-38	Junta de eje	Shaft gasket
6	080510-39	Eje	Shaft
7	080510-40	Cuba	Vat
8	080510-42	Sonda de cuba	Vat probe
9	080510-43	Junta aislante	Insulating seal
10	080510-33	SopORTE motor	Motor support
11	080510-24	Relé 230V	Relay 230 V
12	080490-37	Tuerca M3	Nut M3
13	080510-26	Zócalo relé	Relay support
14	080510-67	Tornillo relé	Relay bolt
15	080460-65	Tornillo motor	Motor bolt
16	080522-09	Arandela motor	Motor washer
17	080511-10	Motor	Motor
18	080510-91	Clixon	Clixon
19	080510-54	Sonda de resistencia	Resistor probe
20	080510-53	Resistencia	Resistor
21	080510-47	Casquillo motor-eje	Motor-shaft bushing
22	080510-49	Espárrago M6x6	Stud M6x6
23	080510-98	Espárrago punta cónica M6x6	Stud conic point M6x6
24	080075-08	Tornillo M3x16	Bolt M3x16
25	080508-21	Regleta de conexiones	Terminal strip
26	080010-05	Pasacables	Cable bushing
27	080510-07	Espárrago M5x10	Stud M5x10
28	080550-23	Tornillo M3x8	Bolt M3x8
29	080490-36	Arandela elástica M3	Spring washer M3
30	080592-26	Fuente de alimentación 24V	Power source 24V
31	080510-41	Junta de grifo	Tap gasket
32	000236	Fusible 7A	Fuse 7A
33	080490-06	Base de conexión	Connection base
34	080400-08B	Ventilador	Fan

Nº	Ref.	DESCRIPCIÓN	DESCRIPTION
35	080510-66	Junta goma tapa	Lid rubber gasket
39	080080-22	Tuerca ciega M5	Cap nut M5
40	080510-56	Bisagra superior	Upper hinge
41	080510-02	Tapa	Lid
42	080510-05	Bisagra parte inferior	Lower part hinge
43	040410-04	Tuerca bisagra	Hinge nut
44	080510-07	Espárrago M5x10	M5x10 stud
45	080510-28	Protector sensor	Sensor protector
46	080510-27	Sensor	Sensor
47	080510-10	Junquillo	Tongue
48	080510-14	Junta tórica grifo	Tap O-ring
49	080510-15	Porta juntas grifo	Tap seal housing
50	080510-16	Tuerca eje de grifo	Tap shaft nut
51	080510-17	Eje del grifo	Tap shaft
52	080510-18	Cuerpo del grifo	Tap body
53	080510-19	Casquillo del grifo	Tap bushing
54	080510-20	Tuerca del grifo	Tap nut
55	080510-21	Pomo del grifo	Tap knob
56	080510-22	Pila CR2032 (3V)	Battery CR2032 (3V)
57	0800510-23	Control	Control
58	080400-10	Interruptor	Switch
59	080080-26	Patas	Legs
60	080251-09	Tornillos de patas	Legs screws
62	000250	Cable	Cable
63	080510-31	Llave para aspas	Blades removing tool
64	080510-03	Tornillo empujador	Pusher screw
65	080510-04	Pulsador	Pusher
66	080510-05	Bisagra parte inferior	Lower part hinge
67	040410-04	Tuerca bisagra	Hinge nut
68	080510-07	Espárrago M5x10	M5x10 stud
69	080510-08	Muelle de empujador	Pusher spring
70	080510-09	Tuerca de empujador	Pusher nut
71	040060-60	Microinterruptor	Micro switch
72	080510-13	Tornillo microinterruptor	Micro switch screw

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The GEL BAT-6 (Ref. 080510) is a gel mixing practise specially designed for use within the mechanical dentistry sector. Safety elements complying to the most rigorous guidelines of EC regulations have been used in its design. The GEL BAT-6 enables the gel placed in the container to a heating and stirring lea defined by the user which consists of three phases:

1. A prior waiting period, during which the unit remains on standby awaiting the heating phase to commence.
2. A heating period, in wich the gel temperature rises to a value of (°C1) set beforehand by the user, between 65 °C and 95 °C. During this phase heating temperature is higher than the subsequent pouring temperature. The idea is to achieve a homogeneous gel mixture and without any clots. Also, this phase the mixer blades are set in motion in order to further favour the correct mixing of the gel.
3. Once the heating phase the gel achieves an optimum texture, but its temperature is too high for pouring. A cooling period then commences forced by a fan which ends when the gel temperature reaches a value of (°C2) between 35 °C and 65 °C preset by the user. Once the third phase has been completed, the gel is then ready for pouring.

### INSTALLING THE UNIT

- Carefully unpack the gel mixer GEL BAT-6.
- Choose a flat surface to place it on, far from any heat or vibration sources. We also recommend choosing a location which is sufficiently ventilated.
- Be sure to leave a space between the surface of the table and the legs (C), given that the cooling air from the container is collected
- Conect the unit to a 230 V, 50/60 Hz electric outlet with an earth connection.
- Proceed to set the control panel on time and programme it according to the indications in the following points.

### SETTING THE CLOCK ON TIME

The GEL BAT-6, is provided with a clock in real time. Every time the main switch (13), the display (1) alternatively shows the real time of the unit and the temperature inside the container. In order to set the time, proceed in the following manner:

1. Activate the main switch (13) pressing the clock key (7) at the same time. For one minute the display (1) will show the real time of the unit and alternatively the word "horA".
2. During this period of time it is possible to modify the time shown. In order to do this, activate the increase (5) and decrease (6) keys until the clock shows the desired time.
3. Press the clock key (7) to save the time.
4. The GEL BAT-6 is delivered from the factory with Spain time.

### PROGRAMMING THE UNIT

The unit can be programmed to execute a cycle made up of the three phases indicated at the beginning of the instruction manual. The (°C1) heating and (°C2) pouring temperatures may be regulated at will by the user between values of 65 °C and 95 °C for (°C1), and 35 °C and 65 °C for (°C2).

The temperature setting is carried out in the following manner:

1. In order to regulate the heating temperature (°C1), press key (8) briefly. You will notice how the heating temperature light on the graph (11) starts to blink, and that display (1) indicates the value set for temperature (°C1).
2. In this situation use the keys to increase (5) decrease (6) until the display shows the desired temperature.
3. Then, briefly press key (°C1) in order to memorize the temperature, you have just selected.
4. In order to regulate the pouring temperature (°C2), press key (9) briefly. You will notice how the heating temperature light on the graph (12) starts blinking, and the display (1) indicates the value set for the temperature (°C2).
5. In this situation use the increase (5) and decrease (6) keys until the display (1) shows the desired temperature.
6. Once the setting process has been carried out, press key (°C2) (9), in order to memorize the temperature, you have just selected.

As regards the periods, the only time wich can be programmed by the user is the time in which he/she wishes to reach heating temperature (°C1). In this sense, the GEL BAT-6 may carry out two types of startups:

On the one hand direct startup- The unit is put into operation directly, in order to reach °C1

temperature in the shortest possible time. In order to carry out a startup of this type it is sufficient to briefly press key (8) °C1. Then press key (8) (°C1) for a few seconds. You will notice that light (11) on the control graph blinks intermittently. If you wish to consult the closest time in which the °C1 temperature will be reached just press the watch key (7) briefly. The display will show the word “horA” and the said hour alternatively.

The other type of startup is delayed- In this case the user inserts the time he/she wishes the °C1 temperature to be reached. Then, the unit remains on standby during a variable amount of time, defined by the control module, and it is automatically turned on at the precised time in order to reach °C1 temperature at the time desired by the user.

There are also two possibilities with the delayed startup procedure:

1. **Specific delayed startup-** This is used when the user wishes the delayed startup to be carried out only once. In subsequent startups, the value of the time to reach °C1 inserted by the user, will be lost.

In order to programme a delayed startup of the type press the clock key (7). You will notice how the display (1) shows the nearest time which according to the microprocessor's calculations the heating temperature (°C1) can be reached. Alternatively, the word “horA” appears. Press the increase (5) and decrease (6) keys until you set the time you wish. When you have finished, press the key (7) before 3 seconds have transpired, in order to set the time and commence the programme.

If when inserting the time, the minimum period the microprocessor deems necessary for heating decreases, the control will interpret that the programmed time is for the following day.

2. **Memory delay start-up** – This option provides two possibilities, as follows: firstly, when we want the value entered by the user so that the machine reaches temperature °C1 to remain for an indefinite period of time in the memory, and to be used by the control function as a default value, i.e. so that temperature °C1 is obtained every day at the same time.

In order to programme a delayed startup of the first type press key (7) twice. You will notice the display (1) alternatively shows a time value, and the word “hPro” (Hora programada). This is the value stored in the unit's memory as a programmed time to reach temperature (°C1). Press the increase

(5) and decrease (6) keys until the display shows the time you wish to reach temperature (°C1). Also, in this case, if the minimum time the microprocessor deems necessary for heating decreases, the control will interpret the programmed time is for the following day. When you have set the time, keep the clock key (7) pressed for several seconds. From this moment on, the unit will startup automatically when required — according to the control calculations— in order to reach temperature °C1 at the programmed time. The time value inserted will be stored in the memory. It is possible to check that the programming has been carried out properly by briefly pressing the clock key (7). The display shows alternately the time programmed and “today”.

During the unit operation, the blades and the resistances in the container will startup automatically each time it is required. In order to indicate the operation of these elements, lights (4) and (2) respectively are kept on.

Once the heating phase is completed, the final cooling phase commences to reach the pouring temperature (°C2). The resistance is disconnected and a forced ventilation commences. The time needed for cooling up to temperature (°C2) is difficult to calculate exactly, because it depends on many factors, among these the specific heat of the gel used, the amount of the gel used, the amount of the gel in the container, the room temperature, the values of (°C1), (°C2), etc.

For this reason, we have decided to limit the scope of the programming time only up to the temperature (°C1), allowing the cooling to proceed without any type of programming. Nevertheless, experimental value show that a reasonable period for cooling the gel from temperature (°C1) to (°C2), may be from 2 to 4 hours.

Once this value is known, it is good practice to programme the time in which we wish to reach the GEL BAT-6 heating temperature (°C1), from 2 to 4 hours before pouring. For example, if we wish to pour at 17:30 PM the unit may be programmed to reach (°C1) at 13:30 PM. In this manner, after four hours cooling, at 17:30 the gel will have reached pouring temperature.

Finally, mention that at any time of the cycle a direct displacement of the (°C1) temperature. In order to do this press key (°C1) (8) for a few seconds. You will then notice that the indicator lamp (11) remains on continuously demonstrating this fact.

It is also possible to carry out a direct displacement up to pouring temperature doing

the same thing on key (°C2) (9). However, a direct startup of the unit at the (°C2) step is not advisable.

### OPERATION

1. Activate the main switch (13). You will notice how the display (1) on the control panel of the unit is activated indicating the GEL BAT-6 is ready to operate.
2. Open the upper lid (A). You will notice that by opening the lid the display shows the message “**PUER**” and the temperature of the gel mixer alternatively.
3. Fill the container with gel. In order to do this, cut the solid block of gel into small pieces of about 4 cm long, and place them in the container (do not cut up the gelatine inside the vat). The final level of the gel should be about 3 cm, above the last blade. When the unit enters into operation, the level of the gel in the container will fall somewhat due to the heating resistance and blades. Then refill the container again until the final level of the gel covers the last blade. (The one closest to the top edge of the container).

**WATCH OUT!! Do not let solidify the jelly into the vessel: when the machine starts up, the temperature sensor could take a wrong value and consequently a wrong machine perform. If this happens, full empty the vessel, cut the jelly in small pieces and start up again the machine.**

4. Close the lid and programme the GEL BAT-6 following the previously indicated instructions. The Ref. 080510 gel mixer will startup automatically if a delayed startup has not been programmed. During the operation the heating resistance light (2) will turn on a stirrer (4) will startup in order to indicate to the user that the programme continues to run without any problems. If a delayed startup were programmed, the clock light (10) would turn on indicate this circumstance. Subsequently, the unit would commence operating at the time foreseen by the microprocessor. (Delayed startup).
5. You can alter the values of the °C1 and °C2 temperatures at any time during the programme execution. On the other hand, if you wish to achieve a gel mix press the (°C1) key for several seconds, because during this step always causes a stirring of the mixture. In order to know what step of the programme is being executed at a given time you must look at the lights on the graph (10), (11), (12).

6. When the heating temperature (°C1) is reached, the GEL BAT-6 will turn on the control panel. The cooling phase will then commence until the pouring temperature (°C2) is reached. You will notice that during this phase light (12) turn on.
7. When the pouring temperature (°C2) is reached, the gel is ready to use. Place a container under the spout (E) and pull the tap knob (F) towards you. The gel will flow out continuously.
8. Once the work with the unit has been completed, empty the unit container by pouring all of the excess gel into another container through the spout (E).
9. In order to finish the programme, press key °C2 (9) continuously for several seconds. You will notice how the three lights on the control graph start to blink. Then disconnect the main switch (13).

During the normal running of the programme the lid of the unit cannot be opened. A safety switch will disconnect the GEL BAT-6 if the upper lid (A) is opened, and will alternatively show an error message on the display (“**PUER**”) along with the container temperature.

Sometimes, although we may have a preset time to reach the °C1 temperature, the control module may detect a specific overheating of the gel, and as a preventive measure alter the control algorithm of the heating resistance, causing a delay of 10 to 15 minutes of the estimated time to reach °C1. This manoeuvre has been foreseen by TALLERES MESTRAITUA, S.L. and should be understood as being completely normal.

If the electric supply is cut off, the gel mixer Ref. 080510 reacts in two different ways: If the cut off is no longer than 10 minutes, the unit understands that it is a once-off failure and keeps the execution step parameters in its memory, so that when the electric power returns, it can continue with its programme from the moment is stopped.

If on the contrary the cut off period is more than 10 minutes, the display (1) shows the message electricity supply failure “**Fred**”. The control then interprets that the gel may solidified and that the stirrer shaft may end up blocked. If the gel has solidified, the unit will carry out a type of special startup (slow) showing the word “**Lent**” on the screen.

During the slow startup, the control executes a series of manoeuvres in the heating resistance and the stirrer motor, trying to return to the normal situation as soon as possible. TALLERES MESTRAITUA S.L. has foreseen this circumstance in the design phase of the

unit, providing the GEL BAT-6 with a high startup torque motor, and a specific control module software. However, in certain occasions the gel may spatter or spill over during a slow startup. Therefore, we recommend being very careful in these situations.

When a programme is being implemented, never turn off the unit directly by means of the main switch (13), given that the unit will consider it a failure in the electric supply grid and when we turn it on again it will show us the word “Fred”, if more than 10 minutes have gone by.

## PROGRAMMING EXAMPLES

### Example 1 (Direct Startup)

Let us suppose we find the unit turned off and we want to carry out a direct startup. As you will remember, it is a unit startup until the temperature (°C1) is reached in the shortest possible time, and a subsequent cooling to temperature (°C2). We also suppose that values (°C1) and (°C2) have been previously inserted by the user:

- Turn on the main switch (13).
- Press key (°C1) (8). Light °C1 (11) will start to blink.
- Press key (°C1) (8) again for several seconds. Light °C1 (11) will turn on permanently.

### Example 2 (Specific delayed startup)

Let's suppose that it is 10:00 AM, and we want to reach the gel the temperature (°C2) of 50 °C by 16:00 PM. We also want temperature (°C1) to reach 85 °C. The programming procedure will be carried out in the following manner:

- Press temperature key (°C1) (8).
- Press the increase (5) and decrease (6) keys until the display shows a temperature of 85 °C.
- Press the temperature key (°C1) again.
- Press temperature key (°C2) (9).
- Press the increase (5) and decrease (6) keys until the displays shows a temperature of 50 °C.
- Press temperature key (°C2) again.
- Press the clock key (7) once.
- As we want the gel to be ready for pouring at 16:00, it is convenient to programme the unit to reach (°C1) 4 hours beforehand, that is, at 12:00. In order to do this press the increase (5) and decrease (6) keys until the display programmes 12:00.
- Press the key (7) once again for a few seconds.

From this moment on, the unit will startup automatically at the precise time in order to

reach temperature °C1 at 12:00. After the cooling period until °C2, the gel will be ready for pouring.

### Example 3 (Memorized delay startup)

Let's suppose we wish to permanently memorize a start up at 8:00 AM, the gel will be ready to pour at a temperature of 55 °C. Furthermore, for a correct mixture a temperature of 80 °C (°C1) has to be reached. Proceed as follows:

- Press temperature key (°C1) (8).
- Press the increase (5) and decrease (6) keys until the display shows a temperature of 80 °C.
- Press the temperature key (°C1) (8) again.
- Press temperature key (°C2) (9).
- Press the increase (5) and decrease (6) keys until the display shows a temperature of 55 °C.
- Press the temperature key (°C1) (9) again.
- Given that you want the gel to always be ready at 8:00 AM, we will have to programme a memorized delayed startup. In order to do this, you press the clock key (7) twice. The display will alternatively show the programme memorized “hPro”.
- We said the pouring temperature was reached 2 to 4 hours after reaching (°C1). Therefore, we recommend programming (°C1) at 4:00 AM in order for the gel to be ready at 8:00. Press the increase (5) and decrease (6) keys until the display shows 4:00.
- Then press the clock (7) for several seconds. The delayed startup light (10) will turn on continuously.

At night, the machine will startup at the time calculated by the microprocessor in order to reach the heating temperature (°C1) by 4:00 AM. Then the cooling phase commences so that the gel is ready for pouring at 8:00 AM.

## MAINTENANCE AND CLEANING

In order to guarantee optimum reliability and long working life of the unit, we recommend you bear in mind the following:

- Cleaning the outer casing of the GEL BAT-6 should always be done with a cloth moistened in soapy water. Do not use solvents or any other inflammable product.
- The control has a lithium battery. It should be changed approximately every two years. When the battery has no further power, the display shows the word “PILA”. The location of the battery appears on the assembly sheet.
- The unit is provided with two thermoprobes. If these were to cause any problems, the

display would show “rot1” or “rot2”. Immediately contact the **MESTRA** technical service.

The motor reducer which drives the blades of the stirrer is protected by an automatic disconnecting switch, and by the microprocessor's own programming software, which disconnects the motor whenever overloads occur.

During the unit's operation, it is perfectly normal for the metal casing to reach temperatures close to 50 °C, a value understood to be correct, with the operations limits foreseen by TALLERES MESTRAITUA, S.L.

All the units are checked before leaving the factory. During the rigorous final quality control carried out in TALLERES MESTRAITUA, S.L. several tests are made using gel, in order to better simulate the true use conditions of the unit. For this reason, it is completely normal to receive your GEL BAT-6 with some remains of gel inside the container.

### PRECAUTIONS

- ◆ **Before carrying out any maintenance on the unit, switch it off and unplug from the power supply.**
- ◆ Before connecting the unit make sure the electric connection is at 230 V, 50/60 Hz with an earth connection.
- ◆ Neither children nor unqualified personnel are allowed to manipulate the unit.
- ◆ Do not open the upper lid when the unit is in operation. Under circumstances must you bridge the safety switch which protects the opening of the lid.
- ◆ Avoid the level of gel rising above 3 cm over the last blade (the closest to the upper edge of the container).
- ◆ Place the unit upon a flat surface, correctly levelled, far from any sources of vibration. Avoid the possibility of it tipping over.
- ◆ Make sure the tap is turned off before starting up the unit.
- ◆ Do not pour the gel at a temperature (°C1), because it may scald you. Wait until the pouring temperatures reaches (°C2).
- ◆ During operation avoid touching the lid. It may have a high temperature. In the gel pouring manoeuvre, wear gloves for protection.
- ◆ Be extremely careful during the slow start ups (the word “LEnt” appears on display (1)). Eventual spattering or spillage of gel may occur.
- ◆ Never leave the gel in liquid state during a delayed start up. It may solidify, preventing a normal start up of the unit. Always insert the

gel in pieces when programming a delayed start up.

### REMINDER

In short, we give you a simple and easy to assimilate manner to carry out the most elementary manoeuvres which can be done on your GEL BAT-6:

#### To switch on the GEL BAT-6:

- ◆ Turn on the main switch (13).
- ◆ Briefly press the key (8) °C1. Light (11) on the control graph will blink intermittently.
- ◆ Press key (8) °(C1) for a few seconds. Light (11) on the control graph will switch on permanently.

#### To switch off the unit (after executing the programme):

- ◆ Press key (9) (°C2) for a few seconds. The lights on the control graph will blink intermittently.
- ◆ Turn off the main switch (13).

#### To switch on the GEL BAT-6 in delayed startup:

- ◆ Turn on the main switch (13).
- ◆ The first possibility is to press key (7) until word “hPro” and alternatively the expected time to reach (°C1) appear. If it is the time you wish to programme for the unit to reach (°C1), press the clock key (7) for several seconds.

#### To turn off the unit (after executing a programme):

- ◆ Press key (9) (°C2) for several seconds. The lights on the control graph will blink intermittently.
- ◆ Turn off the main switch (13).

#### In order to switch the unit off at any time during the programme:

- ◆ Press the key which corresponds to the step being taken for a few seconds, until all the lights are blinking.
- ◆ Disconnect the GEL BAT, acting upon the main switch (13).

### DISPLAY MESSAGES

The GEL BAT-6 is provided with a series of words which appear on the display, in order to facilitate the handling of the unit. These messages are as follows:

- HorA** Indicates the nearest time when the °C1 temperature can be reached.
- HPro** Indicates the time programmed by the user (memorized).
- Fred** Indicates a breakdown in the electric power supply (higher than 10).

- PUEr** Indicates the lid is open.
- Lent** Indicates the unit is working on a slow startup.
- PILA** Indicates the control lithium battery has run out of power.
- rot1** Indicates a fault in the gel temperature sensor.
- rot2** Indicates a fault in the resistance temperature sensor.

#### TECHNICAL CHARACTERISTICS

Height:	500 mm
Width:	325 mm
Depth:	250 mm
Weight:	15 kg
Voltage:	230 V, 50/60 Hz
Power:	700 W
Container capacity:	6 L

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