Ref. 080542

ES.....	MÁQUINA TERMOCONFORMADORA	7
EN	VACUUM FORMING MACHINE.....	9
FR.....	MACHINE A THERMOFORMER	11
PI	MÁQUINA DE TERMO ENFORMAR.....	13
IT	MACCHINA TERMOCONFORMATRICE	15

MESTRA®


TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tel. (+34) 944530388 - Fax (+34) 944711725

mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

1. Resistencia de caldeo / **Heating Resistor** / *Résistance de Chauffe* / Resistência de aquecimento / *Resistenza di riscaldamento*
2. Portaplacas / **Sheet holder** / *Porte Plaques* / Porta placas / *Portapiastre*
3. Interruptor general / **General Switch** / *Interrupteur general* / Interruptor geral / *Interruttore generale*
4. Columna / **Column** / *Colonne* / Coluna / *Colonna*
5. Pomo / **Handel** / *Pommeau* / Puxador / *Pomo*
6. Anillo de cierre / **Ring fastener** / *Anneau de fermeture* / Clip de fecho / *Clip di chiusura*
7. Boca de aspiración / **Suction inlet** / *Bouche d'aspiration* / Boca de aspiração / *Bocca d'aspirazione*



Fig. 1

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Clarifying leaflet to assembly the machina handles

The machine handles are single packed and should be mounting in the vacuum machine before being placed as indicated at the illustrative drawing.

There are two similar smaller handles with M6 threads and two larger handles with M10 threads. **The longest one is at the right side, and first you must insert the small plastic cylinder.**

Mount the handles in the machine at the places indicated at the drawing screwing in the clockwise until feeling the heater at the end of the thread.

Note: The handles cannot be untied and must be pressed only with the hands, do not use any type of setting tool.

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After a thermoforming cycle, let the unit rest for 2 minutes.

Do not make more than 7 cycles; wait 2 hours before continuing.

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The Vacuum Forming Machine (R-080541) is a practical machine with many applications in Dental Laboratory. For optimum performance and a long lasting life, please read the following instructions.

INSTALLATION

1. Carefully unpack the Vacuum Forming Machine. Make sure it was not damaged during transport. In the packing case you should find the following:
 - Vacuum Forming Machine.
 - Two-position model holder (flan – basket).
 - A jar of grit.
2. Choose a firm, rigid and well-levelled surface, far from heat sources, vibrations.
3. Connect the machine to 230 V, 50/60 Hz electric socket with ground wire. The machine is ready to work.

OPERATION

- Select the most adequate model holder position (plane or basket) to the work to be done. We recommend using plane position (A) with sheets of less than 1 mm. In the case of thicker sheets, it is better to use the basket position (B). This basket is also useful when working with models that have retentions so that they avoid keeping a plane seat.
- When using the basket position (B), place the model inside and fill the intermediate space with grit (a jar supplied together with the machine). With this drawing depth is reduced and there is a better finishing.
- Once the suitable model holder position is selected, place it on the suction inlet (6) of the machine.
- Turn and remove the fastener ring (5) on the central column of the machine. You will see that the sheet holder can hold both square (120-131 mm side) and round (134 mm Ø) plastic sheets.
- Place the plastic sheet on the sheet holder and close the fastener again (4). Make sure the sheet is firmly held, there is no looseness and it is well centred.
- Slide the sheet holder down the central column (2) of the machine till it reaches down the closest point to the heating resistor (1). In order to slide it, use the handles (3) located on both sides of the sheet holder.
- Operate the main switch (7). Heating resistor will start and begin heating the plastic sheet softly. The time needed to heat depends on the material and thickness that is going to be vacuum formed. This time does not usually exceed 5 min. A sign that indicates that the sheet is ready is that it should expand between 12 and 15 mm downwards on its central part. In any case, each professional has to make the necessary tests to experimentally determine plasticity conditions of the sheet that best suit work needs.
- Once plastic sheet is under optimal conditions, firmly hold both handles (3) at each side of the sheet holder and slide this by the central column (2) till it reaches and gets to the end of the lower part (model holder).
- Suction motor will start running automatically. You will see that the plastic sheet will adapt to the model irregularities due to vacuum produced. Firmly hold the sheet holder for 10 to 20 seconds.

- **Disconnect machine main switch (7)** and remove model and vacuum formed plastic sheet. The work is ready.

IMPORTANT

Machine heating resistor can be on for a maximum time of approx. 12 minutes. After that, the machine will disconnect automatically. If you want to vacuum form again, you should switch off and then on the main switch (1).

**After a thermoforming cycle, let the unit rest for 2 minutes.
Do not make more than 7 cycles; wait 2 hours before continuing.**

MAINTENANCE AND CLEANING

- Avoid grit coming into the machine. It may cause short circuits or breakages in the suction motor. Handle it carefully.
- The machine is protected with a fuse located on the base of the socket. When needed change the fuse (7 A).
- Cleaning of external body of the machine should be made with a soapy water dampened cloth. Disconnect the machine from the mains to clean it.

PRECAUTIONS

- ◆ Before connecting the machine, make sure it is a 230 V, 50/60 Hz electric point with earth.
- ◆ Do not submerge the machine in any kind of liquid when cleaning it.
- ◆ Do not use the Vacuum Forming Machine in damp places or near flammable liquids (alcohol or solvents).
- ◆ **Heating resistor will start as the main switch (1) is turned on. So, THE MAIN SWITCH HAS TO BE DISCONNECTED WHEN THE MACHINE IS NOT BEING USED.**
- ◆ Do not let children or not qualified people handle the machine.
- ◆ Plastic sheets can reach high temperatures. Do not use fingers or any other part of the body to check their temperature.
- ◆ Avoid grit coming into the machine.

TECHNICAL CHARACTERISTICS

Height:	310 mm
Length:	275 mm
Width:	260 mm
Weight:	4.6 kg
Voltage:	230 V, 50/60 Hz
Heating power:	450 W
Motor power:	1400 W

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