
R-080024

E.....	ELECTROVIBRADOR “TEIDE CIMA”	7
GB	ELECTRO-VIBRATOR “TEIDE CIMA”	9
F.....	VIBRATEUR “TEIDE CIMA”	11
P.....	ELECTROVIBRADOR “TEIDE CIMA”	13

MESTRA®

TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tel. (+34)944530388 - Fax (+34)944711725

E-mail: mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

Rev. 26/05/14

- | | |
|---|--|
| A | Bola de vaciado / Pouring ball / <i>Support-boule pour coulée</i> / Bola de vaziado |
| B | Plataforma de vibración / Vibration platform / <i>Plateau de vibration</i> / Plataforma de vibração |
| C | Regulador de potencia / Power regulador / <i>Potentiomètre</i> / Regulador potencia |
| D | Selector de modo de vibración / Vibration mode selector / <i>Selecteur du mode de vibration</i> / Painel de seleção modo vibração |
| E | Pata ventosa / Suction cap / <i>Pied-ventouse</i> / Pata ventosa |
| F | Interruptor general / Main switch / <i>Emplacement du support-boule</i> / Interruptor geral |
| G | Portafusibles / Fuse holder / <i>Porte-fusibles</i> / Porta fusíveis |

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The R-080024, is a practical electro-vibrator unit specially designed for compacting materials for the jewellery and dental mechanics sectors.

NOTE: Carefully read this manual. It contains important safety information in the installation, use and maintenance. If the equipment is used in a manner not specified by the manufacturer, the security provided by the equipment may be compromised.

MEANING OF THE SYMBOLS ON THE UNIT

Danger. Risk of injury.

Read completely these instructions before operating the unit.

Electric shock risk.

Risk of entrapment.

INSTALLATION

Environmental conditions:

- Use indoors only.
- Environmental temperature: 5-40 °C [41-104 °F]
- Max. relative humidity: 0-95 %
- Degrees of protection: IP20

1. Place the unit on a smooth, firm, perfectly horizontal surface. Try not to install near to other devices, vibrations could be transmitted. If you do not have a surface like that, you can place the unit on a cloth or foamy material.
2. Connect the plug to a 230 V, 50/60 Hz earthed power outlet.
3. Place the device so that it is easy to use the disconnecting device (main switch).
4. Unscrew from the lateral housing and screw it in the upper hole in the vibration platform (B).

PLATFORM AND POURING BALL USE


- Place the cylinder or items to be vibrated on the unit's vibration platform (B). In the case of impression bowls, the mixture can also be compacted by resting the bowl firmly on the edge of the unit and holding down the bowl handle with your hand.

- When the mixture is not fluid, the pouring ball can be useful. Lay the bowl on the ball to transmit vibration to the mixture.

OPERATION

- Press the main switch placed in the back side.
- Pressing upper button, set the vibration time (30'', 1', 1'30'', 2'' or indefinitely).
- Pressing médium button, set the vibration power. We recommend start by the lower power. You can vary the power while unit is running.
- Pressing the lower button, set the vibration mode. There are 4 vibration modes, plus a fifth one that automatically alternates from one to another. You can also vary the vibration mode while the unit is running.
- Press the **START/STOP** Ⓞ button to start the vibration.
- Después de cierto tiempo, la mezcla quedará correctamente compactada.
- After a short time, once the mixture is properly compacted, switch the unit off using the **START/STOP** Ⓞ button.
- When not using the vibrator for a while, disconnect it from the network acting on the rear main switch.
- After putting the main switch off, last used parameters will remain in memory.

MAINTENANCE AND CLEANING

To ensure a long, trouble-free service life and performance of the unit, we recommend as follows:

- After using the unit, always clean away any splashes of mixture before they dry.
- As a precaution, always disconnect the unit from the mains before cleaning it.
- Regularly check the suction caps (E).
- Use only original spare parts supplied by the manufacturer.

PRECAUTIONS

- Before connecting the unit to the power source, make sure the voltage is correct (grounded 230 V, 50/60 Hz).
- The cable may only be replaced by one with the same characteristics as supplied by the manufacturer. Do not use the unit in wet places or close to heat sources.
- Do not place easily broken containers or other objects on the vibration platform.
- Do not allow children or unskilled personnel to handle the unit.
- Do not pull, twist or bend the mains cable excessively.
- Before cleaning the unit, ensure that it is unplugged.
- Never immerse the unit in liquid when cleaning it.
- Do not place hot containers or objects on the vibration platform, as this may damage the rubber.
- Do not remove the vibration platform cover.

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