## R-080060

E.....	LAVAPOLIMERIZADORA .....	6
GB .....	WASHER-POLYMERISER .....	11
F.....	LAVEUSE-POLYMERISEUSE.....	15
P.....	LAVA-POLIMERIZADORA.....	20
I.....	LAVA-POLIMERIZZATRICE .....	25

**MESTRA®**



**TALLERES MESTRAITUA S.L.**

Txori-Erri Etorbidea, 60


Tel. +34 944530388 - Fax +34 944711725

mestra@mestra.es - www.mestra.es


**48150 SONDIKA - BILBAO - ESPAÑA**

A	Tapa / <b>Cover</b> / <i>Couvercle</i> / Tampa / <b>Coperchio</b>
B	Cuadro sinóptico / <b>Synoptic display</b> / <i>Tableau Synoptique</i> / Quadro sinóptico / <b>Quadro sinottico</b>
C	Teclado / <b>Keyboard</b> / <i>Clavier de commandes</i> / Teclado / <b>Tastiera</b>
D	Interruptor general / <b>Main switch</b> / <i>Interrupteur Marche Arrêt</i> / Interruptor geral / <b>Interruttore generale</b>
E	Rejillas de refrigeración / <b>Cooling grille</b> / <i>Grilles de refroidissement</i> / Grelhas de refrigeração / <b>Griglie di raffreddamento</b>
F	Disyuntor / <b>Circuit breaker</b> / <i>Disjoncteur</i> / Disjuntor / <b>Disgiuntore</b>
G	Contacto / <b>Contacto</b> / <i>Contacteur</i> / Contacto / <b>Contattore</b>

1. Display de temperatura / **Temperature display** / *Cadran digital de température* / Display de temperatura / **Display di temperatura**
2. Display de tiempo / **Time display** / *Cadran du temps* / Display de tempo / **Display di tempo**
3. Piloto de control resistencia / **Element control lamp** / *Témoin de la résistance* / Lâmpada controlo resistencia / **Spia di controllo resistenza**
4. Piloto de ducha activa / **Shower on lamp** / *Témoin de la douche* / Lâmpada chuveiro activo / **Spia doccia attiva**
5. Piloto de tiempo cumplido / **Time up lamp** / *Témoin de fin de cycle* / Lâmpada tempo realizado / **Spia tempo concluso**
6. Pilotos selección programa / **Program selection lamps** / *Touches de programme* / Lâmpadas selecção programa / **Spie selezione programma**
7. Piloto indicador marcha / **On lamp** / *Témoin indicateur de Marche* / Lâmpada indicadora funcion / **Spia marcia**
8. Tecla de marcha/paro / **Start-stop key** / *Touche Marche-Arrêt* / Tecla de funcion-para / **Tasto di avvio-arresto**
9. Modo displays / **Display mode** / *Modo displays* / Modo displays / **Modalità display**
10. Tecla de ducha / **Shower key** / *Touche de la douche* / Tecla chuveiro / **Tasto doccia**
11. Tecla decrementar / **Down key** / *Touche de programmation (-)* / Tecla decrementar / **Tasto diminuire**
12. Tecla de función / **Function key** / *Touche de changement de fonction* / Tecla de função / **Tasto di funzione**
13. Tecla de incrementar / **Up key** / *Touche de programmation (+)* / Tecla incrementar / **Tasto incrementare**
14. Pilotos indicadores de paso / **Step indicator lamps** / *Témoins indicateurs des paliers (temps)* / Lâmpadas indicadores de passo / **Spie di passo**
15. Pilotos temperatura / **Temperature lamps** / *Témoins indicateurs de la température* / Lâmpadas de temperatura / **Spia di temperatura**
16. Piloto programación ducha / **Shower programming lamp** / *Témoin de fonctionnement de la douche* / Lâmpada programação chuveiro / **Spia programmazione doccia**

Nº	Ref.	DESCRIPCIÓN	DESCRIPTION
1-A	080060-01A	Carcasa tapa superior	Housing upper lid
1-B	080060-01B	Carcasa base	Housing base
1-C	080060-01C	Carcasa lateral	Housing side
1-D	080060-01D	Bandeja estrecha	Narrow tray
1-E	080060-01E	Bandeja inferior	Lower tray
1-F	080060-01F	Bandeja ancha	Wide tray
1-G	080060-01G	Soporte de carcasa	Tray support
2	080060-02	Bisagras juego	Hinges kit
	080060-02A	Bisagras derecha	Right hinge
	080060-02B	Bisagras izquierda	Left hinge
3	080060-03	Árbol completo	Complete tree
4	080060-04	Boquilla de árbol	Tree nozzle
5	080060-05	Acoplamiento rápido	Fast coupling
6	080060-06	Arandela de cobre	Copper washer
7	080060-07	Tuerca acoplamiento manguera	Hose connecting nut
8	080060-08	Control	Control panel
9	080400-10	Interruptor 230V	Switch 230V
11	080060-11	Placa aislante	Ceramics
12	080060-12	Manguera con abrazaderas (larga)	Hose with clamp (long)
14	080060-14	Manguera con abrazaderas (corta)	Hose with clamp (short)
18	080077-10	Tornillos triac	Triac nut
19	080060-19	Ventilador	Fan
20	080490-37	Tuerca del ventilador	Fan nut
21	080075-08	Tornillo ventilador	Fan screw
22	080077-10	Tornillo de tapa de registro	Inspection cover screw
23	080400-12	Asa	Grip
26	030420-22	Tuerca de termopar	Termocouple nut
27	030420-14	Arandela de termopar	Termocouple Washer
28	080400-23	Termopar	Termocouple
29	080060-29	Resistencia	Resistor
30	080400-27	Tuerca de resistencia	Resistor nut
31	080400-26	Arandela de cobre de resistencia	Resistor copper washer
32	080400-25	Junta tórica de resistencia	Resistor O-ring
33	080400-19	Grifo de desagüe	Drain cock
34	000212	Cable 230V	Cable 230V
	000261	Cable 230V (clavija americana)	Cable 230V (American plug)
35	080070-12	Pasa cables	Cable bushing
37	080060-37	Condensador	Capacitor
38	080060-24	Triac	Triac
39	080060-41	Disyuntor	Circuit breaker
40	080060-36	Motobomba (50 Hz)	Motor pump (50 Hz)
	080060-61	Motobomba (60 Hz)	Motor pump (60 Hz)
41	080400-29	Contactador	Thermal relay

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The R-080060 is a practical washer and polymeriser with an easy-to-use control panel on which sequences of polymerising and shower times and temperatures can be programmed as required by the user. The polymeriser works in the following 6 steps.

1. Waiting time: a time in which the water in the tank is at ambient temperature.
2. 1st ramp time: a time when the temperature of the water in the tank rises steadily to a final temperature set by the user.
3. 1st plateau time: a time when the temperature of the water in the tank remains steady at the level reached at the ends of the previous step. At the end of this time a warning buzzer is set off.
4. 2nd ramp time: a time when the temperature of the water in the tank rises steadily again, to a final level set by the user.
5. 2nd plateau time: a time when the temperature of the water in the tank remains steady at the level reached at the end of the previous step. After these first 5 steps the unit disconnects automatically, a warning buzzer sounds and lamp 4 lights up. The time of each step can be set between 0 and 99 hours 59 minutes. The end temperatures at each plateau can be regulated by the user between 0° and boiling.
6. Shower time: a time when the shower comes on to clean the muffles on the shower tray. This function is independent of the steps programmed previously, and must be activated manually by pressing the shower key (10).

To enable the user to see the actual program status at all times, the R-080060 polymeriser has a synoptic display. LEDs (14) show which step is ongoing and LEDs (15) show the temperatures reached.

The unit also has two digital displays: the top one (1) shows the present temperature of the water in the tank, and the bottom one (2) shows the time remaining to the end of treatment.

If W/T (9) is pressed, display (1) shows the temperature setting expected by the microprocessor, and display (2) the time remaining to the end of the step ongoing.

## INSTALLATION

- Carefully unpack the R-080060 washer & polymeriser. Inside the tank you will find two large muffle trays and two small ones.

- Set the unit up on a flat, horizontal, firm surface away from sources of heat or vibration. Ensure that there is a gap of at least 15 cm between the ventilation grilles and the nearest wall or other object.
- Connect the unit to a 230 V, 50/60 Hz earthed power socket which is protected by a suitable thermal magnetic cut-out. Consult your installer on this.

## SHORT PROGRAM (SHOWER ONLY)

The R-080060 washer-polymeriser can be set for a short shower program in case you wish merely to clean muffles without polymerising them first. This program has two steps:

1. Shower temperature: a time when the temperature of the water in the tank rises to the level set by the user.
2. Shower time: a time set by the user for shower operation.

To choose between the long (polymerisation + shower) program and the short (shower only) program, press key (6 L) long or (6 S) short: the lamp on the relevant key will light up to indicate which program has been selected.

To program shower time with no prior polymerisation, proceed as follows:

1. Press short program key (6 S). The lamp will light up to indicate that short program mode (shower only) has been selected.
2. Press function key (12) to enter programming mode. The top display will read "U1", thus indicating that this step can be programmed. The LED for U1 (shower temperature) on the synoptic display will also begin to flash.
3. Use the Up/Down keys (11 & 13) to bring up the desired first step temperature on the lower display (2). Temperatures of less than 55 and more than 95 °C may not be set.
4. Press function key (12). The top display will read "A5", thus indicating that this step can be programmed. The LED for step A5 (shower time) on the synoptic display will also come on.
5. Use the Up/Down keys (11 & 13) to bring up the desired shower operating time on the lower display (2).
6. When programming is complete, or at any time during programming, programming mode can be exited by pressing W/T (9), which acts as an

escape key. The short (shower only) program is memorised automatically.

7. After exiting programming mode, hold down Stop/Start (8) for approximately one second to start up the shower program. The program can also be stopped at any time by holding down Stop/Start (8) for around 5 seconds.

For the shower to operate the following conditions must prevail:

- The water temperature in the tank must not be less than 55 °C: otherwise the wax dissolved in the water may solidify and clog the pump shaft.
- The water temperature must not be more than 95 °C, otherwise bubbles caused by boiling will appear, leading to cavitation in the pump. If the temperature exceeds 95 °C the machine disconnects the heating elements and attempts to drop to this level.

When the shower is running lamp (4) is lit. The time shown on display (2) is the time remaining until the shower stops, while the figure shown on display (1) is the shower water temperature.

### **LONG PROGRAM (POLYMERISATION + SHOWER)**

It is extremely easy to program a heating and shower sequence for the water in the tank. Merely proceed as follows:

1. Press the long program key (6 L). The lamp for the long program (polymerisation + shower) will light up.
2. Press function key (12) to enter programming mode. Check that the top display reads "A0", thus indicating that this step can now be programmed. The LED for A0 (waiting time) on the synoptic display should also light up.
3. Use the Up/Down keys (11 & 13) to bring up the desired first step time on the lower display (2).
4. Press function key (12). Check that the top display reads "A1", thus indicating that this step can now be programmed. The LED for A1 (1st ramp time) on the synoptic display should also light up.
5. Use the Up/Down keys (11 & 13) to bring up the desired second step time on the lower display (2).
6. Press the function key (12) again. This time the display will read "U1" to indicate that the end temperature for the first ramp can be programmed.
7. Proceed in this way until the times and temperatures for all steps are programmed in.
8. The sequence in which the steps appear on the display when the unit is in programming mode is as follows:

A0- waiting time

A1- 1st ramp time  
U1- 1st plateau temperature  
A2- 1st plateau time  
A3- 2nd ramp time  
U2- 2nd plateau temperature  
A4- 2nd plateau time  
A5- Shower time

9. When programming is complete, or at any time during programming, programming mode can be exited by pressing W/T (9), which acts as an escape key. The program is memorised automatically.

10. After exiting programming mode, hold down Stop/Start (8) for approximately one second to start up the polymerisation and shower program. The program can also be stopped at any time by holding down Stop/Start (8) for around 5 seconds.

Should you wish the program to end at step A2 and miss out the second ramp and plateau, merely allocate a time of 00.00 for steps A3 and A4 and the same temperature for U2 as for U1.

The temperature may not exceed 110 °C at the plateaux U1 and U2. If you try to program a higher temperature, display (2) will show "A" at step U1. If you accept this, the water in the tank will be held at boiling point throughout plateau 2.

Water may boil at slightly above or below 100 °C, depending on altitude above sea level, water composition and atmospheric pressure. This means that if the boiling point programmed into the unit is 100 °C the polymerisation process may not proceed to the next step because the set temperature may never be reached.

Likewise, the polymeriser may go on to the next step without the water having boiled properly. This problem can be solved by selecting "A" at step U1. The R-080060 polymeriser will then automatically detect when the water boils and hold it at that temperature throughout the step in question.

Another possibility is to detect the exact boiling point in a laboratory. The polymeriser is programmed with the following functional parameters: A0=00.00, A1=00.00, U1=110°C. The temperature of the water in the tank will rise slowly to a level at which it holds steady. This is the laboratory boiling point, and this temperature should be programmed whenever you wish to polymerise at boiling point.

To prevent excessive moisture due to steam given off by the washer & polymeriser, Talleres Mestraitua, S.L recommends that you do not program the unit for more than 96 °C. Remember also that the fan of the unit will remain on until main switch D is switched off, to prevent the control module from overheating.

**OPERATION (LONG PROGRAM)**

The R-080060 washer & polymeriser performs both polymerisation and shower functions. It should be used as follows:

- Place the large muffle trays into the bottom of the tank.
- Fill the tank with water up to 3-4 cm below the shower tray support (the small cross-member used to hold muffle trays for washing). There should always be at least 12 cm of water in the tank.
- Place the muffles in the tank, ensuring that they are completely covered with water. Close the cover.
- Switch on the main switch (D). The lamps on the control panel will light up to indicate that the R-080060 washer & polymeriser is running. Lamp (3) comes on when the heating element is running. When the set temperature is reached the control panel automatically controls the heating element to hold that temperature steady in the tank.
- Press (6 L) Long and then program the washer & polymeriser as indicated above. After programming, hold down Start/Stop (8) for one second and the unit will start up automatically.
- When polymerisation time is completed the R-080060 washer & polymeriser disconnects automatically and buzzes for one minute. The buzzer can be cancelled by pressing W/T (9). Lamp (5) also lights up. Open the cover carefully, as outrushing steam from the tank could cause burns. Always use long protective gauntlets, a mask and tongs for muffles. Talleres Mestraitúa, S.L. recommends MESTRA tongs. Close the cover.
- Ensure that the shower is disconnected before removing the muffles. If you are not completely sure how the unit works, switch off at the main switch (D) before removing polymerised muffles. If the shower comes on unexpectedly it can cause major burns.
- The shower can be turned on at any time during the long program (polymerisation + shower) by pressing key (10), and turned off by pressing (10) again. The shower works only when the temperature requirements (between 55 and 95 °C) are met.
- If the shower is stopped via (10) and re-started in less than one minute (e.g. after checking how muffle cleaning is going) the program will re-start from the point where it left off. If the stoppage is longer than one minute the program restarts from the beginning.
- If the water temperature is higher than 95 °C when the shower key (10) is pressed the unit

will disconnect the heating element until the water temperature drops to 95 °C.

- When starting and stopping the shower ensure that the unit cover is closed.

**SHORT PROGRAM (SHOWER ONLY) OPERATION**

The R-080060 washer & polymeriser has a shower function especially designed to scald muffles. To use it correctly, proceed as follows:

- Fit the muffle trays into the top support inside the tank.
- Place the muffles on the tray, using long gauntlets, a mask and muffle tongs. Place them so that the jet of water from the injector sprays the muffle thoroughly. Close the cover.
- The water used for the shower is water from inside the tank. A closed circuit enables this water to be recycled. This means that the scalding temperature is that shown on the display
- Press (6 S) Short (shower only) to enter programming mode. Program the shower temperature and time on the control panel as indicated above.
- Hold down the Start/Stop button for one second. The short shower program can then be run when required by pressing (10). When the shower is running control display (2) will show the time remaining until the shower stops automatically, and display (1) will show the water temperature in the tank. The shower can also be stopped manually at any time by pressing (10).
- Once the set shower time is completed, open the cover carefully, as outrushing steam from the tank could cause burns. Use long gauntlets and a face mask. Remove muffles with special tongs. Close the cover.

**MAINTENANCE & CLEANING**

- Clean the inside of the tank after every 60 hours' operation approximately, by heating the water to 75 °C, then emptying the tank through the drainage tap at the back. Then unplug the R-080060 washer & polymeriser and clean the inside of the tank with a fine scourer moistened in soapy water, taking care not to damage the heating element or temperature probe. Special products for cleaning polymerisers may also be used. Rinse the inside of the tank with plenty of water and then close the drainage tap.
- For easy access to the tank of the unit, the shower spray tube can be dismantled by turning the quick-lock collar to the left. Ensure

that the unit is re-assembled correctly after cleaning.

- At intervals which will depend on how hard the water is, clean away any limescale which may have been deposited on the shower nozzles or electric heating element, using a scale remover. Anti-limescale products may also be used.
- To clean the outside housing of the unit use a cloth moistened in soapy water. Do not use solvents or other inflammable products.
- To get the best results from the washer & polymeriser additives available on the market to dissolve wax may be added to the water. The manufacturers' recommendations should always be followed in this case.
- For optimum performance, Talleres Mestraitúa, S.L. recommends the use of MESTRA muffles and flanges.
- Clean the unit's control panel regularly with a damp cloth. Always disconnect the machine before doing this.
- After using the washer & polymeriser, remove wax, foam, remains of resin and other particles which may be left floating on the water.
- The R-080060 has a circuit breaker that protects the machine against electric overloads. Once the problem is solved, press the circuit breaker located on the side of the machine (F) to get it into work again.
- Message “bro” means that the thermal probe is broken and must be replaced.
- Message “HHHH” means that the triac is broken and must be replaced.
- The R-080060 has a contactor that protects the machine against accidental overheating as it detects excessive temperature and cuts off machine operation. Once the machine gets back to a normal temperature, the security contactor has to be reset. In order to get to it, do as explained below:
  1. Remove the side panel unscrewing the fixation screws.
  2. Disconnect the electric fan.
  3. Press the contactor (G).
  4. To assemble the machine do the other way round.

Try and check the Wahser-Polymeriser functioning.

**PRECAUTIONS**

- ◆ Before connecting the unit ensure that the power supply is 230 V, 50/60 Hz and earthed. It should also be protected by a thermal magnetic cut-out. Consult your installer on this.
- ◆ Do not allow children or unskilled personnel to handle the unit.

- ◆ While the shower is running, keep the cover of the unit closed.
- ◆ After removing the shower spray tube, check that the connector is properly connected. Never run the shower without the spray tube fitted.
- ◆ Whenever you put muffles into the tank or take them out, use long gauntlets and a muffle tongs.
- ◆ If the shower stops suddenly during normal operation, disconnect the unit from the power supply before manipulating it in any way.
- ◆ Open the cover carefully: outrushing steam from the tank could cause burns.
- ◆ Never place more than 12 muffles on the polymerising tray, or more than 6 on the washing tray.
- ◆ Always keep at least 15 cm of water in the bottom of the tank, as the heating elements may be damaged if the level is lower than this. Do not allow water level to reach less than 3 cm from the muffle washing tray.
- ◆ Always leave a gap of at least 15 cm between the unit's ventilation grill and the nearest wall or other object. Place the unit on a flat, firm, vibration-free surface where it cannot overturn.
- ◆ 95 °C is a recommendable polymerisation temperature: at higher temperatures too much steam may be produced.
- ◆ A typical shower temperature would be 85 °C. Higher temperatures may lead to cavitation in the pump.

**TECHNICAL CHARACTERISTICS**

Height:	450 mm
Width:	635 mm
Depth:	405 mm
Weight:	34 kg
Power supply:	230 V, 50/60 Hz
Consumption:	2200 + 400 W

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