
**R-080262 / R-080263 / R-080264**  
**R-200262 / R-200263 / R-200264**

<b>ES ... CHORREADORAS DE ARENA CONSTELLATION.....</b>	<b>5</b>
<b>EN ... CONSTELLATION SANDBLASTERS .....</b>	<b>8</b>
<b>FR ... SABLEUSES CONSTELLATION .....</b>	<b>10</b>
<b>PO ... JACTO DE AREIA CONSTELLATION .....</b>	<b>13</b>
<b>IT ..... SABBIATRICI CONSTELLATION.....</b>	<b>15</b>
<b>DE ... CONSTELLATION SANDSTRAHLER .....</b>	<b>18</b>

**MESTRA®**

TALLERES MESTRAITUA S.L.  
Txori-Erri Etorbidea, 60  
48150 SONDIKA - BILBAO - ESPAÑA  
Tel. +34 944530388 - Fax +34 944711725  
mestra@mestra.es - www.mestra.es

Nº	Ref.	DESCRIPCIÓN	DESCRIPTION
1A	080262-56	Carcasa	Housing
1B	080262-57	Carcasa de módulo	Module housing
1C	080262-58	Bandeja	Tray
1D	080262-59	Soporte cristal	Glass support
2	080262-02	Cristal	Glass
3	080262-03	Burlete cristal	Rubber strip
4	080235-60	Tira de LED	LED bar
5	080235-62	Cubierta LED	LED shield
6	080235-63	Tapas LED	LED lids
7	080242-28	Manguitos	Wool sleeves
8	080270-30	Bridas	Flanges
9	080262-07	Pomo de llave	Key knob
10	080262-08	Retén pequeño	Small seal
11	080262-09	Retén grande	Large seal
12	080270-33	Tapón vaciado	Drain plug
13	080010-05	Pasa cables	Cable bushing
14	080262-31	Cuerpo lapicero	Gun body
15	080262-32	Pieza intermedia lapicero	Gun intermediate part
16	080240-15	Boquilla 0,75 mm 50 micras	0,75 mm 50-micron nozzle
17	080240-12	Boquilla 1,00 mm 150 micras	1,00 mm 150-micron nozzle
18	080240-13	Boquilla 1,25 mm 250 micras	1,25 mm 250-micron nozzle
19	050221	Soplador	Blower
20	030350-48	Junta tórica	O-ring
21	080262-33	Punta de lapicero rojo	Gun point (red)
	080262-38	Punta de lapicero azul	Gun point (blue)
	080262-39	Punta de lapicero verde	Gun point (green)
	080262-45	Punta de lapicero amarillo	Gun point (yellow)
22	080490-31	Pata	Leg
23	080242-33	Pedal	Pedal
24	080242-25	Codo de plástico	Plastic elbow
25	080240-31	Filtro electroválvula	Electro valve filter
26	080242-23	Electroválvula 230V	Electro valve 230V
	080242-23B	Electroválvula 115V	Electro valve 115V

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
The Constellation sandblaster family is composed of several sandblasting chambers specially designed for precision-blasting small parts used in the dental prosthesis industry. For optimum output and performance of this unit, we advise you to read the following instructions with great care:

### INSTALLING THE UNIT

- Carefully unpack the unit. You will find the following fittings inside the chamber:
  1. An air connection tube, diameter 6 mm.
  2. A 90 ° PVC elbow fitting to connect the outgoing tube to suction if required.
  3. Pedal drive mechanism.
- A rigid flat surface must be chosen for installation of the unit, at a safe distance from any sources of vibrations or heat, and the surrounding atmosphere must not contain excess moisture, since this could cause compacting of the emery particles inside the chamber, and lead to obstructions in the microblasting nozzles. The unit must also be positioned on a level which allows objects inside the chamber to be handled with ease.
- Hook up the air connection tube supplied with the unit to the connection at the rear of the unit. The other end must be fitted to a compressed air facility, with at least 7 bars of pressure and minimum flow rate 100 litres/minute.
- Hook up the electric plug to 230 V, 50/60 Hz, which must be earthed.
- Connect the pedal drive tube to the coupling nozzle on the right-hand side of the unit. Then position the pedal for comfortable foot operation. Depending on operational requirements and/or handling comfort, it is possible to position the pedal for knee operation.
- Unscrew the sandblasting module covers (A), and fill the pans with emery. During this operation, **it is most important not to fill the pans to more than ¾ of their capacity.** Put the covers back on, and ensure they have been screwed on tightly.
- In relation to the calibre of the nozzles, the machine is supplied ex-works to the following criterion of nozzle diameters and colours:
 

Yellow nozzle	Ø 0.75 mm	
Blue nozzle	Ø 1 mm	
Green nozzle	Ø 1 mm +	
Red nozzle	Ø 1.25 mm	

### REGULATING THE PRESSURE

The blowing pressure of the blasting guns can be regulated depending on the nature and characteristics of the operation to be carried out.

Use the pressure-adjustment knob (B) on the top right of the unit – pull the knob upwards, and rotate it to the required pressure. The unit has a pressure gauge (C) for easy selection of operating pressure. When the pressure has been adjusted, lock out the regulator knob by pushing it down to the initial position.

Pressure gauge (C) shows the blasting module gun blow pressure at all times, even when the gun is non-operational.

### CONNECTION TO THE SUCTION UNIT

The microblasting chamber must be permanently connected to forced suction via an outgoing tube on the top left of the unit (G). Suction must be able to produce a slight depression inside the chamber for continuous gentle circulation of the powder generated towards the outgoing tube.

Depending on the requirements and geometry of the external suction unit, it is possible to use the PVC elbow fitting supplied for specific needs.

### USING THE UNIT

The following operations are recommended for proper precision-blasting:

1. Open the lid of the chamber, and put in the part you wish to blast. Parts may also be placed inside the chamber via sleeve (F). Check there is enough emery in the pans.
2. Close the lid of the chamber, and ensure it is tightly shut (the glass should fit snug against the housing).
3. Suitable protection gear must be used: goggles, gauntlets and a mask.
4. Start up the forced suction facility. The light in the chamber may be activated for a better view of the part, if required, using switch (D).
5. Use the main switch (E) to select the pan. The colour coding system used makes it easy to see which is the active pan at all times.
6. Pick up the mobile gun selected inside the chamber. The colour of the gun's nozzle-holder matches the colour on the pan covers (red, blue, yellow, green). When the gun has been selected, press the drive pedal to begin sandblasting.
7. Direct the blast carefully to prevent the stream hitting the sides of the chamber directly, since this could damage the chamber. If you wish to test the blasting pressure or strength, direct the stream at the rubber dolly inside the chamber.
8. After the blasting operation, release the pedal, stopping the sand flow and activating immediate air discharge to relieve pressure in the pans. Then switch off the suction unit, open the lid of the chamber, and take out the part.

## MAINTENANCE AND CLEANING

To ensure a lengthy service life, we recommend the following:

- Always use approved quality emery - **Talleres Mestraitua, S.L.** recommends **MESTRA** emery. The mobile blasting guns with nozzle diameters 1.25 mm and 1 mm can respectively blow emery of maximum diameter 250 and 150 microns. The mobile blasting gun with nozzle diameter 0.75 mm can blow a maximum diameter of 75 microns (see the table).
- It is preferable to use dry-piston air compressor units, since oil in the air could cause the sand to compact and obstruct the gun nozzles.
- Never fill the modules with emery deposited on the floor of the chamber, since particles with large diameters could obstruct the nozzles, and this emery is now contaminated with particles eroded off the parts blasted.
- The silencer filters inside the pans must be replaced every 5,000 or 10,000 operations (around 1 or 2 years). The incoming filter fittings on each pan must be replaced every 500 or 1,000 hours (much depends on the quality and metrage of the sand used). Original MESTRA spare parts must always be used.
- To prevent any abrasion wear on the walls inside the chamber, take care not to direct the stream of emery directly onto the walls.
- Clean out the emery inside the chamber on a periodic basis, using the practical emptying plug inside.
- Any water which has condensed inside the incoming air filter must be bled out regularly. The chamber has a window to view the amount of water in the filter.
- Use a dry cloth or a cloth moistened with soapy water for general cleaning of the unit.
- If the guns become obstructed, remove the nozzle-holder and take out the nozzle. Unblock the nozzle by blowing in the opposite direction to the flow of air, or use a pin or a thin wire to unblock it, and then refit the attachment.
- The blasting module guns on the 2-pan Constellation model are supplied ex-works with nozzles 1 and 1.25 mm in diameter. 0.75 nozzles, however, may also be used with this model (consult your **MESTRA** distributor).
- The maximum diameter of emery particles in relation to nozzle diameter is shown in the table below:

MAXIMUM PARTICLE DIAMETER / NOZZLE	
Nozzle	Diameter
1.25 mm	Up to 250 microns
1 mm	Up to 150 microns
0.75 mm	Up to 75 microns

## PRECAUTIONARY MEASURES

- ◆ Never exceed 5 bars of pressure. Optimum working pressure is between 2 and 4 bars.
- ◆ Operations should always be carried out with the lid of the chamber closed, with suitable protection. Masks and goggles should be worn.
- ◆ Connect the chamber's outgoing tube to a suction facility.
- ◆ The protection sleeves (F) must be replaced immediately following any tearing or damage.
- ◆ Do not allow children or unqualified staff to operate the unit.
- ◆ Use only guaranteed-quality materials as emery (we recommend **MESTRA** emery). Poor quality emery may be toxic when inhaled, and may also cause obstructions in the blow nozzles.
- ◆ Hook up the power supply cable to 230 V, 50/60 Hz, earthed.
- ◆ Consult your fitter immediately if any leakages are observed in the compressed air circuit.
- ◆ Before carrying out any operation, first ensure that the micro blasting module covers are closed tightly.

## TECHNICAL CHARACTERISTICS

Height:	320 mm
Length:	390 - 510 mm
Width:	440 mm
Weight:	12.5 - 14 - 15.5 kg
Working pressure:	2-6 bar
Voltage:	230 V, 50/60 Hz
Power consumption:	4 W

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**GUARANTEE**

All **MESTRA®** units are guaranteed against every defect in manufacturing for a period of one year from date of purchase. Prior to honouring a guarantee, **MESTRA®** will require that the attached coupon be filled in with date of sale, duly stamped by the distributor.

The guarantee does not cover breakdowns caused by incorrect installation or use.

**MESTRA®** will not effect any repair outside its own workshops and will require transport costs.

Repairs carried out on the premises of the customer, be they with or without guarantee, are entrusted to the technical assistance of the distributor who has sold the unit or alternatively to the technicians he has contracted

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*TALLERES MESTRAITUA S.L. is a company manufacturing tools and appliances for the dental prosthesis sector, from the business address shown at the bottom of the page.*

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*I HEREBY STATE: That the “CONSTELLATION” SANDBLASTERS R-080262, R-080263, R-080264, R-200262, R-200263, R-200264, with the following characteristics:*

Alto / Height:	320 mm
Largo / Length:	390 - 510 mm
Ancho / Width:	440 mm
Peso / Weight:	12.5 - 14 - 15.5 kg
Presión trabajo / Working pressure:	2-6 bar
Tensión / Voltage:	230 V, 50/60 Hz
Consumo / Power consumption:	4 W

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*Meet the low voltage requirements set out in Directive 2014/35/UE and in Directive 2014/30/UE concerning electromagnetic compatibility.*

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*In witness whereof, the company issues this declaration of conformity in Sondika, on 09 February 2022.*

Firma / Signed:

Ignacio Mestraitua,  
Director Gerente / *Managing Director*

### MESTRA®

TALLERES MESTRAITUA S.L.  
Txori-Erri Etorbidea, 60  
48150 SONDIKA - BILBAO - ESPAÑA  
Tel. +34 944530388 - Fax +34 944711725  
mestra@mestra.es - www.mestra.es