R-080460 POLY-PLUS

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The R-080460 POLI-PLUS unit is a practical pressurized polymerizer which incorporates a control panel in order to programme a sequence of polymerization temperatures and time, according to what is preset by the user. Therefore, any polymerization procedure is made up of the following 5 steps:

- A period when the water inside the container is at room temperature and atmospheric pressure (initial period).
- A period when the water temperature inside the container increases linearly, until it reaches a final temperature set by the user (period of 1st slope). At the beginning of this phase, the inside of the container is pressurized at a value set by the user, which is maintained throughout the entire programme.
- A period when the water temperature inside the container remains constant, taking the final value used in the previous step (temperature reaches 1st plateau).
- A period when the water temperature inside the container varies linearly, until it reaches the final temperature set by the user (period of 2nd slope).
- A period when the water temperature inside the container remains constant taking the final value acquired in the previous step (temperature reaches 2nd plateau).

Once this sequence of steps is finished, the unit will disconnect automatically and let off a warning beep for one minute, and the water inside the container is discharged into a drainage bottle, also re-establishing the atmospheric pressure. The duration of each step may be set from 0 to 99 hours and 59 minutes. The final temperatures of each of the two plateaus may be regulated between 0 °C and 135 °C.

In order to show the user the true status and evolution of the programme, the control panel on the R-080460 pressurized polymerizer comes with a synoptic graph (11). The LED lights indicate at all times the step being taken and the temperatures reached. The unit is also equipped with two digital Displays. The upper one (1) indicating the instantaneous water temperature inside the container. The lower Display (2) indicating the time remaining to the end of the programme.

Both displays have a double function, which is activated by pressing (W/T) (7). By pressing this key, the upper display (1) shows the temperature the user has programmed as the final value of the step (assigned temperature). The lower display

shows the time left to finish the step being executed.

PROGRAMMING

Programming a water heating law inside the container is very simple. In order to do this, we recommend you carry out the following steps:

1. Activate the main switch (12).
2. Press function key (10). This operation enables you to enter into the programming mode. You will notice the upper display reads "A0" indicating that the programming of this step has been activated. The LED light on the synoptic graph will turn on which corresponds to "A0" (starting time).
3. Press the increase/decrease keys (8), (9), until the lower display (2) shows the time value you wish to use for the first step. (Standby time).
4. Press function key (10) again. You will notice the upper display reads "A1", indicating that the programming for this step has been activated. The LED light on the synoptic graph will turn on which corresponds to "A1" (period for the 1st slope).
5. Press the increase/decrease keys (8), (9) until the lower display (2) shows the time value you wish to use for the second step.
6. Press function key (10) again. This time you will read U1 on the display, asking you to programme the final temperature of the first slope.
7. Proceed the same above to complete the programming of the duration and temperature of the different steps in the programme.
8. The ordered sequence of the steps which will appear on the display when the unit is in programme mode, is as follows:
 - (A0)- Starting time
 - (A1)- 1st slope period
 - (U1)- Temperature of 1st plateau
 - (A2)- Time of 1st plateau
 - (A3)- 2nd slope period
 - (U2)- Temperature of 2nd plateau
 - (A4)- Time of 2nd plateau
9. You have three ways of leaving the programming mode: Go through any of the function menus; press key (7) (W/T) at any time, which acts as an escape key; or remain in visualization position without pressing any of the keys for a period of more than 30 seconds, the programme will be automatically stored in memory. We should point out that when key (W/T) is used as an escape key, the control

jumps into start situation, saving the selected parameters in memory.

10. If for some reason you wish the programme to finish in step "A2" (1st plateau), without running the rest of the programme, just go to step "A3" and then press key (9) decrease until the display (2) shows the word 'Fin'.
11. After leaving the programming mode, in order for the unit to run the programme you need to open and close the lid, this is required in order to fill the container with water and to insert the elements you wish to polymerize. In order to facilitate handling and programming the unit, the control display (2) gives the user assistance messages. The lid opening message is represented by the word "AbrA". Close the lid with the message on "Cerr". Once the opening and closing manoeuvre has been completed, the word "StAr" will appear. As of this moment, the programme will start-up automatically when the (P/M) (5) key is pressed.

INSTALLING THE UNIT

- Carefully unpack the R-080460 pressurized polymerizer. The accessories will include a drainage bottle (which does not appear in the figure) with its connection pipe, and a connection pipe for the pressurized air network. You will also find two trays inside the container. One of these has a handle and is used to easily take out the polymerized elements. The other has no handle, and is used as a protection for the heating resistance. Place the tray without a handle on the bottom of the container just above the resistance.
- Place the unit upon a horizontal, flat and rigid base, far from heat or vibrating sources. Make sure that a connection to the pressurized air network is close by. The sides of the unit are provided with ventilation grooves. In order to guarantee appropriate ventilation, make sure there is at least 15 cm between the grooves and the nearest wall.
- Connect the device to a 220 V, 50/60 Hz electricity outlet, provided with an earth connection.
- Connect one of the pneumatic tube connection ends to the manoregulator connector (F). The other end should be connected to the pressurized air network.
- When locating the drainage bottle be sure to choose a place at least 20 cm below the lowest part of the unit. The water discharged from the container is at boiling temperature, therefore if the bottle were to fall it may cause an accident. Therefore, make sure the drainage bottle is

placed in a location which has optimum stability and safety. Once the bottle has been placed, connect the drainage tube to the connector located on the back of the unit, and insert the other end through the upper bale of the bottle.

HOW TO PROGRAMME

Let's suppose we want to programme the polymerizer in order to carry out the following steps:

1. Start time- We want the polymerizer to start operating without a waiting time.
2. Period of the 1st slope- We want the water inside the container to reach a temperature of 130 in 1 hour.
3. Period of the 1st plateau- We want the water inside the container to continue at 130 °C for 5 minutes.
4. Period of the 2nd slope- Subsequently, the water temperature inside the container will descend to 100 °C in 1 hour.
5. Period of the 2nd plateau- When the 100 °C are reached, the water will be kept at this temperature for 30 minutes. At this time, the programme will be completed and the polymerizer will let off a signal and discharge the water from the container. Light (4) will also turn on.

The programming steps are as follows:

1. Press key (10) to enter into programme mode. "A0" will start blinking.
2. Press the decrease key (8), to place the start time at zero.
3. Press key (10), to enter into the next step. The "A1" lamp will start to blink.
4. Press keys (8) and (9) increase/decrease, to select 1 hour as the 1st slope period.
5. Press key (10) to go onto the next step. "U1" will start to blink.
6. Press keys (8) and (9) increase/decrease to select the value of 130 °C, as the temperature value during the 1st plateau.
7. Press key (10) to go onto the next step. The "A2" light will turn on.
8. Press keys (8) and (9) increase/decrease, to select 5 minutes as the period for the 1st plateau.
9. Press key (10) to go onto the next step. Light "A3" will turn on.
10. Press keys (8) and (9) increase/decrease, to select 1 hour as the period for the 2nd slope.
11. Press key (10) to go onto the next step. The "U2" light will turn.
12. Press keys (8) and (9) increase/decrease to select the value of 100 °C, as the temperature value during the 2nd plateau.

13. Press key (10) to go onto the next step. Light "A4" will turn on.
14. Press keys (8) and (9) increase/decrease, to select 30 minutes as the period for the 2nd plateau.
15. Press the function key (10) or (7) (W/T) to record the programme you've just carried out in the memory.

From now on when you open and close the lid, the display (2) will show the word "STAr". Then press key (5) (P/M), and the R-080460 pressurized polymerizer will start up automatically running the programme.

OPERATION

In order for the polymerization programme to commence it is always necessary to open and close the lid, enabling you to fill the container with water and insert the elements you wish to polymerize. The steps taken in order to execute a polymerization cycle is as follows:

1. Activate the main switch (12). You will notice how part of the unit's control panel elements are turned on indicating that the R-080460 polymerizer is in operating mode. In this initial situation, you will notice that the display (2) shows the word "AbrA", indicating that the unit is ready to commence the polymerization cycle.
2. Open the lid of the unit first of all by turning the safety handle (B) and then moving the sealing slide (A). You will notice that the word on the display has changed to "Cerr", indicating that the unit is waiting for the lid to be closed before commencing the programme.
3. Place the elements you wish to polymerize in the tray with the handle located inside the container.
4. Fill the container with water the final level of the water should come up to 3 cm below the top edge of the container. Make sure the water level in no case reaches the pressurized air inlet connector level which is located inside the container (upper part). We do not recommend water levels below 5 cm of the container tray, as they may leave the heating resistances uncovered.
5. Check to see if the drainage bottle is empty. If the bottle is full, water may overflow at boiling temperature during the discharge phase.
6. Programme the pressurized polymerizer by following the instructions indicated in the programming section.
7. Close the lid and make sure it is fastened by using the sealing slide and safety handle. Sometimes, the rubber seal placed on the upper edge of the container, which guarantees the water tightness of the lid, may be displaced. Pay

special attention to this, because if the seal is incorrectly placed this may cause steam to leak out. If you detect the seal is incorrectly placed, proceed to fit it correctly in the container groove.

8. Once the lid is closed, the display message will read "StAr", indicating the unit is ready to commence the polymerization cycle.
9. When you press the (P/M) (5) key the POLY-PLUS unit will automatically start-up executing the programme which has been defined. Light (6) will stay turned on, which indicates the programme is being executed without any problems. Light (3) will also remain turned on when the heating resistance is implemented. The programme's graph lights (11) will gradually turn on indicating the phase of the programme you are in.

During the normal running of the programme it is possible to interrupt the unit by pressing the (P/M) (5) key for 5 seconds.

The display will show the message "dESC" meaning you can open the cock to unload the water from the tank.

You may also, at anytime enter into the programming mode and alter the parameters of the programme being executed. In order to do this press the programming key (10) and proceed according to the programming section. During a power supply failure, the POLY-PLUS goes on to off status, the display reading "FrEd", when the power supply is returned.

The pressure gauge (13) indicates at all times the pressure inside the container. This value may be changed at will between 0 and 6.5 bars (maximum pressure) by using the manoregulator (F). In order to do this, pull the manoregulator knob upwards. Then by turning the knob clockwise you will increase the pressure. The pressure is reduced by turning the knob anti-clockwise. Once you have completed the pressure regulation procedure, lock the manoregulator knob into position by pressing downwards. The unit is delivered with a regulated working pressure of 5 bars. From the experimental results achieved, TALLERES MESTRAITUA, S.L. recommends you always maximum working pressures of between 5 to 5.5 bars.

10. Once the polymerization programme has been completed, the R-080460 will warn you by letting off a beep. Light (4) on the control panel will also turn on. Press any key to stop the beep and close the air intake. The machine is ready for download. The display will show the message "dESC". The discharge time takes approximately 3.5 minutes. It is of utmost importance not to open the lid of the unit until

the discharging procedure has been entirely completed. As the water steam released may cause burns.

11. Once the discharge procedure has been completed (approximately 2 minutes), the display (2) will intermittently read "Fin" and "AbrA", indicating that you may open the lid.
12. In special circumstances, for example, when the container is full of water and at a very high temperature, the 2 minute discharge time may be insufficient to completely empty the container. Therefore, before opening the lid of the POLY-PLUS unit always make sure the machine has stopped discharging water or steam, and that the pressure on the pressure gauge (13) is 0.
13. When the programme is completed without any user interruptions, the display (2) will continue to intermittently read "Fin" and "AbrA" even after the lid has been opened. In order to return to regular control status, just press key (7) (W/T).
14. Take out the elements you have just polymerized. the tray supplied with the unit has a practical handle which is of great use when taking it out of the container. When carrying out this operation, wear gloves and a mask, because the inside of the container is at a high temperature.
15. Finally, close the lid of the POLY-PLUS unit and disconnect it by pressing the main switch (12). After a period of time, when the drainage bottle temperature has decreased, you may empty it.

MAINTENANCE AND CLEANING

In order to guarantee the optimum reliability and a long working life of the unit, we recommend the following:

- In order to avoid wax obstructions, regularly clean the interior of the container (approx. after every 5 uses). In order to do this heat the water to 120 °C. Then, empty the container by means of the drainage bottle.
- Then unplug the POLY-PLUS and clean the inside of the container with a fine swab dipped in soapy water, taking care not to damage the

resistance nor the temperature probe. Specially prepared products may also be used to clean polymerizers. Rinse the inside of the container with plenty of water.

- Periodically and depending upon the hardness of the water, clean the lime remains which may be deposited on the electric resistance. A decalcifying product may be used for this work.
- Do not use solvents or any other inflammable product.
- We recommend lubricating the shaft of the lid seal slide (A), and the rotation shaft of the safety handle (B) regularly. Use a thick oil.
- After using the polymerizer get rid of the wax, foam, resin remains or other particles which may be found on the container walls
- If too much water is placed in the container, it's possible that part of this water may recirculate towards the manoregulator (F), being discharged through some holes located under the knob, these act as spillways. This situation - although not at all recommended- has been foreseen in the unit design phase, it is therefore not at all dangerous.
- Regularly purge the water accumulated in the manoregulator filter. In order to do this, slowly loosen the screw located on the lower part of the filter until water commences to flow through the screw. When the water has been purged, proceed to firmly re-tighten the purging screw.
- If the temperature probe on the POLY-PLUS breaks, display (2) will show the word "rot". Immediately get in touch with a MESTRA authorized dealer.
- It is convenient to regularly replace the rubber seal on the crown of the unit's container, in order to achieve a good lid tightness. In order to do this use MESTRA original spare parts, given that these are the only ones able to guarantee an optimum performance of the unit.

PRECAUTIONS

- ◆ Before connecting the unit make sure the connection plug runs on 220 V, 50/60 Hz including an earth connection.
- ◆ Place the drainage bottle in a stable location which prevents it from overturning by mistake. Before commencing a new cycle, make sure the drainage bottle is empty. Do not connect the drainage pipe to another container different from the bottle supplied with the unit.
- ◆ Before starting up the unit, make sure the drainage bottle cap is tightly screwed on, and that the connection between the pipe and the cap is tight.
- ◆ Do not allow children or unqualified personnel to handle the unit nor the drainage bottle.

- ◆ Protect yourself appropriately when placing or taking out muffles from inside the container. Always use gloves and a protection mask.
- ◆ Open the lid of the unit carefully. Re sudden release of steam from the accumulated water in the container may bum.
- ◆ Always close the lid of the unit with the safety components (slide (A) and safety handle (B)). Never open the lid when a programme is running.
- ◆ Do not exceed 6.5 bars of pressure. The recommended maximum operating pressures are approximately 5 or 5.5 bars.
- ◆ When a programme has been completed, wait until the display alternatively shows the word "Fin" "AbrA", in order to open the lid. Also take care and check that the unit does not continue to discharge water or steam from the bottle, and that the pressure indicated on the pressure gauge (13) is 0.
- ◆ Do not open the lid during the discharge phase ("dESC") or while the unit is running a programme.
- ◆ Make sure the water level of the container is never less that 5 cm from the bottom of the tray, or over 3 cm from the upper edge of the container.

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