G E L B A T - 6

G E L B A T - 6

G E L B A T - 6

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The MESTRA gel mixer is designed for a long, efficient, maintenance-free working lifetime. Its efficiency in cutting, heating, beating and pouring makes it an essential element for duplication. It is easy to handle, and offers denture makers a wide range of possibilities.

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Easy-to-use tap made of non-stick thermoplastic.

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
Vane inclination and distribution designed for even gel mixing.

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Its high-torque gear-motor enables the unit to start up with solidified gel. A heat relay safety system protects the motor against overheating.

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A safety switch prevents the motor and heating element from working while the lid is open.

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
Forced-draught tank cooling on reaching setting C1.

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LATEST GENERATION MICROPROCESSOR

- Temperature control
- Forced draught cooling control
- Last program memory
- Buzzer to indicate completion of cycle C1
- Power off warning
- Tank thermocouple breakage warning
- Door open warning
- Low battery warning
- Protection against voltage variations
- Protection against overvoltage
- Protection against motor blockage
- Protection against element overheating

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ERGONOMICS

Easy-to-read display screen
Graphic with lights showing which process is ongoing.
Smooth-action touch-controls
EASY TO HANDLE

- [C1] C1 temperature selector
- [C2] C2 temperature selector
- [⌚] Timer key
- [▲] Key to increase settings
- [▼] Key to decrease settings

FEATURES

With the microprocessor the user can set:

- gel melting temperature from 65 to 95°C
- gel pouring temperature from 35 to 65°C
- delayed start up times (nocturnal program)

FEATURED SYSTEMS TROUBLESHOOTING CHECK

In normal operation the machine runs a check once a minute to ensure that motor resistance is not too high. This can be detected visually during heating.

The motor cuts out for a tenth of a second every minute, and when this occurs the pilot lamp for the vanes on the microprocessor controls goes out.

MESTRA has incorporated this system into the microprocessor software in case the gel solidifies in the tank for any reason (power cut, etc.). The control detects this and attempts to start up the motor and cut the gel. If it cannot do this, it cuts off power to the motor and only allows the elements to heat on slow until it considers that the resistance met by the vanes is low enough for them to start up without difficulty.

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