
R-080560


ES.....	5	Manual original
EN.....	11	
FR.....	16	

Nº	Ref.	Descripción	Description
1A	080561-33	Tapa trasera	Back lid
1B	080561-34	Tapa lateral izquierda	Left side lid
1C	080561-35	Puerta	Door
1D	080561-36	Chapa protección control	Control protection sheet
1E	080561-37	Escuadra fijación pared	Wall fixing square
1F	080561-38	Carcasa	Housing
1G	080561-39	Tapa lateral derecha	Rigth side lid
1H	080561-40	Tapa delantera	Front lid
2	080561-02	Varilla bisagra	Finge stick
3	080561-03	Chapa protección aislante	Insulation protection sheet
4	080561-04	Junta	Gasket
5	080561-05	Cámara de calentamiento	Heating chamber
6	080561-10	Placas aislantes	Insulations boards
7	080561-07	Resistencias	Resistence
8	080561-08	Sonda	Termal probe
9	080561-09	Inyector	Injector
11	080561-11	Racor inyector	Injector connector
12	080561-12	Tuerca inyector	Injector nut
13	080561-13	Distribuidor	Distribuitor
14	080561-14	Racor distribuidor codo	Elbow distribuitor connector
15	080561-15	Racor distribuidor recto	Straight distribuitor connector
16	080561-16	Regulador de aire	Air regulator
17	080561-17	Racor en codo	Elbow connector
18	080561-18	Racor en T	T connector
19	080561-19	Arandela inyector	Injector washer
20	080561-20	Cilindro	Pneumatic cilinder
21	080561-21	Racor en codo cilindro	Cilinder elbow connector
22	000333	Racor recto hembra	Straigth feme connector
23	080561-22	Racor en codo manómetro	Manometer elbow connector
24	000308	Racor en codo Ø4	Ø4 elbow connector
25	000312	Racor en codo Ø6	Ø6 elbow connector
26	080242-34	Manorreductor	Pressure control
27	080550-06	Microinterruptor	Microswitch
28	080561-01	Husillo con pomo	Spindle with pome
29	080561-24	Porta husillo	Spindle holder
30	080561-25	Arandela pata	Leg washer
31	080080-26	Pata	Leg
32	080280-25	Tornillos alineación	Alignment screws
33	080561-27	Tornillos sujeción cilindro	Cilinder cleamping screws
34	080561-28	Manómetro	Manometer
35	080561-29	Control	Control
36	000200	Cable	Cable
37	080251-07	Racor entrada aire	Air intel connector
38	080400-07	Interruptor	Switch
39	080490-06	Base con porta fusible	Base with fuse holder
40	000260	Fusible	Fuse
41	080400-08	Ventilador	Fan
42	080561-31	Imán puerta	Door magnet
	020210	Mufla inyectora para nylon	Injector flask for nylon

- A Manubrio de husillo / **Screw handle** / *Axe de serrage du moufle*
- B Cámara de inyección / **Injection Chamber** / *Chambre d'injection*
- C Control electrónico / **Electronic Control** / *Contrôle électronique*
- D Manómetro / **Pressure Gauge** / *Manomètre*
- E Manorregulador / **Pressure Control** / *Manorégulateur*
- F Interruptor general / **Main switch** / *Interrupteur general*
- G Entrada de aire / **Air inlet** / *Entrée d'air*

- 1. Display / **Display** / *Affichage*
- 2. Tecla de Start-Stop / **Start-Stop Key** / *Touche Start-Stop*
- 3. Tecla subir / **Up key** / *Touche +*
- 4. Tecla bajar / **Down key** / *Touche -*
- 5. Tecla de acceso al modo de programación / **Program mode access key** / *Touche d'accès au mode de programmation*
- 6. Tecla Pistón abajo / **Piston down key** / *Touche piston en bas*
- 7. Tecla Pistón arriba / **Piston up key** / *Touche piston en haut*

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The Resin Injection Machine R-080560 has been designed for Dental Lab use. It is useful to make dental prosthesis with polyamide base (Nylon) by injection of material inside 56 calibre cartridges (56 sixty-fourth of an inch = 7/8 inch = 22,25 mm.). For an optimal working of the machine and an adequate performance we advise you read carefully this Instruction Manual.

INSTALLATION OF THE RESIN INJECTION MACHINE

1. Carefully unpack the Resin Injection Machine. Inside you will find two flasks, a pipe to connect it to the air supply system, a silicone anti-seizing up spray tin and an Allen spanner to tighten the flask screws.
2. Choose a flat and firm surface far from heat sources or vibrations. Make sure that the ventilation grilles at the back of the machine are free. Keep a distance of at least 15 cm from the wall or any other object.
3. Connect the machine to 230 V, 50/60 Hz power point equipped with earth return. It is also necessary to be protected with a magneto thermal element. Consult your installer.
4. Connect the machine air inlet to the compressor outlet with at least 6 or 7 working press Bars. Minimum volume of flow should be 80 litres/minute and it is advised to use dry filtered air.
5. Now the machine is ready to work.

INSTALLATION ON A WALL (page 23)

A kit for the installation on a wall is included with the machine. It consists of 3 rubber bands, one support, 4 rawls and 6 screws.

1. Stick the three rubber bands on the back of the machine.
 2. Draw two marks on the wall with an awl. Distance between centres should be 116 mm. Make 2 initial holes with a 3 mm bit.
 3. Drill the holes again with a 6 mm bit. Put the 10 mm rawls.
 4. Screw the two Ø6 screws. Distance between head screw and wall should be 9 mm.
 5. Put the support in the lower part of the machine using screws M4x10.
 6. Hang the machine on the wall and draw the position of the support two holes.
 7. Remove the machine and drill two holes with a 6 mm bit.
 8. Put the 6 mm rawls.
- Hang the machine and fix it to the wall by screwing the two Ø4.2x25 mm screws. Distance between the machine and the wall should be 10 mm.

STARTING UP

1. Once the machine is connected to the air supply, pull the pressure control knob and turn it till pressure reaches 6/10 Bars. (We advise using maximum pressure that the compressor can supply).
2. Make sure that the door of the injection chamber is closed and switch on the machine. After a few seconds the display will light up and the caption “TALLERES MESTRAITUA” will be shown.
3. Then, the machine will make a test of piston movement. After this, the injection machine will remain in stand by position.

BASIC OPERATION WITH DEFLEX TYPE CARTRIDGES

Any operation cycle of the Resin injection machine has four different steps:

1. Warming up time that the machine needs to reach working temperature.
2. Warming up time that the material inside cartridge needs to reach injection temperature
3. Injection step, time that the piston needs to insert the resin inside the flask with enough pressure.
4. Cooling down time so that the material solidifies copying the casting silhouette.

Therefore, a working cycle is defined by three variables that are defined by the user:

1. Working temperature. It generally ranges from 260 to 280 °C.
2. Cartridge warming up time. Time needed for the cartridge to reach injection temperature.
3. Cooling down time needed so that material inside the flask solidifies and we get the piece required.

Resin Injection Machine R-080560 has 10 programmes that can be edited by the user. The first five programmes are the ones advised by the manufacturer for using standard cartridges of these types: “M10 XR”, “CLASSIC SR”, “SUPRA SF”, “ACRILATO”, “FLUENCE” (see parameters in section “4 Edition of Programmes”). These programmes can be edited easily by the user. You may find hereafter how to carry out a basic working cycle using standard programmes advised by the manufacturer:

1. Starting from Stand-by position, press up/down keys to select the programmes appropriated for the cartridge we are going to use. There will always appear working temperature, warm up time and cool down time for the selected programme on the lower line of the display in Stand-by position.
2. Press Start/Stop Key for a few seconds. The display will show “Calentando” (warming up) and also the machine temperature.

- When the machine reaches the desired temperature, it emits an acoustic signal and it will appear the caption “Ponga Cartucho” (put cartridge). Open the door and place the cartridge on the injection chamber. In order to avoid possible seizures, we recommend sprinkle the surface of the cartridge with silicone lubricating spray tin (a tin supplied with the

(Fig. 1) Detail of how to place the flask

- machine). Take care not to sprinkle the part of the cartridge that is in contact with the flask in order not to contaminate the material when it is injected.
- Close the door and press Start/Stop key. “Calentando Cartucho” (warming up cartridge) and the time remaining to end this phase will appear on the display.

(Fig. 2) The screw handle must be completely up when removing the flask

- When this time finishes, the machine emits an acoustic signal and the caption “Inyección lista” (injection ready) will appear. Open the door of the injection chamber and put the flask inside: place the opening of the flask in the heating chamber orifice and tighten the screw handle till it is firmly secured. The correct position to place the flask is as in Fig. 1. In order to avoid negative effects of the thermal shock during injection, we recommend preheating the flask before putting it inside the injection chamber. Optimal preheating temperature is between 70 and 80 °C.
- Once the flask is correctly place, close the door and press Start/Stop key. “Inyectando” (injecting) will be shown on the display and the piston will raise. The machine will have the piston up while injecting for the time programmed. Once the time elapses, and acoustic signal will be emitted and the caption “Inyección finalizada” (injection finished) will appear, and the piston will return to low position. Then press Start/stop key and the machine will return to Stand-by mode, and it will be ready for another cycle.
- Open the door of the injection chamber, loosen the screw that holds the flask, wait 10 minutes and remove the flask. The screw must remain completely up so that it is hold by the second thread to the body of the machine (See Fig. 2) Sometimes the flask could get stuck and so it cannot be removed easily by hand. In this case, you should go to manual operation. Make sure that the **LID OF THE INJECTION CHAMBER IS CLOSED**, then press piston up key so that the force that makes the pneumatic cylinder releases the flask. **MAKE THIS**

CAREFULLY. (See chapter “5. Manual Operation Mode”).

While injecting, an air-cooling pneumatic system operates so it is normal to hear a slight air leak inside the machine. The purpose of this device is to cool down piston rod to protect it against overheating.

Once the cycle starts it can be cancelled at any step by pressing Start/Stop key for a few seconds. Then the caption “Ciclo Abortado” (cycle cancelled) appears and the machine returns to Stand-by position.

If there is an electrical cut while running a cycle, this can affect the machine in two different ways:

- If it is a short cut (it lasts less than 10 to 15 seconds), the machine will keep on running the cycle as if no power cut occurred.
- If the power cut lasts more than 15 seconds, then when electricity comes back, the machine emits an acoustic signal and the display will show an message about the error. Then press the Start/stop key to return to stand-by position.

PROGRAMME EDITION

Parameters of programme cycle (working temperature; heating time; injection time) can be easily changed to adapt these to the characteristics of a new material, or to adapt to needs and tastes of the user. In order to edit information about a programme, follow these steps starting from stand by position:

- Select with up and down keys the program that you want.

2. Keep key “key” till “P” of the selected programme is shown intermittent (about 2 seconds) Then the programme can be edited.
3. The next step is choosing editing the name of the programme (this is, add an alphanumeric label that identifies that programme), or editing the three numerical parameters that define the programme. In order to do this, press “key” key again as indicated below:
 1. If the keystroke is simple, numerical parameters that define a programme are edited one by one.
 - a. With the first simple keystroke of the “key” key, temperature value blinks. Press up and down keys to select temperature.
 - b. Press again “key” key and value of cartridge heating time starts blinking. Press un and down keys to select value for that parameter.
 - c. Press once again “key” key and injection time will blink. Do as before to change that parameter.
 2. If you press the “key” key for 2 seconds, you choose editing the Name of the programme. You go from one letter to the next one with simple

keystrokes. After entering the name, press “key” key for 2 seconds to get out of editing mode.

- a. Once the letter to edit is selected, use up/down keys to modify these values.
- b. If you want to get out of editing mode without validating the value or letter, press any of the other keys that are not used in editing. The machine will also return to stand-by position if no key is pressed in 30 seconds.

The five programmes of the machine are protected. Protected programmes are identified with a dot at the right of the number of the programme, (for example: “P2.”). If the selected programme is protected, up/down keys do not modify the value in editing mode, so it cannot be edited in a simple way. If you want to modify these protected values, keep “up” or “down” key pressed and then press another key. You press “up” or “down” depending if the value is to be increased or decreased.

The machine is factory programmed with the following parameters for cycle:

<i>Programme number</i>	<i>Programme name</i> (maximum 12 characters)	<i>Set temperature</i> 50°C to max limit of set temperature	<i>Heating time of cartridge</i> 1 to 60 minutes	<i>Injection time</i> 1 to 60 minutes
0	'D. M10 XR'	300	15	2
1	'D.CLASSIC SR'	280	15	2
2	'D. SUPRA SF'	260	15	2
3	'D. ACRILATO'	265	15	2
4	'D. FLUENCE'	275	15	2
5	'Ciclo 5'	280	13	2
6	'Ciclo 6'	280	13	2
7	'Ciclo 7'	280	13	2
8	'Ciclo 8'	280	13	2
9	'Ciclo 9'	280	13	2

MANUAL OPERATION MODE

This operation mode allows carrying out piston movement manually at any step, starting from stand-by position or while running a cycle. Besides, it allows going to the next step of a cycle even if the condition of step change does not take place. This mode is normally used when the piston gets blocked for any circumstance, when testing new materials without knowing necessary heating or injection time or when the flask is stuck after injection and cannot be taken out from the injection chamber easily using both hands.

IMPORTANT: in manual operation mode, injection chamber door must be closed for the piston to move. WITH OPEN DOOR, THE PISTON WILL NOT WORK.

In order to enter this mode, press for a second any of the piston keys: “Piston up” o “Piston down”, (the machine

can be in Stand-by position or running a cycle, but you cannot enter this mode while injecting). When entering or getting out of this operation mode, the machine emits an acoustic signal and the display will show 'Piston manual' (manual piston) and '!! PRECAUCION !!' (precaution) alternately.

In manual operation the piston is moved by pressing “Piston up” or “Piston down”. A cycle can be carried out in manual operation. In this case, the cycle is carried out as in automatic operation mode but with the following differences:

- You can change to the next step of the programme by pressing Start/stop key. For example, when carrying out a cycle in manual operation, if the machine is heating the cartridge, after pressing Start/stop key, it will do directly to the next step that is injection.
- While injecting, injection timer remains stop if the piston is down. When the piston is up, injection time is deducted.

Press again “key” key in order to get out of manual operation. If the machine is running a cycle when “Key” is pressed, the cycle will be aborted (a beep will be emitted and the message 'Ciclo abortado' (aborted cycle) will be shown for 5 seconds), except if the cycle is in the last step (end of injection shown on the display). Manual operation also ends when there is a long power cut.

EXAMPLE OF PROGRAMME

An example of how to program the Resin Injection Machine R-080560 and then how to run a cycle is described below: Suppose you want to programme program 5 with a cycle that has the following parameters:

- Temperature 265 °C.
- Heating time 25 minutes
- Injection time 5 minutes.

1. From stand-by situation, press up/down keys till the display shows programme 5 (P5)
2. Then press “key” key for a few seconds till “P” letter of the programme blinks.
3. Press again “key” key. Temperature value will start blinking. Use up and down keys to select 265 °C.
4. Press again “key” key. This time heating time will blink. Use up and down keys to select 25 minutes.
5. Press once more “key” key and injection time will blink. Use up and down keys to select 5 minutes.
6. After pressing “key” key again, the machine will return to Stand-by mode.

The programme is now ready to work doing the following:

1. Starting from Stand-by mode, press up/down keys till programme P5 is shown on the display.
2. Press Start/Stop key for a few seconds. On the display the caption “Calentando” (heating) will appear.
3. When the machine reaches set temperature, it will emit an acoustic signal and the caption “Ponga Cartucho” (put cartridge) will appear.
4. Close the door and press Start/Stop key again. The display will show “Calentando Cartucho” (heating cartridge) and the time that remains till that step ends.
5. When time elapses, the machine will emit an acoustic signal and the display will show “Inyección lista” (injection ready). Open the door of the injection chamber and place the flask inside.
6. Then close the door and press Start/Stop key. The display will show “Inyectando” (injecting) and the piston will raise.
7. Once injection time elapses, the piston will return to down position and the display will show “Inyec. Finalizada” (injection finished). Press any key to go back to Stand-by mode.

WARNING OR ERROR MESSAGES

These messages of warning or failure are usually shown while running a cycle. As they appear the machine stops and remains in warning situation till the user notices it. These failure situations are identified by: corresponding message on the display; flashing light on the display and intermittent beep. As soon as any key is pressed, warning situation is taken as notified and warning message disappears, light stops flashing, the beep silences and the machine returns to stand-by situation. There are some situations in which the machine does not get in warning situation as described above:

- Warning of aborted cycle. There is no need of recognition of the user. It emits a warning beep for 5 seconds.
- Failure of sensor or Failure of detectors, when they fail out of the cycle. In these cases, there is no need of recognition of the user and there is no light and no acoustic signal, it is just indicated on the display till the failure is solved.

The following are some of the different messages of warning and failure:

1. FAILURE OF POWER SUPPLY “Fallo de red”- It appears when there is a power cut for more than 15 to 20 seconds while a cycle is running except in the step of “inyección finalizada” (injection finished) warning.
2. ABORTED CYCLE “Ciclo abortado”- These warning is notified to the user and needs to be recognized. It is a beep for 5 seconds that takes place if:
 - A cycle is cancelled by pressing Star/stop key for 2 seconds, except in the step of “inyección finalizada” (injection finished)
 - A cycle is cancelled pressing any key for 2 seconds while injecting.
 - Manual Piston Mode is cancelled during a cycle pressing the key “key”, except if the cycle is in the last step (Aviso inyección finalizada: end of injection warning).
 - 30 minutes passes from “ponga cartucho” (put cartridge) warning o “inyección lista” (injection ready) warning.
3. SENSOR FAILURE “A. Sonda”- It is produced when the machine detects failure in temperature sensor.
4. HIGH TEMPERATURE FAILURE “A. temperat.alta”- It is produced if:
 - Temperature exceeds for 20 seconds temperature maximum limit.
5. LOW TEMPERATURE FAILURE “A. temperat.baja”- It is produced if:

- It is heating and Vacuum Heating is not passed in 45 minutes, this is, the machine does not reach set temperature in that time.
 - Temperature goes down temperature limit for 20 seconds in the following steps: “Ponga cartucho” (put cartridge) warning, “Calentamiento cartucho” (heating cartridge) warning or “Inyección lista” (injection ready) warning. (In “Manual piston” mode this situation is not taken into account).
6. **BLOCKED PISTON FAILURE** “A.bloqueo piston”- It is produced if:
- After 4 seconds in Vacuum heating step, the machine does not detects piston down (In “Manual piston” mode this situation is not taken into account).
 - After 4 seconds of the Injection step, the machine detects that the Piston is down when in normal conditions it must be up (In “Manual piston” mode this situation is not taken into account).
 - After 4 seconds at the end of injection step, the machine detects piston up, when it has to be down (In “Manual piston” mode this situation is not taken into account).
 - In the Piston Test, the correspondent detector (piston up, piston down) is not activated after 4 seconds from pressing the correspondent key.
7. **DETECTOR FAILURE** “A. detectores “ The machine has electronic detectors that control at any step the position of the pneumatic cylinder (Piston up; Piston down). There is a failure warning if:
- Piston down is not detected except while injecting or during the first 4 seconds of the Vacuum heating step (In “Manual piston” mode this situation is not taken into account).
 - Piston down is detected while injecting except the first 4 seconds of this step (In “Manual piston” mode this situation is not taken into account).
 - Piston up and Piston down are detected at the same time
4. Lubricate all surface around the cartridge with the silicone spray before putting it inside the machine. This avoids possible piston blockings. However, make sure you do not lubricate the part of the cartridge that gets in contact with the flask to avoid contamination.
5. In case the piston gets blocked, check working pressure (it should be 6 or 7 bars at least). Then in Manual Mode move the piston alternatively. **MAKE THESE OPERATIONS WITH INJECTION CHAMBER DOOR CLOSE.** If the piston is still blocked, ask an authorized distributor.
6. Outer body can be cleaned with a soft polish dampened cloth. A good option is using products for car dashboard cleaning.

PRECAUTIONS

1. Before connecting the machine, make sure it is a 230 V, 50/60 Hz electric point with earth. It is also necessary to be protected with a magneto-thermal element. Ask your installer.
2. Do not let children or not qualified people handle the machine.
3. When the piston is working, keep closed the door of the injection chamber.
4. Make sure that the flask is firmly secured with the screw handle once it is put inside the injection chamber.
5. Use security gloves to place and get out flasks from the injection chamber.
6. If the piston gets blocked in normal operation, disconnect the machine from the mains and air supply beforehand.
7. If the flask gets blocked inside the injection chamber, start moving the piston manually **WITH THE CHAMBER DOOR CLOSED.** Make sure the screw handle is completely up (lower screw, see Fig. 2) If further doubts, ask your MESTRA distributor.
8. Air pressure should not be lower than 6 bars or above 10.
9. Do not operate the piston (manual mode) with injection chamber door open.
10. Leave ventilation grilles free and at least at 15 cm from the wall or any other object. Place the machine on a firm well levelled surface so it does not turn over.

MAINTENANCE AND CLEANING

To ensure an optimum reliability and long life of the machine, we recommend::

1. Clean the heating chamber inside periodically. Before doing that, unplug the machine form the mains and air supply. Clean all plastic remains, aluminium casings and all dirt. You can use tweezers but do it carefully.
2. Clean all plastic remains, aluminium casings and all dirt that may have fallen inside the machine. In order to do this, unplug the machine, unscrew the back lid of the machine.
3. We recommend using dry and filtered air. If any doubt, ask your pneumatic installer.

TECHNICAL DATA

Height: 840 mm (screw handle up)
 Width:240 mm
 Depth:255 mm
 Weight:29 kg
 Pressure: 6/10 bars
 Voltage: 230 V, 50/60 Hz
 Consumption: 1100 W

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