# SOPLETE DE LABORATORIO Ref. 100002

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Rev. 14/02/2023

# MICRO TORCH Ref. 100002

Please read the following instructions thoroughly before using the torch.

## A. Ignition:


- The air adjustment switch must aim along with the holes before igniting the torch.
- Turn the ignitor button counter clockwise until you can hear the sound of gas flow.
- Press the ignitor button to light the torch.
- To turn off the torch simply just tighten the ignitor button in clockwise direction.

## B. Flame Adjustment:

- The flame of the torch or the temperature of the soldering head may be adjusted in two ways:
  1. Turn the ignitor button clockwise to reduce the gas flow, thus, decrease the flame strength. Turn the ignitor button counter clockwise would give contrary result.
  2. Rotate the air adjustment switch to reduce the flame strength.

## C. Refilling Gas Tank:

- Turn off the torch and make sure the fire is off. Wait until the torch temperature is cool down before refilling gas.
- Hold the butane can nozzle into the torch refilling valve.
- Stop refilling once gas starts overflowing. Please allow a few minutes for gas to stabilize.
- When refilling, keep away from fire and make sure that room is well ventilated.

## PRECAUTIONS

1. The operation of recharging the lighter must be carried out with special caution: use suitable recharging gas for these devices; recharge in a well-ventilated place AWAY FROM ANY SOURCE OF HEAT; and insert the refill nozzle into the lighter tank so that there is no gas leak.
  2. Sometimes (especially depending on the lighting), the flame of the lighter can take on a light blue hue that is very difficult to perceive. The user must remain alert and make sure the burner is lit, even if the flame cannot be seen.
- Improper use may cause a fire or explosion, resulting in personal injury or damage to objects and facilities.
  - Make sure there are no flammable residues on the outside of the appliance. Wipe the exterior with a rag before use.
  - Use only in well-ventilated areas.
  - The use of fire retardant gloves, protective clothing, goggles and face protection is recommended as appropriate. Anyone who is near the lighter while it is in use must use the same personal protection elements that are appropriate for the user.
  - The lighter must always be kept and handled in an upright position.
  - Before use, check that there are no gas leaks and that the ignition nozzle is in good condition.
  - After using the lighter, make sure it is completely extinguished. Never leave the lighter burning or leave it alone on the table or any other surface.
  - The lighter should only be handled by qualified personnel in a dental laboratory environment. It is recommended to have a nearby water source and, if possible, a fire extinguisher prepared for type "B" fires (flammable liquids).
  - In case of inflammation, the best way to put out the fire is by covering the flames with a blanket or resistant fabric. Never work near flammable substances.
  - When operating with the lighter, be careful with your clothing (ties, handkerchiefs, etc.) and with long, loose hair.
  - In case of burns, wet the burned area with cold water (tap stream) and apply water-soaked compresses. Go to a medical center.

## TECHNICAL DATA

Dimensions	36 mm x 142 mm x 82 mm
Temperature	1300 °C (2370 °F)
Soldering iron	200~450 °C (390~840 °F)
Capacity	Butane 64 ml
Burning time	2-3 h / each full filled

**MESTRA**®

*Su marca de confianza*

Talleres Mestraitua S.L.  
Txori-erri Etorbidea, 60  
48150 SONDIKA - BILBAO - SPAIN  
T: (+34)944530388 - F: (+34)944711725  
info@mestra.es - www.mestra.es