R-100006

E.....	ESPÁTULA ELECTRONICA PARA PARALELÓMETRO.....	4
GB.....	ELECTRONIC SPATULA FOR PARALLELOMETER.....	6

MESTRA®



TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

48150 SONDIKA - BILBAO - ESPAÑA

Tfno. + 34 944530388 - Fax + 34 944711725

E-mail: mestra@mestra.es - www.mestra.es

A	Adaptador / Adaptor
B	Caja de control / Control box
C	Clavija / Plug
D	Teclas ▼▲ / ▼▲ keys
E	Display / Display
F	Interruptor / Main switch
G	Eje de fijación / Fixing shaft
H	Puntas / Tips

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- Carefully read this manual. It contains important information on safety during the installation, operation and maintenance of the unit.
- Keep this manual in a suitable place to consult it in the future.
- Check the overall state of the unit just upon unpacking it.
- Before connecting the unit, make sure the voltage shown in the identification plate matches that of the power source to be used. The installation must fulfil all the applicable safety standards and the instructions of the manufacturer.
- The manufacturer will not be responsible of the damage caused by an incorrect ground connection. The damage caused by a faulty installation or usage of the unit will not be covered by the warranty.
- In case of doubt, do not start the unit. Contact a qualified and authorised technician.

- Do not touch the unit or the power cord with wet hands.
- Unplug the unit in case of storm and when you are not going to use it for a long time.
- Place the unit in a dry place, away from dust and water splashes.
- Do not pull, twist or fold the power cord excessively.
- Do not wet the unit; do not use chemical products to clean it.

- **WARNING:** to prevent electric shock hazards, connect the unit to a proper power source (grounded 230 V, 50/60 Hz).

- Do not touch the spatula heads while the unit is in operation, since they reach high temperatures.
- After turning the main switch off, wait a time before touching the spatula heads.

- The machine features a precision circuit. Prevent fortuitous impacts and stresses that could cause mechanical damage.

- Remove any dirt or water remains from the cord connection.
- Unplug the unit immediately whenever you notice an unusual noise, smell or smoke.
- Do not insert plastic objects into the bowls or any other material that melts at low temperatures.
- Do not use the unit in damp rooms or near inflammable materials (alcohol, solvents,...).

- Do not allow children or non-qualified personnel handling the unit.
- Choose a flat and firm surface to place the unit, away from heat or vibration sources.

- Contact the supplier for questions or detailed information about problems with the unit.

INTRODUCTION

The R-100006 is an electronic spatula for wax modelling specially designed to be used by professionals in the field of dental prosthesis. Its ergonomics, easy handling and small size make it an ideal element for this kind of work. The high accuracy of the electronic controls allows to maintain the temperature of the spatula tip within a narrow error range. The working temperature is shown in the digital display in the control box.

It is an ideal device to be attached to a parallelometer or similar device for wax modelling. The device is factory supplied with a $\varnothing 2$, $\varnothing 2.35$ or $\varnothing 3$ mm fixing shaft. It includes a set of 4 different tips.

INSTALLATION

- Place the control box (B) on a flat, horizontal and firm surface, away from vibration or heat sources. Make sure it is away from flammable materials.
- Select one of the tips (H) and screw it into the adaptor (A).
- Fit the adaptor to a parallelometer, micro-milling support, etc, and fix it to the shaft (G).
- Connect the plug (C) to the control box.
- Connect the unit to a grounded 230 V, 50/60 Hz power source.

OPERATION

- Turn on the main switch (F) in the right side of the control box (B).
- Press the ▼▲ keys (D) to select in the display (E) the temperature you wish for the tip to reach.
- After finishing the work, disconnect the device by turning off the main switch (F), and handle the tip very cautiously for several minutes until it has cooled down completely.

MAINTENANCE AND CLEANING

To ensure a long working life of the device and an optimum performance, we strongly recommend:

- Remove the wax remains from the device periodically. To do this, you can use a small spatula. Do not forget to disconnect the device beforehand.
- Do not spill melted wax over the casing of the device. Handle the spatula with caution.

SPECIFICATIONS

Voltage:	230 V, 50/60 Hz
Power:	12 W
Fuse:	1 A
Temperature range:	50 ~ 200 °C
Height x width x depth (control):	65 x 75 x 120 mm
Fixing shaft:	$\varnothing 3 / 2.35 / 2$ mm
Height of adaptor:	23 mm
Weight of control box:	650 g