R-100200

E..... ELECTROVIBRADOR	4
GB..... ELECTRO-VIBRATOR.....	5

MESTRA®

TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tfno. + 34 94 453 03 88 Fax + 34 94 471 17 25

E-mail: mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

Rev. 26/02/14

A

B

C

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The R-100200 is a practical electro-vibrator unit specially designed for compacting materials for the jewellery and dental mechanics sectors. To obtain best performance from the unit, we recommend you read the following instructions carefully.

- Place the unit on a smooth, firm, perfectly horizontal surface. The electro-vibrator is fitted with rubber suction caps (A) to ensure firm attachment.
- Connect the plug to a 230 V, 50/60 Hz earthed power outlet.
- Place the cylinder or items to be vibrated on the unit's vibration platform. In the case of impression bowls, the mixture can also be compacted by resting the bowl firmly on the edge of the unit and holding down the bowl handle with your hand.
- Adjust the vibration regulator (B) to its minimum setting (turn the knob anti-clockwise).
- Slowly move the unit's regulator clockwise and check that the oscillation rate gradually increases. Set the regulator to a suitable amplitude for the operation.
- After a short time, once the mixture is properly compacted, switch the unit off by turning the regulator to the left.

MAINTENANCE AND CLEANING

To ensure a long, trouble-free service life and performance of the unit, we recommend as follows:

- After using the unit, always clean away any splashes of mixture before they dry.
- Clean the vibration platform and outer housing of the electro-vibrator regularly with a damp cloth. The rubber cover protecting the platform can be easily removed to facilitate cleaning.
- As a precaution, always disconnect the unit from the mains before cleaning it.

PRECAUTIONS

- ◆ Before connecting the unit ensure that the power supply is 230 V, 50/60 Hz and earthed.
- ◆ Do not use the unit in wet places or close to heat sources.
- ◆ Do not place easily broken containers or other objects on the vibration platform.
- ◆ Do not allow children or unskilled personnel to handle the unit.
- ◆ Do not pull, twist or bend the mains cable excessively.
- ◆ Before cleaning the unit, ensure that it is unplugged.
- ◆ Never immerse the unit in liquid when cleaning it.
- ◆ Do not place hot containers or objects on the vibration platform, as this may damage the rubber.

TECHNICAL CHARACTERISTICS

Voltage:	230 V, 50/60 Hz
Dimensions:	Ø127 x 90 mm
Power:	110 W
Weight:	1.6 kg