
Ref. 100347

ES.....	3	Manual original
EN.....	7	

Uso previsto

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
Ignacio Mestraitua
Director Gerente
Sondika, 23.10.2023

Intended use

Position the mouthpiece next to the patient's mouth and prevent the spread of all types of aerosols, suspended particles, vapors and odors in the clinic. Greatly reduce the risk of cross-contamination between professionals and patients and protect respiratory health.

To achieve optimal performance and performance of the device, we recommend that you read the following instructions carefully:

Description

- Two power levels.
- 2 filter cartridges composed of:
 - a) Coarse filter for small particles.
 - b) Active carbon filter.
 - c) H13 HEPA filter.
- UV light.
- Plasma generator.
- Ergonomic shape and easy to move to have a good working position.
- Semi-rigid suction arm, self-supporting and easy to clean.
- Easy-to-clean housing.
- Hours of use counter.

Installation

- Carefully unpack the device.
- Choose a rigid, flat surface for the location of the device, away from sources of vibration or heat.
- Connect the electrical plug to a 230 V, 50/60 Hz power outlet with a ground connection.
- Install the arm into the plastic connection on the top of the vacuum body. Position the articulated tube so that the mouthpiece is over the patient's mouth.
- Handle the arms gently to avoid damage.

Operation

- Carefully unpack the device.
- Choose a rigid, flat surface for the location of the device, away from sources of vibration or heat.
- Connect the electrical plug to a 230 V, 50/60 Hz power outlet with a ground connection.
- Install the arm into the plastic connection on the top of the vacuum body. Position the articulated tube so that the mouthpiece is over the patient's mouth.
- Handle the arms gently to avoid damage.

UV lamp maintenance

The expected life of the UV lamp is approximately 10 000 hours.

To replace it, remove the back panel and twist both ends.

Place the new lamp (Ref. 100347-03) and close the machine again.

Filter maintenance

Periodically, when you notice a decrease in suction capacity, proceed to clean/replace the filters. Replace them at least once a year.

Before opening the filter compartment, protect yourself adequately with gloves, a mask and safety glasses.

To do this, remove the rear panel and the next panel that closes the filter chamber.

Remove the two filter cartridges, being careful not to allow debris to fall on the motor.

Place the new filters (Ref. 100347-01) and close the machine again.

Plasma sterilizer maintenance

The expected life of the plasma sterilizer is approximately 20 000 hours.

Replacing the fuses

The suction has two 10 A fuses.

Use a Phillips screwdriver to open the fuse holders.

Precautions

- Keep this manual in an appropriate place for future reference.
- Check the general condition of the device immediately after removing the packaging.
- Failures caused by poor installation or misuse of the device are excluded from the warranty.
- To avoid the risk of electric shock, this equipment should only be connected to a suitable grounded power source. Before connecting the appliance, check that the identification label indicates the correct voltage. Installation must be in accordance with all applicable safety standards.
- Remove any dirt or water from the cable connection.
- Avoid pulling, twisting or excessively bending the network connection cable.
- Do not touch the appliance or the cable with wet hands.
- Install the product away from heat sources, dust, high humidity, water splashes and flammable materials (alcohol, solvents, etc.).
- Unplug the appliance immediately if you detect noise, odor or smoke.
- Prevent children or non-specialized personnel from manipulating the device.
- If in doubt, do not operate the device. Contact a qualified and authorized repairer.
- Do not disassemble, repair or modify the device beyond maintenance operations on filters, lamp, plasma generator and fuses.
- Unplug the appliance in the event of a storm and when not in use for a long time.
- Unplug the appliance before cleaning it, do not wet it directly, do not use chemicals.
- Contact your dealer for questions or detailed troubleshooting information.
- Before using the device, make sure that it is sitting on a firm and solid base.
- The cable can only be replaced by one with the same characteristics as that supplied by the manufacturer.
- Protect yourself adequately, use a mask, gloves and safety glasses.
- Do not use the machine outside its intended use.

Technical characteristics

Tension: 220 V, 50 Hz
Power: 1000 W
Two power levels: (50-100 %).
Noise: ≤ 60 dB.
Flow rate: 4000 l/min.
Máximum vacuum: 25 kPa
UV lamp: 254 nm (lifespan 10 000 h)
Plasma generator: 5000 V (lifespan 20 000 h)
HEPA grade: H13 (0.3 µm, > 99.97 %)
Dimensions: 32 x 33 x 75.5 cm (arm not included)
Weight: 30 kg (arm included)
Arm length: 1.3 m
More than 12 hours of continuous work without overheating

Declaration of conformity UE

The firm **Talleres Mestraitua S.L.**
with social address: Txori-erri Etorbidea, 60
48150 Sondika - Bilbao - Spain

Hereby declares that the product

Extraoral suction unit Ref. 100347

Meets the requirements set out in the following directives:

2014/35/UE (low tension)

2014/30/UE (electro magnetic compatibility)

Ignacio Mestraitua
Managing Director
Sondika, 2023-10-23

MESTRA[®]

Your trusted brand

Talleres Mestraitua S.L.

Txori-erri Etorbidea, 60
T: (+34)944530388 - F: (+34)944711725
mestra@mestra.es - www.mestra.es
48150 SONDIKA - BILBAO - ESPAÑA