R-050000/20 ; R-050030/50

ES	MECHERO DE GAS	4
EN	GAS BURNERS	5
FR	BRULEURS D' ETABLI A GAZ	6
PO	QUEIMADOR DE GÁS	7

MESTRA®


TALLERES MESTRAITUA S.L.


Txori-Erri Etorbidea, 60

Tel. +34 944530388 - Fax +34 944711725

mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

A	Boquilla de conexión a la red de gas / Connection nozzle to the gas manis / <i>Raccord d'admission du gaz</i> / Boquilha de ligação à rede de gás
B	Llave de paso de gas / Gas shutoff cock / <i>Vanne</i> / Torneira de gás
C	Válvula de regulación de paso de aire / Airflow adjustment valve / <i>Valve de régulation d'air</i> / Válvula de regulação de passagem de ar
D	Boquilla de llama / Flame nozzle / <i>Embout de flamme</i> / Boquilha da chama
E	Recipiente de cera / Wax tray / <i>Réipient pour la cire</i> / Recipiente de cera

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The R-050030/50 burner range is supplied with wax trays. The R-050000/20 range is identical, although this component is not supplied.

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EN

INSTALLATION

1. Remove the gas burner packaging materials with care. You will find a label on the rear side of the burner base describing the type of fuel to be used. Depending on the burner reference number, fuel options are as follows: Butane gas (R-050000/30) or natural gas (R-050020/50). **Make sure that the burner fuel coincides with your mains supply.**
2. Position the burner on a flat, solid surface in the laboratory, well away from inflammable materials.
3. Check the gas mains. Make sure that all the necessary inspections have been performed by licensed technicians.

4. Use a standard rubber gas tube to connect the burner to the gas input nozzle (A) and standard metallic clasps to secure this component firmly. Make sure that all clasps are completely tight.

HOW TO USE THE BURNER

1. If you need to work with the burner, place the wax tray (E) on the burner's flame nozzle.
2. Rotate the airflow adjustment valve (C) until the circular window is completely open. Adjust by turning the knurled cylinder on its axis.
3. Open the gas mains stopcock. Apply a ¼ turn to the burner's shutoff cock (E) until it is parallel to the central axis of the burner.
4. Bring up a flame to the nozzle (D). The burner should light up immediately without any problem.
5. Use regulator (C) to adjust the airflow. By opening the regulator completely you will obtain a flame with a blue base and a bright yellow tip, with a constant length and brightness. A flame of this kind means that combustion is perfect and that the gas supply is optimum. When regulator (C) is used to shut off the airflow, the blue color eventually becomes reddish, meaning that combustion is lacking oxygen. Under extreme circumstances, and when the airflow is too low, black (unburned) smoke may appear. **This situation should always be avoided, since the resulting gases are extremely dangerous.**
6. Close the gas airflow stopcock (B) by applying a ¼ turn until it is perpendicular to the burner's axis. Close the gas mains stopcock.

MAINTENANCE & CLEANING

1. Clean the burner externally on a regular basis. First make sure that the gas mains stopcock has been shut off completely. Use a cloth soaked in soapy water to clean the unit. Never use alcohol or inflammable products.
2. Gas facilities must be inspected regularly by licensed technicians.
3. During normal operation, molten wax may accidentally fall through the orifice in the flame nozzle (D). Wax accumulating inside the nozzle elbow may cause obstructions and impair burner performance. Should this occur, proceed to remove the wax as follows:
 - Close the mains stopcock.

- Release the metallic clasp on the connection nozzle (A). Disengage the burner from the rubber mains connection tube.
- Once the burner has been disengaged from the mains, apply a ¼ anticlockwise turn to the flame nozzle (D).
- Apply a flame to the nozzle elbow to melt the wax accumulated inside. Tilt the burner to allow the wax to flow out freely through the flame nozzle.
- After removing the wax, apply a ¼ turn to the flame nozzle to return to the initial position, and connect the burner to the gas mains as described above, in reverse order.

PRECAUTIONS

- ◆ Before connecting the burner, ensure that the gas supply is correct (butane, natural gas).
- ◆ Only connect the unit to a gas supply which has undergone official inspections performed by licensed technicians.
- ◆ Make sure that the rubber hose and connection clasp of the burner to the gas mains are solid and secure.
- ◆ Never allow children or unqualified staff to use the burner.
- ◆ Make sure the flame has a bright blue base and a strong yellow tip. Reddish colors, and black smoke in particular, indicate poor combustion. Always avoid this situation using the airflow valve.
- ◆ After using the burner, always close the shutoff cock and gas mains stopcock.
- ◆ Always position the burner in your laboratory well away from inflammable products.
- ◆ Never leave the burner operating unattended.
- ◆ If you smell gas, close the mains stopcock immediately and inform a licensed technician.

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