Ref. 080133, 080135, 080139

VULCANO

ES..... 6 Manual original
EN..... 9

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References Ref. 080133, Ref. 080135 and Ref. 080139 are three practical ovens specially designed for heating cylinders and crucibles within the dental prosthesis sector. The devices have a control panel to program a sequence of temperatures and times pre-established by the user. In addition, it is possible to store up to a total of 20 different programs in memory, so that it is simple and practical to recover previously made programs, significantly simplifying programming operations.

INSTALLATION

- Carefully unpack the oven. Inside the box you will find a SiO₂ hearth. Also check that the muffle that lines the interior of the furnace has not been damaged during transport.
- Choose a horizontal, flat and rigid base for the location of the oven, away from sources of vibrations. Make sure there is also enough space for ventilation.
- The closest element must be at least 15 cm from the appliance. Take special care that the ventilation tube is free and more than 20 cm from the nearest object or wall.
- During normal operation of the oven, smoke leaks may occur through the ventilation tube. For this reason, it is necessary to place the oven under an extractor hood.
- Connect the device to a 230 V, 50/60 Hz power outlet with a ground connection. It is also important that the outlet be protected by a circuit breaker of the appropriate amperage (consult your installer).
- Before using the oven, it is necessary to preheat it in order to evaporate any moisture that the muffle may contain. To do this, program the oven to remain at 300 °C for 1 hour.

USER SETTINGS

- From the home screen, press and hold the **P** key.

CUSTOMER CONFIG.

1.0 User

- Briefly press the **P** key again to enter option *1.0 User*.

1.0 User

TWARN = 10 "

LANGUA = English

UNITS = °C

SLOP= min

- By briefly pressing the **P** key you can select:

TWARN: End of cycle warning time

LANGUAGE: Spanish English, German, French

UNITS: Celsius, Fahrenheit

SLOP: Rising times in *min* o *°C/min*

- Use the keys **▲▼** to modify the values and the **P** key to save.

SMOKE EXTRACTOR

On the lower side of the appliance there is a 230 V power outlet, very useful when a smoke extractor is attached to the oven and which serves to power the extractor.

The use of Mestra extractors Ref. 080124 and Ref. 080125 is recommended, as they are specially designed for that function.

By default, the extractor connects at 75 °C and disconnects at 650 °C. To modify these temperatures, consult the manufacturer.

EDITING A PROGRAM

- Press the **P** key to access the program data.

P1 PROGRAMA-1

t0 : 01h50m

C2: 45°C

r1: 1°/min

t2: 10'

C1: 95°C

r3: 1°/min

t1: 30'

C3: 95°C

r2: 2 °/min

t3: 55'

- Use the keys **▲▼** to modify the values and the **P** key to save and go to another parameter:

- t0: Delayed start time (so that the machine starts working after a while).
- r1: Time to reach the 1st temperature
- C1: 1st temperature
- t1: Maintenance time of the 2nd temperature
- r2: Time to reach the 2nd temperature
- C2: 2nd temperature
- t2: Maintenance time of the 2nd temperature
- r3: Time to reach the 3rd temperature
- C3: 3rd temperature
- t3: Maintenance time of the 3rd temperatura

SELECTING A PROGRAM

Being on the program number we can select another one by pressing the keys **▲▼**.

OPERATION

- Place the SiO₂ hearth inside the furnace, taking special care not to damage the muffle.
- Activate the main switch.
- Introduce the cylinders and crucibles to be heated. To do this, always use gloves and tongs when the oven is in operation. TALLERES MESTRAITUA, S.L. recommends the use of **MESTRA** tongs.
- The oven is supplied from the factory with a series of standard programs in its memory. However, it is advisable that you program the oven following the instructions outlined above, to enter different programs that adapt to your own needs. In the future you can use or modify these programs whenever you want.
- Select the program you want to run.
- Press the Start/Stop button to start or stop the cycle.

- If necessary, a long press of the key **▲** allows you to jump to the next stage.

- In the event of a possible power outage, the electronic control keeps the parameters of the step in execution in memory, so when the power is restored again, it will continue its program from the stopping point. Also, if we disconnect the main switch before the program ends, the oven will consider it as a power failure.
- If the alarm is active, it will sound at the beginning of the third heating (if any) and at the end of the program.
- The oven will maintain the last temperature setpoint until the user stops the operation by pressing and holding the Start/Stop button
- When finished, remove the cylinders and crucibles with caution. During operation of the device, high temperatures can reach the outer casing. Furthermore, the temperatures of the cylinders and crucibles reach very high values. Always use gloves and tongs to carry out this type of manipulation.
- If after finishing a program you want to run another one, you must first stop the oven by holding down the Start/Stop button.
- To disconnect the oven from the mains, action the main switch.

ERROR MESSAGES

INTERNAL TEMP.	Overheating of the electronics has been detected. Get in contact with your distributor.
HIGH TEMP.	The safety temperature has been exceeded.
CHAMBER PROBE	Indicates failure of the temperature probe. If the problem persists, contact your dealer.
NET FAILURE	The machine warns that there has been a prolonged power outage. Although it will try to recover the program, it may have affected work.
TRIAC FAILURE	Contact your dealer.

MAINTENANCE AND CLEANING

To ensure optimal reliability and long life of the device, we recommend the following tips:

- Periodically clean the outside of the oven.
- Take care not to damage the muffle that covers the interior walls of the furnace hearth when inserting or removing cylinders and crucibles.
- Always use the SiO₂ hearth to avoid leaks that could damage the muffle.
- Prevent liquids or solid substances from dirtying or contaminating the furnace muffle.
- To achieve optimal performance of the device, TALLERES MESTRAITUA, S.L. recommends the use of cylinders and crucibles from the **MESTRA** series.
- The two circuit breakers located next to the cable protect the oven in case of electrical overload. To restore operation, once the problem has been solved, it will be necessary to reset the circuit breaker.
- If the equipment detects that the triac has broken, the upper display will show the message “*TRIAC FAIRULE*”. When this alarm signal appears, turn off the main switch and immediately contact technical service.

IMPORTANT

If during normal operation of the oven you notice that the upper display shows the message “F. CAMERA PROBE”, disconnect it from the network and consult an authorized MESTRA distributor. This symptom indicates that the thermocouple has suffered some damage. The thermocouple is an element subject to wear and its durability depends on the working temperature. This is an item not covered by the warranty.

PRECAUTIONS

- Before connecting the device, make sure that it is a 230 V, 50/60 Hz power outlet with a ground connection. It is also necessary that it be adequately protected with a thermal magnet (consult your installer).
- Do not allow children or unqualified personnel to handle the appliance.
- During oven operation, always keep the door closed.
- Open the oven door with caution. The sudden release of vapors could cause burns.
- Never overfill the oven, nor place objects that could fall when opening the door.
- Do not put flammable materials or liquids inside the oven. Also do not introduce substances that can produce toxic fumes.
- Place the oven on a flat, firm and vibration-free surface. The location must also have an extractor hood. Take special care that the oven is away from flammable materials.
- Place the appliance leaving the sides of the oven free and more than 15 cm from the wall or any other object. Take special care that the ventilation tube is free and more than 20 cm from the nearest object.
- Whenever you remove or insert cylinders and crucibles, use gloves and tongs. You should also avoid touching the metal casing of the oven. Its temperature can be high.
- Handle the materials you remove from the oven after heating, with extreme caution. Its extremely high temperature (up to 1100 °C) can cause accidents.
- Before using the oven for the first time, keep it for 1 hour at 300 °C, so that any moisture that may be contained in the muffle evaporates.
- Always use the oven under an extractor hood.

TECHNICAL DATA

	080133	080135	080139
Height x width x depth:	465 x 335 x 340 mm	500 x 375 x 405 mm	525 x 440 x 490 mm
Weight:	30 kg	35 kg	42 kg
Power supply:	230 V, 50/60 Hz	230 V, 50/60 Hz	230 V, 50/60 Hz
Power:	1500 W	2000 W	3000 W
Maximum temperature:	1100 °C	1100 °C	1100 °C
Muffle height x width:	100 x 160 mm	117 x 180 mm	149 x 234 mm
Muffle length:	165 mm	227 mm	305 mm
Capacity:	2,7 L	4,8 L	10,6 L