**Ref. 080142 / Ref. 080144 / Ref. 080148**

|          |                                  |    |
|----------|----------------------------------|----|
| ES ..... | HORNOS DE PRECALENTAMIENTO ..... | 6  |
| EN ..... | PRE-HEATING FURNACES .....       | 10 |
| FR ..... | FOURS DE CHAUFFE .....           | 14 |
| PO ..... | FORNOS DE PRÉAQUECIMENTO .....   | 18 |
| IT ..... | FORNI DI PRERISCALDAMENTO .....  | 22 |
| DE ..... | VORWÄRMÖFEN .....                | 26 |
| NO ..... | FORVARMINGSOVNER .....           | 30 |

**MESTRA®**


TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tel. +34 944530388 - Fax +34 944711725

mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA


**Ref. 080142 / Ref. 080144 / Ref. 080148**


|      |   |
|------|---|
| A    | Interruptor general / <b>On-Off switch</b> / <i>Interrupteur general</i> / Interruptor geral / <b>Interruttore generale</b> / <i>Hauptschalter</i> / På-Av-bryter                         |
| B    | Panel de control / <b>Control panel</b> / <i>Panneau de contrôle</i> / Painel de controlo / <b>Pannello di controllo</b> / <i>Bedientafel</i> / Kontrollpanel                             |
| C    | Display de temperatures / <b>Temperature display</b> / <i>Cadran de température</i> / Display de temperaturas / <b>Display temperature</b> / <i>Temperaturdisplay</i> / Temperaturdisplay |
| D    | Display de tiempos / <b>Time display</b> / <i>Cadran du temps</i> / Display de tempos / <b>Display tempi</b> / <i>Zeitdisplay</i> / Tidsdisplay   |
| E    | Puerta / <b>Door</b> / <i>Porte</i> / Porta / <b>Porta</b> / <i>Tür</i> / Dør   |
| F    | Hogar / <b>Hearth</b> / <i>Âtre</i> / Fornalha / <b>Focolare</b> / <i>Ofenkammer</i> / Brennkammer  |
| G    | Mufla / <b>Muffle</b> / <i>Moufle</i> / Mufla / <b>Muffola</b> / <i>Muffel</i> / Muffel   |
| H    | Tubo de ventilación / <b>Ventilation tube</b> / <i>Conduit d'aération</i> / Tubo de ventilação / <b>Tubo di ventilazione</b> / <i>Abzugsrohr</i> / Ventilasjonsrør                        |
| I    | Enchufe schuko / <b>Current pickup</b> / <i>Prise de courant Schuco</i> / Tomada Schuko / <b>Spina Schuko</b> / <i>Schukosteckdose</i> / Stikkontakt                                      |
| J, K | Disyuntor / <b>Circuit breaker</b> / <i>Disjoncteur</i> / Disjuntor / <b>Disgiuntore</b> / <i>Ausschalter</i> / Overlastbryter  |

**Fig. 1**

1. Display de temperaturas / **Temperature display** / *Cadran des températures* / Display de temperaturas / **Display temperature** / *Temperaturdispla* / Temperaturdisplay
2. Display de tiempos / **Time display** / *Cadran des temps* / Display de tempos / **Display tempi** / *Zeitdisplay* / Tidsdisplay
3. Funcionamiento de resistencia / **Resistor operation** / *Témoin de fonctionnement de la résistance* / Funcionamento da resistencia / **Funzionamento resistenza** / *Widerstandsbetrieb* / Indikatorlampe for v.element
4. Aviso fin de programa / **Program end warning** / *Témoin du signal acoustique* / Aviso de fim de programa / **Avviso fine programma** / *Hinweis auf Programmende* / Varsellampe f. prog.slutt
5. Tecla stop-start / **Start-stop key** / *Touche Marche-Arrêt* / Tecla stop-start / **Tasto stop-start** / *Start-Stopp-Taste* / Start-stopp tast
6. Tecla escape - conmutación displays / **Escape - display switch key** / *Touche déconnection-commutation des cadrans* / Tecla escape - comutação displays / **Tasto escape - commutazione display** / *Escape-Taste – Displayumschaltung* / Escape - display tast
7. Aviso funcionamiento de aparato / **Unit-on lamp** / *Témoin de fonctionnement du four* / Aviso de funcionamento do aparéelo / **Avviso funzionamento apparecchio** / *Hinweis auf Gerätebetrieb* / Driftsindikatorlampe
8. Decrementar / **Down** / *Touche –* / Decrementar / **Diminuire** / *Verringern* / Ned
9. Incrementar / **Up** / *Touche +* / Incrementar / **Aumentare** / *Erhöhen* / Opp
10. Tecla de programación / **Programming key** / *Touche de programmation* / Tecla de programação / **Tasto di programmazione** / *Eingabetaste* / Programmeringstast
11. Tiempo de pasos / **Step time** / *Temps de chaque étape* / Tempo dos passos / **Tempo passi** / *Schrittzeit* / Platåtid
12. Temperatura de pasos / **Step temperature** / *Température de chaque palier* / Temperatura dos passos / **Temperature passi** / *Schritttemperatur* / Platåtemperatur
13. Interruptor general / **On-off switch** / *Interrupteur général* / Interruptor geral / **Interruttore generale** / *Hauptschalter* / På-Av-bryter

| Nº | Denominación                | Description               | HP-25     | HP-50     | HP-100    |
|----|-----------------------------|---------------------------|-----------|-----------|-----------|
| 1A | Carcasa tapa trasera        | Back lid housing          | 080142-31 | 080144-29 | 080148-23 |
| 1B | Carcasa tapa trasera mufla  | Muffle back lid housing   | 080142-32 | 080144-30 | 080148-24 |
| 1C | Carcasa mufla               | Muffle housing            | 080142-33 | 080144-31 | 080148-25 |
| 1D | Carcasa cableado            | Hardwired housing         | 080142-39 | 080144-02 | 080148-26 |
| 1E | Tope de puerta              | Door stop                 | 080142-35 | 080142-35 | 080142-35 |
| 1F | Carcasa marco               | Frame housing             | 080142-37 | 080144-35 | 080148-29 |
| 1G | Carcasa puerta              | Door housing              | 080142-36 | 080144-34 | 080148-28 |
| 1H | Carcasa laterales           | Sides housing             | 080142-38 | 080144-36 | 080148-30 |
| 2  | Tornillos M4x8              | M4x8 screws               | 080422-29 | 080422-29 | 080422-29 |
| 3  | Mufla cerámica              | Ceramic muffle            | 080142-03 | 080144-03 | 080148-03 |
| 4  | Placas aislantes            | Insulation boards         | 080142-04 | 080144-04 | 080148-04 |
| 5  | Solera                      | Hearth plate              | 080142-05 | 080144-05 | 080148-05 |
| 6  | Separadores aluminio        | Aluminium separator       | 080142-06 | 080142-06 | 080142-06 |
| 7  | Arandela de imán            | Magnet washer             | 080142-07 | 080142-07 | 080142-07 |
| 8  | Tornillo arandela imán      | Magnet washer screw       | 080523-07 | 080523-07 | 080523-07 |
| 9  | Tornillo sujeta imán        | Magnet supporting screw   | 080142-09 | 080142-09 | 080142-09 |
| 10 | Caja conexiones             | Connection box            | -         | 080144-38 | 080144-38 |
| 11 | Imán                        | Magnet                    | 080142-11 | 080142-11 | 080142-11 |
| 12 | Muelle imán                 | Magnet spring             | 080142-12 | 080142-12 | 080142-12 |
| 13 | Tapa cerámica               | Ceramic lid               | 080142-13 | 080144-13 | 080148-13 |
| 14 | Fibra de puerta cerámica    | Ceramic door fiber        | 080142-14 | 080144-14 | 080148-14 |
| 15 | Tornillo varilla de bisagra | Hinge stick screw         | 080251-09 | 080251-09 | 080251-09 |
| 17 | Bisagra                     | Hinge                     | 080142-17 | 080142-17 | 080142-17 |
| 18 | Varilla de bisagra          | Hinge stick               | 080142-18 | 080144-18 | 080148-18 |
| 20 | Tornillo hexagonal puerta   | Door hexagonal screw      | 080142-20 | 080142-20 | 080142-20 |
| 21 | Tuerca hexagonal puerta     | Door hexagonal nut        | 030230-04 | 030230-04 | 030230-04 |
| 22 | Manilla de baquelita        | Bakelite handle           | 080130-31 | 080130-31 | 080130-31 |
| 23 | Tornillos pata              | Legs screws               | 080142-22 | 080142-22 | 080142-22 |
| 24 | Pata                        | Legs                      | 080142-23 | 080142-23 | 080142-23 |
| 25 | Interruptor                 | Switch                    | 080400-10 | 080400-10 | 080400-10 |
| 26 | Control                     | Control                   | 080142-50 | 080142-50 | 080148-50 |
| 27 | Microinterruptor            | Microswitch               | 040060-60 | 040060-60 | 040060-60 |
| 30 | Cable                       | Cable                     | 000250    | 000212    | 000212    |
|    | Cable (América)             | Cable (America)           | 000252    | 000214    | 000214    |
| 31 | Disyuntor horno             | Furnace circuit breaker   | 080142-46 | 080144-37 | 080148-31 |
| 32 | Termopar                    | Thermocouple              | 080130-04 | 080130-04 | 080130-04 |
| 33 | Regleta                     | Terminal board            | 080130-06 | 080130-06 | 080130-06 |
| 34 | Base empotrable             | Base schuko               | 080084-01 | 080084-01 | 080084-01 |
| 35 | Triac                       | Triac                     | 080060-24 | 080060-24 | 080060-24 |
| 36 | Disyuntor extractor         | Extractor circuit breaker | 080142-45 | 080142-45 | 080142-45 |

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The R-080142, R-080144 and R-080148 are three practical furnaces especially designed to heat casting rings and crucibles for the dental technicians. They have a control panel on which a pre-set sequence of temperatures and times can be programmed by the user. Furthermore, up to 20 different programmes can be stored. This simplifies programming, as previous programmes can be recovered when required. Both furnaces work in the following 7 steps:

- Waiting Time A0: a time when the unit is at room temperature and is not running.
- 1st Ramp Time A1: a time when the temperature inside the furnace rises gradually to a level pre-set by the user.
- 1st Plateau Time A2: a time when the temperature inside the furnace remains steady at the level reached at the end of the previous step.
- 2nd Ramp Time A3: a time when the temperature inside the furnace rises gradually again, to a level pre-set by the user.
- 2nd Plateau Time A4: a time when the temperature inside the furnace remains steady at the level reached at the end of the previous step. When this time is up, the unit emits a warning buzz.
- 3rd Ramp Time A5: a time when the temperature inside the furnace rises gradually again to a level pre-set by the user.
- 3rd Plateau Time A6: a time when the temperature inside the furnace remains steady at the level reached in the previous step.

At the end of 2nd Plateau Time A4 the unit emits a warning buzz. Another buzz is emitted at the end of 3rd Plateau Time A6 and the furnace remains at this level until it is disconnected. Pilot lamp 4 also comes on.

**Heating ramps of the programmes can be edited either in grades/minute or just in time. In order to choose editing mode press key (10) for 8 seconds. The display on the left will show “-6-” and in the display on the right you enter the value from 0 (indicating minutes) or 1 (indicating °C/min). In order to save configuration, press key (10). Whenever this parameter is modified, all programmes stored in memory will erase.**

The time for each step can be adjusted from 0 to 99 hours 59 minutes. It is also important to point out that if we chose a “0” time for a given step and then we press the decrease key again, the display will show “UU.UU”. This means that when this step ends, the

control will understand that the program has also finished.

**Heating ramps of each step can also be adjusted between 1 and 15 °C/min. If we choose ramp “1” for a step and then we press key decrease, the display will show “AA.AA”. This means that the furnace will work at highest power till it reaches temperature set.**

The final temperatures at each ramp can be adjusted from 0 to 1100 °C.

To enable the user to see the actual program status at all times, the furnace control panel has a synoptic display. Lamps (11) indicate which step is ongoing, while lamps (12) show the temperatures reached.

The furnace also has two digital displays: the left hand one (1) indicates the present temperature in the furnace. The right hand display shows the time remaining to the end of A4 (when at steps A0, A1, A2 or A3) or the time remaining to the end of A6 (when at steps A5 or A6). These displays have two functions. To switch between them press W/T (6). When this key is pressed the left-hand display (1) shows the theoretical temperature setting expected by the microprocessor, and the right-hand display shows the time remaining to the end of the step in progress. When programming key (14) is pressed the right-hand display (2) shows the number of the program enabled.

Finally, on the left-hand side of the machine there is a current pickup which is extremely useful for connecting up a suction system. We recommend using the Smoke Extractor R-080125. The extractor’s control output operates when steps “A1” and “A2” are being carried out, and furnace temperature is between 75 °C and 850 °C. When temperature exceeds 850 °C, the control output disactivates, and will not reactivate even if the temperature falls below 850 °C. When steps “A3” and “A4” are being carried out, and furnace temperature is between 75 °C and 650 °C, the control output activates. When temperature exceeds 650 °C, the control output disactivates, and will not reactivate even if the temperature falls below 650 °C.

## PROGRAMMING

It is extremely easy to program a heating sequence into the furnace hearth and then save it as one of the 20 programs which can be stored. Proceed as follows:

1. Press the programming key (14). Check that the display (2) shows a number between 1 and 20: this is the number of the present program.
  2. Use the Up/Down keys (8 & 9) to select the number under which the program to be entered is to be saved.
  3. Press the function key (10). This accesses programming mode. Check that the left-hand display reads "A0" - this indicates that this step can now be programmed. The lamp for step "A0" (waiting time) on the synoptic display should also come on.
  4. Use the Up/Down keys (8 & 9) to bring up the time desired for the first step on the right-hand display (2).
  5. Press the function key (10). Check that the left-hand display shows "A1" - this indicates that this step can now be programmed. The lamp for step "A1" (1st ramp time **or** °C/min) on the synoptic display should also come on.
  6. Use the Up/Down keys (8 & 9) to bring up the time **or** °C/min desired for the second step on the right-hand display (2).
  7. Press the function key (10) again. This time the left-hand display should show "U1", indicating that the final temperature for the first ramp can be programmed.
  8. Proceed in the same way to program the time **or** °C/min and temperature of each step in the program.
  9. The order in which the steps appear on display (1) when the furnace is in programming mode is as follows:
    - (A0) - waiting time.
    - (A1) - first ramp time **or** °C/min.
    - (U1) - first plateau temperature.
    - (A2) - first plateau time.
    - (A3) - second ramp time **or** °C/min.
    - (U2) - second plateau temperature.
    - (A4) - second plateau time.
    - (A5) - third ramp time **or** °C/min.
    - (U3) - third plateau temperature.
    - (A6) - third plateau time.
  10. When the program is complete, or at any time during programming, you may exit programming mode by pressing W/T (6), which acts as an Escape key. The program will automatically be saved in the memory under the program number allocated to it at the commencement of programming.
- If you wish to eliminate any of the ramps or plateaux in the program, simply allocate zero time and the same temperature as for the previous plateau: e.g. to eliminate the 2nd ramp and 2nd plateau, set a figure of 00.00 for steps A3 and A4, and the same temperature at U2 as at U1. **When working in**

**°C/min, temperature for "U2" should be the same as for "U1" and any heating speed.**

After exiting programming mode, you can run the program by holding down Stop/Start (5) for approximately 5 seconds. To select another program, press the programming key (14) and select the number of the program desired by means of the Up/Down keys.

## INSTALLATION

- Carefully unpack the unit. Inside the box there should be an SiO<sub>2</sub> hearth plate. Check also that the muffle lining the inside of the furnace has suffered no damage during transport.
- Set the furnace up on a smooth, firm, horizontal base away from any source of vibration. Allow enough space for ventilation. There should be no other objects less than 15 cm from the unit. Ensure especially that the ventilation tube is free and at least 20 cm from the closest object or wall.
- During normal furnace operation fumes may filter through the ventilation tube. This tube should therefore be placed below an extractor hood.
- Connect the unit to a 230 V, 50/60 Hz, earthed power outlet which protected by a thermal magnetic cut-out with a suitable amp rating (consult your installer).
- Before using the furnace heat it up to 300 °C and hold it there for one hour to evaporate off any moisture in the muffle.

## OPERATION

- Place the SiO<sub>2</sub> hearth-plate on the hearth, taking care not to damage the muffle.
- Switch on the On/off switch (13). Some of the controls on the panel should flash to indicate that the furnace is ready to operate.
- Put in the cylinders and crucibles to be heated. When the furnace is operating, always use gauntlets and tongs. Talleres Mestraitúa, S.L. recommends MESTRA tongs.
- The furnace is factory-supplied with a series of standard programs in its memory. However, it is advisable to program the furnace as indicated above to ensure that the programs stored meet your own requirements. These programs can be used and modified as required in future.
- Press the programming key (14). Display (2) should show the number of the present program. Use the Up/Down keys (8 & 9) to select the desired program number. After programming, hold down Stop/Start (5) for approximately 5 seconds, and the furnace will run the program selected. Lamp (7) also comes on to indicate that the program is running with no problems.

Lamp (3) also remains on whenever the resistor is working.

- In the normal course of the program the furnace can be stopped at any time by holding down Stop/Start (5) for approximately 5 seconds. You can also enter programming mode and modify the parameters of the program which is running if you wish.
- If there is a power cut the furnace retains in its memory the parameters of the step in progress, and when power is restored the program will continue from the point at which it stopped. If the switch (13) is turned off before the program is completed, the furnace takes this as a power cut, and will therefore use maximum power when it is re-connected to attempt to retrieve the temperature level before disconnection.
- Once program step A4 is completed the furnace emits a warning buzz. A further buzz is emitted at the end of A6. This buzzer can be disconnected by pressing any key on the control panel.
- At the end of the process, remove the cylinders and crucibles carefully. While the unit is running the outer housing may heat up. The cylinders and crucibles also heat up considerably, and so gauntlets and tongs should be used when handling them.
- To run another program, first stop the furnace by holding down Stop/Start (5) for 5 seconds, then selecting another program via the programming key (14) and the Up/Down keys (8 & 9). Finally, restart the furnace by pressing Stop/Start (5) for one second.
- Use the On/Off switch to disconnect the furnace from the mains (13).

### MAINTENANCE AND CLEANING

To ensure long, reliable operation of the unit, proceed as follows:

- Clean the outside of the furnace regularly with a cloth moistened in soapy water.
- Take care not to damage the muffle lining the inside walls of the furnace hearth when you put in or take out cylinders and crucibles.
- Always use the SiO<sub>2</sub> hearth plate to prevent filtration, which may damage the muffle.
- Prevent liquids or solids from contaminating or dirtying the muffle of the furnace.
- To make best use of the unit, TALLERES MESTRAITUA, S.L. recommends the use of **MESTRA** cylinders and crucibles.
- If the machine detects that the temperature is higher than 1150 °C, a special display will appear, and also a characteristic acoustic signal indicating a dangerous situation (Fig. 1).  
These characteristic symbols will appear alternatively in the two displays, and disappear

after 3 seconds. The display will then show the usual signal. This situation will repeat every 5 seconds.

This alarm appears also when the triac breaks.

- When this happens and the alarm signal appears, disconnect the main switch and call for technical assistance.
- The furnace has two circuit breakers located near the cable (J and K) to protect the furnace against electric overloads. To get the furnace working again, reset the circuit breaker.

### IMPORTANT

**If the temperature indicator display (1) flashes during normal furnace operation and reads 1200 °C, disconnect the unit from the mains and consult your authorised MESTRA dealer: the thermocouple is damaged. The thermocouple is subject to wear, and its durability depends on working temperature. This element is not covered by guarantee.**

### PRECAUTIONS

- ◆ Before connecting the unit ensure that the power outlet is for 230 V, 50/60 Hz, and that it is earthed. It must also be properly protected by a thermal magnetic cut-out (consult your installer).
- ◆ Do not allow children or unskilled personnel to handle the unit.
- ◆ Always keep the door closed while the furnace is operating.
- ◆ Take care when opening the door: the hot fumes pouring out can cause burns.
- ◆ Do not overfill the furnace or place objects in it which might fall when the door is opened.
- ◆ Do not put inflammable materials or liquids in the furnace, or substances which may give off toxic fumes.
- ◆ Place the furnace on a flat, firm, vibration-free surface close to an extractor hood. Ensure that the furnace is well away from any flammable materials.
- ◆ Place the unit so that its sides are at least 15 cm from any wall or other object. Take particular care to place the ventilation tube in a free area at least 20 cm from the closest object.
- ◆ Always use gloves and tongs to remove and put in cylinders and crucibles. Do not touch the metal housing of the furnace when it may be hot.
- ◆ Handle heated materials removed from the furnace with great care: their temperature may be as high as 1100 °C and may therefore cause accidents.
- ◆ Before using the furnace for the first time, hold it at 300 °C for 1 hour to eliminate any moisture in the muffle.
- ◆ Always use the furnace under an extractor hood.

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