R-080270 / R-080280 / R-080290 T1 - T2 - T3

| | | |
|----------|------------------------|----|
| ES..... | CHORROS DE ARENA..... | 6 |
| EN..... | SANDBLASTERS..... | 9 |
| FR..... | SABLEUSES..... | 12 |
| PO | JACTOS DE AREIA | 15 |
| II | SABBIATRICI | 18 |
| DE..... | SANDSTRAHLGERÄTE | 21 |

MESTRA[®]


TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tel. +34 944530388 - Fax +34 944711725


mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

R-080270 / R-080280 / R-080290

| | |
|---|--|
| A | Grapa Pestillo / Bolt clip / <i>Manette-fermoir</i> / Grampo fecho / Morsetto nottolino / <i>Verriegelungsklammer</i> |
| B | Llave pistola fija / Fixed gun valve / <i>Anneau déclencheur</i> / Chave pistola fixa / Chiave pistola fissa / <i>Fixierte Umlaufdüse</i> |
| C | Sufridera de goma / Rubber pad / <i>Protecteur en caoutchouc</i> / Placa de choque de borracha / Tassello di gomma / <i>Gummieinlage</i> |
| D | Guante de manipulación / Handling glove / <i>Gants de travail</i> / Luva de manipulação / Guanto di manipolazione / <i>Strahlhandschuh</i> |
| E | Pedal microchorreado (T2, T3) / Micro-blasting pedal (T2,T3) / <i>Crayon (T2, T3)</i> / Pedal micro jacto (T2, T3) / Pedale microsabbatura (T2, T3) / <i>Fußschalter Feinstrahlen (T2, T3)</i> |
| F | Tapón depósito (T2,T3) / Tank stopper (T2, T3) / <i>Bouchons des silos (T2,T3)</i> / Tampa depósito (T2,T3) / Tappo serbatoio (T2,T3) / <i>Deckel Strahlmitteltank (T2,T3)</i> |
| G | Pomo manorreductor (T2,T3) / Pressure reducer knob (T2,T3) / <i>Manoréducteur (T2,T3)</i> / Puxador regulador de pressão (T2,T3) / Pomo riduttore di pressione (T2,T3) / <i>Drehknopf Druckminderer (T2, T3)</i> |
| H | Manómetro (T2,T3) / Pressure gauge (T2,T3) / <i>Manomètre (T2,T3)</i> / Manómetro (T2,T3) / Manometro (T2,T3) / <i>Manometer (T2, T3)</i> |
| I | Tapones aspiración / Exhauster stoppers / <i>Bouchons pour l'aspiration</i> / Tampas aspiração / Tappi aspirazione / <i>Zuluftstopfen</i> |
| J | Racor de entrada (T2,T3) / Inlet connector (T2,T3) / <i>Raccord d'admission d'air (T2,T3)</i> / Junção de entrada (T2,T3) / Raccordo di mandata (T2,T3) / <i>Anschlussstutzen (T2, T3)</i> |

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The R-080270 (T1), R-080280 (T2) and R-080290 (T3) are three practical sandblasting booths specially designed for micro-blasting small workpieces for the dental prosthesis sector. To obtain optimum performance from these units, proceed as follows:

INSTALLATION

- Carefully unpack the unit. Inside the booth you will find the following accessories:
 - 6 / 8 mm compressor connection tube (T2, T3).
 - 4 / 6 mm compressor connection tube (T1).
 - Bag filter (T1, T2, T3).
 - Two protective coverings for the sandblaster glass (R-080270-50).
- Set up the unit on a flat, rigid surface well away from any sources of vibration or heat, in an environment which is not too humid. Too much humidity could lead aluminium oxide particles to compact in the booth, and clog the micro-blasting nozzles. Locate the unit at a height where the objects in the booth can be handled comfortably.
- Connect the tube provided with the unit to the air inlet (J) on the right hand side of the booth. Connect the other end to an air outlet at a pressure of no more than 5 bar.
- Connect the inside light to 230 V, 50/60 Hz with earth.
- Position the pedal (E) of the micro-blasting gun within reach of your foot, where it can be worked comfortably (T2 and T3 only).
- Cover the floor of the booth with aluminium oxide (approx. 3 kg) with particle diameter less than 250 microns. While spreading the aluminium oxide, take care to cover the area around the filter on the left hand side inside the booth.
- Remove the top threaded stopper (F) from the tank of the micro-blasting module, and fill the vessel with a grinding medium with particle diameter less than 250 microns, leaving at least 30 mm clearance at the top of the brass tube (T2 and T3 only).

Grinding medium heights of less than 35 microns are not recommended. Screw the tank stopper (F) back in place, ensuring that the seal is tight (T2 and T3 only).

PRESSURE REGULATION

The blowing pressure of the fixed micro-blasting gun is the same as the pressure at the air outlet (compressor). This pressure may be regulated between 2 and 5 bars to suit the nature and characteristics of the work in hand. **It should never exceed 5 bars.**

To adjust the pressure use the pressure regulator, which limits the output pressure from the compressed air unit.

The pressure of the moving gun micro-blasting module (T2 and T3) can also be regulated between 2 and 5 bars, bearing in mind that the upper limit depends on the blowing pressure of the air unit used. **In any event, it should never exceed 5 bars.**

To adjust the pressure in the micro-blasting unit, pull up the pressure regulator knob and turn it to the desired setting.

The micro-blasting module pressure is shown on the pressure gauge (H) built into the regulator. Once pressure has been selected, press the regulating knob back down to lock it.

EXHAUSTER ADJUSTMENT

The blasting cabin should always be used connected to a pressured aspiration unit, via the outlet tube at the back of the units. For proper exhausting, a slight vacuum pressure must be created inside the booth so that there is a smooth, continuous flow of dust particles towards the outlet tube.

The booth vacuum pressure depends on the blowing pressure used. It can be adjusted as follows:

- Remove the three exhaustor stoppers (I) on the top right-hand side of the booth.
- Connect the exhaustor to the outlet tube at the back of the booth and start it up.
- Start up the fixed micro-blasting gun by activating the valve (B).
- Place your thumb over the hole of one of the exhaustor stoppers. You will feel a slight suction effect which tends to stop up the hole.
- If the suction level is low, insert one of the stoppers and repeat this process.

6. Do this until air circulates smoothly through the holes into the booth.

The unit is supplied with a bag filter for exhausting. This should only be used sporadically (temporary installation, exhauster malfunction, etc.). In general, motorised exhausting should be used.

NOTE: make sure that the aspiration unit used can suck abrasive particles. We strongly recommend the MESTRA Eolo Plus R-080531 aspiration unit.

USING THE UNIT

To micro-blast workpieces correctly, proceed as follows:

1. Open the booth lid and put the workpiece inside. Check that the aluminium oxide level inside is sufficient.
2. Check that the level of grinding medium in the micro-blasting module vessels are within the limits for work (T2 and T3).
3. Close the booth lid and ensure that the bolt clip (A) is properly in place.
4. Use suitable protective goggles and a face-mask.
5. Start up the exhauster. If you need to see the workpiece clearly, switch on the inside light.
6. To use the fixed gun, move the slide on the valve (B). Move it forward to make the gun blow at maximum pressure, half-way along its stroke for half pressure, or fully back to cut off the air flow.
7. Using the rubber gloves (D) built into the booth, move the workpiece so that it is pointed properly towards the blasting jet.
8. If you wish to use the micro-blasting module (T2 and T3), pick up the moving gun in the booth and press the pedal. In the case of T3, first choose the module with the selection valve. When a module starts up, point the sand jet towards the rubber pad (C) on the rear wall of the unit. If the jet blasts against any other area, the booth walls will be subject to abrasion over time (T2 and T3).
9. When work is completed, move the slide on the fixed gun valve (B) fully back. Keep the exhauster working for a few seconds to remove all dust particles from inside the booth.
10. Turn off the exhauster and open the booth lid to remove the workpiece.

MAINTENANCE AND CLEANING

To ensure a long service life for the unit, proceed as follows:

- Always use high-quality grinding medium. Talleres Mestraitua, S.L. recommends MESTRA medium. The fixed micro-blasting gun (nozzle diameter 3 mm) can blow material with a particle diameter of less than 300 micron. The moving micro-blasting module gun of T2 has a nozzle diameter of 1.25 mm and can handle grinding medium up to 250 microns in diameter. The T3 has a gun with a diameter of 1.25 mm, and the other gun has a diameter of 1 mm and can handle grinding medium up to 150 microns in diameter. If grinding media with diameters greater than 200 microns are used in the fixed gun, remove the filter at the end of the gun tube in the booth.
- Use dry piston air compressors whenever possible. Oil particles in the blowing air can cause lumps of sand which will clog the gun nozzles.
- The grinding medium discarded by the moving micro-blasting gun (T2, T3) can be used in the fixed gun. However, never fill the micro-blasting module vessel with grinding medium from the booth floor, since large-diameter particles may clog the nozzle of the moving gun.
- To prevent abrasion and wear on the inner booth walls, always point the jet at the rubber pad on the rear wall when using the micro-blasting module (T2, T3).
- Stir up the grinding medium inside the booth periodically, and at the same time clean the filter hood on the left-hand side of the booth.
- For an overall cleaning of the unit, use a dry cloth or moistened with a proper product.

If a moving gun (T2, T3) becomes clogged, remove the nozzle-holder nut with an M-12 fixed wrench and remove the nozzle from inside. Unclog the nozzle with a pin or a thin wire, and refit it on the moving gun.

If the fixed gun becomes clogged, unscrew the nozzle-holder manually and remove the blowing nozzle from inside. Unclog with a pin or a thin wire, and refit.

The moving gun on the micro-blasting module T2 is factory-supplied with a 1.25 mm diameter nozzle. However, 1 mm and 0.75 mm nozzles can also be fitted: consult your **MESTRA** dealer. The moving guns on the micro-blasting modules T3 are factory-supplied with a 1.25 mm and 1 mm. diameter nozzle. The table below shows the grinding medium particle diameters for the different nozzle diameters:

| MAXIMUM PARTICLE DIAMETER FOR EACH NOZZLE | |
|--|--------------------------|
| Nozzle | Particle diameter |
| 1.25 mm | 250 microns |
| 1 mm | 150 microns |
| 0.75 mm | 75 microns |

PRECAUTIONS

- ◆ Never exceed a pressure of 5 bars.
- ◆ Always work with the booth lid closed, and use suitable protective goggles and a face-mask.
- ◆ Connect the booth outlet tube to an exhauster, and adjust the pressure inside the booth via the stoppers on the right-hand side.
- ◆ Do not use the bag filter supplied with the unit on a permanent basis: it is for short-term use only in unforeseen or temporary circumstances.
- ◆ If the rubber gloves become damaged or torn, replace them immediately.
- ◆ Do not allow children or unskilled personnel to handle the unit.
- ◆ Use only guaranteed high-quality grinding media (we recommend MESTRA media). Low quality media can be toxic if dust is inhaled.
- ◆ Connect the inside light to 230 V 50/60 Hz with earth.
- ◆ Ensure that the compressed air system used to power the unit is properly protected by at least one pneumatic discharge valve.
- ◆ If leaks are detected in the compressed air system, consult your installer immediately.

TECHNICAL CHARACTERISTICS

| | T1 | T2 | T3 |
|---------|-----------|-----------|-----------|
| Height: | 300 mm | 345 mm | 345 mm |
| Length: | 380 mm | 435mm | 435mm |
| Width: | 325 mm | 360 mm | 360 mm |
| Weight: | 6.5 kg | 8 kg | 9.5 kg |

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