
Ref. 080517

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The IRIS 2 EVOLUTION vacuum mixer is used to make investmets, plaster and similar mixtures. The machine has been designed as a simple and reliable tool for Dental Labs. For optimum performance and a long service life, please read the instructions for use carefully.

INSTALLING THE UNIT

1. Carefully unpack your IRIS vacuum mixer, and check for these fittings in the packaging:
 - Power cable.
 - Spare filter.
 - 550 cc. mixing bowl (D).
 - Kit to secure the mixer to the wall.
2. Choose a wall well away from any sources of heat, vibrations or moisture to set up the mixer unit. It must also be positioned on a firm, rigid, level surface.
3. Make two marks with a marker pen exactly 80 mm. apart at the chosen location. Ensure that the imaginary line linking the two marks is perfectly horizontal. A level may be used if required.
4. Drill out the holes marked by hand using a Ø 6 mm widia bur.
5. Insert the nylon studs included in the vacuum mixer kit in the drill holes. Then screw in the cheese-head shed bolts, allowing the heads to protrude around 12 mm.
6. Use the bolts to secure the IRIS vacuum mixer firmly to the wall.
7. Use the pen to mark the position of the fixing eye on the lower section of the unit. Remove the mixer and drill out the hole manually. Insert a nylon stud in the orifice and put the unit back on.
8. Use a wrench to put the bolt through the lower hole. Ensure the unit has been firmly secured.
9. Hook up the mixer's power cable to an adecuate power supply with earthing.

OPERATING THE MIXER

The mixer is now ready for use:

1. Use the main switch (E) to turn on the machine.
2. The upper display will show the spin speed in rpm and the lower display will show the working time in seconds.
3. You can change both of them by pressing the ▲/▼ keys.
4. To shift from a display to the other press the ☐ key.
5. Put the mixing ingredients into the bowl (D). **Pre-mix the products in the bowl manually with a spatula.**
6. Use gentle pressure to position the bowl and the mixer connection (F) with a good sealing.
12. Press the **START/STOP** Ⓞ key again for the unit to start a cicle.
13. Use the vacuum meter (B) to check the vacuum in the bowl. If there is no vacuum in the bowl, move it gently to improve its seating. Keep hold of the bowl until the meter needle moves and the cup keeps holded. You may now release the bowl.
14. During mixing operations, the needle must remain in the "VACUUM" section. When the pre-set time has elapsed, the unit will stop.
15. Use the vacuum drain valve (C) to remove the bowl from the unit. The mix is now ready to be used.
16. To stop the operation of the unit at any moment keep the **START/STOP** Ⓞ key pressed.

IMPORTANT

This machine has a "diaphragm" vacuum pump (maintenance-free). If the unit is started at very low room temperature, the membrane rubber may not be sufficiently elastic for best operation. Thus we recommend that the mixer be run for approximately 1 minute before starting operations in cold environments.

MODIFYING THE PROGRAMS

- The electronic controls can store in memory four programs. To store a program, use the **▲/▼** **↻** keys as in the previous section.
- After setting the data you want to store, press the **P** key as many times as required to select the P1, P2, P3 or P4 program.
- Then, press the **P** key again for 3 seconds. The data will be stored in memory.

USING THE PROGRAMS

- To use the data stored in any of the four programs, press the **P** key as many times as required to select the P1, P2, P3 or P4 program.
- The data stored in the program will show directly on the displays until you press the **START/STOP** **⊙** key.

MAINTENANCE AND CLEANING

For best reliability and a long service life, we recommend the following:

- Clean the mixing bowl carefully each time it is used. This will prevent dirt building up and seizure of the blade shaft. Make sure the bowl seat and joint are cleaned for a complete sealing.
- Keep the unit and the surroundings clean. Do not use solvents or any other inflammable products.
- The filter inside the bowl connection (F) must be replaced periodically. Unscrew the filterRef. silencer and replace it with a new one. You can use a 10 mm box spanner.

- We recommend using always original MESTRA spare parts for mixer maintenance.

PRECAUTIONARY MEASURES

- ◆ Before connecting up the unit, check the power supply is according to the label data and with earthing.
- ◆ Do not fit the bowl to the unit when the mixer motor is on. Make sure the blades-bowl set is properly assembled.

- ◆ Do not let children or unqualified staff handle the unit.
- ◆ Choose a solid well-leveled wall to place the mixer, and make sure that fixing bolts are firmly screwed into the wall.
- ◆ Wait for the motor to stop completely after finishing a cycle before removing the bowl.
- ◆ Disconnect the unit from the power supply during cleaning operations.

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