
Ref. 080533

ES.....	6	Manual original
EN.....	10	
FR.....	14	

MESTRA®
UNIDAD DE ASPIRACIÓN
ASPIRATION UNIT
UNITÉ D'ASPIRATION
EOLO INDUCTION
Ref. 080533

Nº	Referencia	Descripción	Description	Description
1	080530-04	Tubo de aspiración	Aspiration tube	Tube d'aspiration
2	080532-14	Microfiltro de motor	Motor micro filter	Micro-filtre moteur
3	080531-30	Filtro de bolsa	Bag filter	Sac d'aspiration
4	000232	Fusible 5A	Fuse 5A	Fusible 5A
5	080533-01	Motor de aspiración	Aspiration motor	Moteur d'aspiration
6	000250	Cable de conexión	Connection cable	Câble de connexion
7*	080530-26	Bifurcación dos puestos - Y	Two places Y-connector	Bifurcation d'aspiration - Y

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The Ref. 080533 EOLO INDUCTION is a practical aspiration unit specially designed to be used within the field of dental prosthesis. It allows the aspiration of aluminium oxide, silicon oxide, resin dust, plasters, and investments and, in general, any material used in the dental lab. The unit's most important features are the following:


- Three-stage filtering circuit, including a filter-bag, a cartridge-filter and a plate-filter in the clean chamber.
- Indicator of the filter's filling level. "Obstruction/Filter Full" alarm. Long operation alarm, which disconnects the unit after operating continuously for more than three hours.
- Electronic regulation of aspiration power.
- Possibility of "automatic" mode operation, synchronizing the machine's start up with 2 micromotors and a sand blast.
- Possibility of simultaneous operation in two working places.

In order to achieve the unit's optimal performance and a long working life, we strongly recommend to read carefully the following operating instructions.

INSTALLATION


1. Carefully unpack your EOLO PLUS aspiration unit. You will find a Ø 35 mm flexible connecting hose inside the packing.
2. The machine has been designed to operate under the workbench, but it can be placed anywhere else. Choose a firm surface to place it, far from heat sources, damp or vibrations.
3. Open the filter lid (D) and make sure the filter-bag is correctly placed. Close the lid again.
4. Connect the unit's aspiration hose to the aspiration nozzle (A). To do this, simply press it firmly. The other end of the hose can be connected to a micromotor, sand blast, dry trimmer or any other similar device. In case you want to use the aspiration unit in two working places simultaneously, you need a "T junction for two working places (ref. 080530-26)", supplied as an accessory.
5. Connect the power cord to a grounded 230 V, 50/60 Hz power source.
6. The unit is ready to be operated.

MAIN DISPLAY


USER CONFIGURATION

1. On the home screen, hold down the key .

2. Briefly press the key again to enter option 1.0 User.

3. Select the language by pressing the keys y .
4. Save your selection by pressing the key .


OPERATION

1. Activate the main switch (3).
2. Press the key to select *Automatic* or *Manual* mode.
 - **Automatic mode:** It is indicated by the green light on. In this case, the suction unit works synchronized with a micromotor, so that when you step on the pedal, the suction starts automatically and remains running until about five seconds have elapsed after the pedal is released.
 - **Manual mode:** In this case, the suction works continuously governed only by the key .


- The display will show the suction power (from 30 to 100%). Adjust the power acting on the keys  and .

CALIBRATION FOR THE DETECTION OF MICROMOTORS

- Connect the power cable of the micromotor to one of the rear power sockets, which are marked with numbers 1, 2, 3.
- Turn on the micromotor and set it to spin at 5000 rpm but don't start it yet.
- From the home screen, press and hold the key .


- Press the key  to select the menu 2.0 Devices and briefly press the key  again.

Select device 1, 2 or 3 with the key . Press the key .

- The “HENA” parameter allows enabling (YES) or disabling (NOT) the automatic detection in one of the sockets. It is recommended to disable the sockets that are not going to be used, to eliminate possible interference.

- To calibrate a micromotor,
 - access the parameter “C.Nmot” with the key ,
 - select the value YES with the key 
 - and press the key .
- The display will show the message “ACTIVATE”. Now run the micromotor at 5000 rpm.
- After a few seconds, the display will show the message “RELEASE” while emitting an acoustic signal. Now stop the micromotor.

- The aspiration will start for a few seconds and one of the following messages will appear:

“CAL OK”: Calibration is successful.

“ERROR!”: The calibrated utensil will probably not be correctly detected automatically. Electrical noise can be excessive due to a bad ground or because the level to be calibrated is too low.

- Select another micromotor to calibrate or abandon the process by pressing the key .

FILTER FILL DETECTOR

An indicator of the filter fill level is displayed at all times on the main screen.

PROTECTIONS

If an obstruction occurs in the suction tube, or if the filter is full, the suction will stop.

IMPORTANT

The full or clogged filter detection system has been conceived as a protection for the motor of the suction unit against possible eventualities. In no case should you wait for the device to act before proceeding to change the filter. Periodically check the fill level of the filter-bag. (See section “Replacement of filters” of this manual).

The suction unit also has a motor protection system that disconnects it if it works continuously for more than 2 hours. This device is intended to prevent damage to the motor in the event that, due to forgetfulness, the unit remains in operation for too long a time. The equipment will only allow vacuuming if we keep the machine on for 15 minutes.



The rear power sockets allow devices to be connected up to a maximum combined power of 400 W. If this power is exceeded, the thermal switch (5) will act, disconnecting the micromotor. If this circumstance occurs, wait about 10 minutes before resetting the thermal switch by pressing the button.

REPLACEMENT OF FILTERS

The EOLO PLUS suction unit has a filter circuit in three stages: the first consists of a filter-bag; a second stage with a cartridge-type absolute filter; and a third stage assigned to a sheet filter located in the clean chamber. This last filter has an unlimited duration, so it should only be changed occasionally, and its replacement is reserved for the Technical Service. To clean and replace the filters we recommend:

Replacing the cartridge filter:

1. Open the upper cover (D) of the suction unit by acting on the latch (C), and remove it, together with the filter-bag.
2. Loosen the screw holding the filter holder frame to access the cartridge filter. Remove it carefully so that no dust particles fall on the inlet of the vacuum motor

3. Place the new cartridge filter, and proceed to assemble the rest of the components.

The cartridge filter cannot be recycled, so it always requires its replacement. **It is recommended to proceed to the same every 1500 or 2000 hours of operation of the machine.**

Replacing the filter-bag:

1. Open the upper cover (D) of the suction unit by acting on the latch (C). Remove the filter-bag together with the cover of the suction chamber of the machine. Carry out this operation very carefully to prevent dust particles from falling into the suction chamber of the appliance.
2. Present the new filter inside the suction chamber. Make sure the filter fits over the squeegee ring (A), and the filter base rests on the bottom tray of the vac chamber.
3. Close the top cover of the machine again.

MAINTENANCE AND CLEANING

To achieve optimal performance and a long machine life, we recommend:

- Cleaning the outer casing of the appliance should always be done with a cloth dampened in soapy water. Do not use solvents or any other flammable product.
- Check from time to time the filling level of the filter-bag. **Do not wait for the Obstruction / Filter Full detection device to activate.**
- Periodically proceed to replace the filters. The cartridge filter has an estimated life of about 1500/2000 operating hours. Always use original MESTRA spare parts.
- Carry out the filter change operations very carefully, trying to dirty the inside of the suction chamber as little as possible. Make sure that the filter-bag is well seated in the suction nozzle (A) and that it rests on the lower tray of the suction chamber.
- In the event of an overload in a Schuko socket (4) the thermal switch (5) will actuate. If this circumstance occurs, wait a few seconds and reset it by pressing its button.
- Check that the suction tubes in your laboratory installation do not have leaks due to cracks, breaks, or poor connections.
- It is not convenient for the suction motor for the machine to work continuously for long periods of time (over 2 hours). In any case, the machine has a protection device against this eventuality.

PRECAUTIONS

- Before connecting the appliance, make sure that it is a 230 V, 50/60 Hz socket with earth.
- Do not connect micromotors or other devices with a consumption greater than 400 W to a Schuko socket.
- Do not allow children or unqualified personnel to handle the appliance.
- Choose a rigid and well-leveled surface for the location of the suction unit.
- Do not open the top cover of the suction unit when the machine is running.
- During the tasks of cleaning the machine and replacing filters, brushes or motor, disconnect it from the network.

TECHNICAL SPECIFICATIONS

Voltage:	230 V, 50/60 Hz
Power:	750 W
Height:	585 mm
Width:	200 mm
Depth:	485 mm
Weight:	17 kg

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