
Ref. 080533

ES	7	Manual original
EN	13	
FR	19	



MESTRA®
UNIDAD DE ASPIRACIÓN
ASPIRATION UNIT
UNITÉ D'ASPIRATION
EOLO INDUCTION
Ref. 080533

Nº	Referencia	Descripción	Description	Description
1	080530-04	Tubo de aspiración	Aspiration tube	Tube d'aspiration
2	080532-14	Microfiltro de motor	Motor micro filter	Micro-filtre moteur
3	080531-30	Filtro de bolsa	Bag filter	Sac d'aspiration
4	000232	Fusible 5A	Fuse 5A	Fusible 5A
5	080533-01	Motor de aspiración	Aspiration motor	Moteur d'aspiration
6	000250	Cable de conexión	Connection cable	Câble de connexion
7*	080530-26	Bifurcación dos puestos - Y (accesorio opcional)	Two places Y-connector (optional fitting)	Bifurcation d'aspiration - Y (accessoire en option)

- | | |
|---|---|
| A | Boquilla de aspiración / Suction hole / <i>Entrée d'aspiration</i> |
| B | Asa / Handle / <i>Poignée</i> |
| C | Cierre / Lock / <i>Fermeture</i> |
| D | Tapa del filtro / Filter lid / <i>Couvercle du filtre</i> |
| E | Ruedas / Wheels / <i>Roulettes</i> |

1. Display
2. Teclado / **Keyboard** / *Clavier*
3. Interruptor general / **Main switch** / *Interrupteur général*
4. Tomas de corriente / **Power sockets** / *Prises de courant*
5. Interruptor térmico / **Thermal switch** / *Interrupteur thermique*

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The Ref. 080533 EOLO INDUCTION is a practical aspiration unit specially designed to be used within the field of dental prosthesis. It allows the aspiration of aluminium oxide, silicon oxide, resin dust, plasters, and investments and, in general, any material used in the dental lab. The unit's most important features are the following:


- Three-stage filtering circuit, including a filter-bag, a cartridge-filter and a plate-filter in the clean chamber.
- Indicator of the filter's filling level. "Obstruction/Filter Full" alarm. Long operation alarm, which disconnects the unit after operating continuously for more than three hours.
- Electronic regulation of aspiration power.
- Possibility of "automatic" mode operation, synchronizing the machine's start up with 2 micromotors and a sand blast.
- Possibility of simultaneous operation in two working places.

In order to achieve the unit's optimal performance and a long working life, we strongly recommend to read carefully the following operating instructions.

INSTALLATION

1. Carefully unpack your EOLO INDUCTION aspiration unit. You will find a Ø 35 mm flexible connecting hose inside the packing.
2. The machine has been designed to operate under the workbench, but it can be placed anywhere else. Choose a firm surface to place it, far from heat sources, damp or vibrations.
3. Open the filter lid (D) and make sure the filter-bag is correctly placed. Close the lid again.
4. Connect the unit's aspiration hose to the aspiration nozzle (A). To do this, simply press it firmly. The other end of the hose can be connected to a micromotor, sand blast, dry trimmer or any other similar device. In case you want to use the aspiration unit in two working places simultaneously, you need a "T junction for two working places (ref. 080530-26)", supplied as an accessory.
5. If you want to synchronize the suction with micromotors or other devices, connect them to the rear power sockets (maximum 400 W).
6. Connect the power cord to a grounded 230 V, 50/60 Hz power source.
7. The unit is ready to be operated.

MAIN DISPLAY

USER CONFIGURATION

1. On the home screen, hold down the key .

CUSTOMER CONFIG.
 1.0 User
 2.0 Devices
2. Briefly press the key again to enter option 1.0 User.

3.0 User
 LINGUA = English
 HRMOTOR = 0h
 HTMOTOR = 0h
3. Select the language by pressing the keys y .
4. Save your selection by pressing the key .


OPERATION

1. Activate the main switch (3).
2. Press the key to select *Automatic* or *Manual* mode.
 - **Automatic mode:** It is indicated by the green light on. In this case, the suction unit works synchronized with a micromotor, so that when you step on the pedal, the suction starts automatically and remains running until about five seconds have elapsed after the pedal is released.
 - **Manual mode:** In this case, the suction works continuously governed only by the key .
3. The display will show the suction power (from 30 to 100%). Adjust the power acting on the keys and .


SYNCHRONIZATION

If the automatic mode does not work properly, follow this calibration procedure:

1. Connect the power cable of the micromotor to one of the rear power sockets, which are marked with numbers 1, 2, 3.
2. Turn on the micromotor and set it to spin at 10 000 rpm but don't start it yet.
3. From the home screen, press and hold the key .

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
Select device 1, 2 or 3 with the key . Press the key .


5. **HBIL**: Enable device.
 - **YES**: Automatic detection enabled.
 - **NO**: Automatic detection disabled.

It is recommended to disable sockets that are not going to be used, to eliminate possible interference.


6. **GAIN**: Configurable sensitivity.
 - **LO**: Low sensitivity. It is recommended when connecting a utensil for the first time.
 - **HI**: High sensitivity. Use only if satisfactory results are not achieved with the **LO** value.

When changing the sensitivity, the device will need 10 seconds to reconfigure. The **GAIN ADJ.** message will flash.


7. To calibrate a micromotor,
 - access the parameter **CNmota** with the key ,
 - select the value **YES** with the key 
 - and press the key .

8. The display will show the message “**ACTIVATE**”. **Immediately** run the micromotor at 10 000 rpm.


9. After a few seconds, the display will show the message “**RELEASE**” while emitting an acoustic signal. Now stop the micromotor. **Immediately stop** the micromotor.

10. The aspiration will start for a few seconds and one of the following messages will appear:

“**CAL OK**”: Calibration is successful.

“**ERROR!**”: The calibrated utensil will probably not be correctly detected automatically. Electrical noise can be excessive due to a bad ground or because the level to be calibrated is too low.

11. If the device detection after calibration is not satisfactory, it is recommended:

- Repeat the process several times.
- Increase the speed of the device in steps of 2000 rpm and repeat the calibration.
- Set the sensitivity **GAIN = HI**, set the device to 10 000 rpm and repeat the calibration.

12. Select another micromotor to calibrate or abandon the process by pressing the key .

13. Finally, for exceptional cases, the calibration values obtained could be adjusted manually:

- Set **V_{mot}** = 0.
- Adjust **VN_{mot}** based on the behavior of the machine:
 - If aspiration is activated unexpectedly, increase **VN_{mot}** in steps of 10 at a time.
 - If the aspiration does not detect the utensil, decrease **VN_{mot}** in steps of 10 at a time.

It is recommended to periodically recalibrate the utensils because their operation may be affected by changes in environmental and electrical conditions.

Electrical noise may be excessive due to poor ground connection.

FILTER DETECTION

The main screen continuously shows an indicator of the filter fill level.

If an error appears for no reason, scan the code to access the calibration tutorial

PROTECTIONS

If an obstruction occurs in the suction tube, or if the filter is full, the suction will stop.

IMPORTANT

The full or clogged filter detection system has been conceived as a protection for the motor of the suction unit against possible eventualities. In no case should you wait for the device to act before proceeding to change the filter. Periodically check the fill level of the filter-bag. (See section “Replacement of filters” of this manual).

The suction unit also has a motor protection system that disconnects it if it works continuously for more than 2 hours. This device is intended to prevent damage to the motor in the event that, due to forgetfulness, the unit remains in operation for too long a time. The equipment will only allow vacuuming if we keep the machine on for 15 minutes.

The rear power sockets allow devices to be connected up to a maximum combined power of 400 W. If this power is exceeded, the thermal switch (5) will act, disconnecting the micromotor. If this circumstance occurs, wait about 10 minutes before resetting the thermal switch by pressing the button.

REPLACEMENT OF FILTERS



The EOLO INDUCTION suction unit has a filter circuit in three stages: the first consists of a filter-bag; a second stage with a cartridge-type absolute filter; and a third stage assigned to a sheet filter located in the clean chamber. This last filter has an unlimited duration, so it should only be changed occasionally, and its replacement is reserved for the Technical Service. To clean and replace the filters we recommend:

Replacing the filter-bag:

1. Open the upper cover (D) of the suction unit by acting on the latch (C). Remove the filter-bag together with the cover of the suction chamber of the machine. Carry out this operation very carefully to prevent dust particles from falling into the suction chamber of the appliance.
2. Present the new filter inside the suction chamber. Make sure the filter fits over the squeegee ring (A), and the filter base rests on the bottom tray of the vac chamber.
3. Close the top cover of the machine again.

Replacing the cartridge filter:

1. Open the upper cover (D) of the suction unit by acting on the latch (C), and remove it, together with the filter-bag.
2. Loosen the screw holding the filter holder frame to access the cartridge filter. Remove it carefully so that no dust particles fall on the inlet of the vacuum motor

3. Place the new cartridge filter, and proceed to assemble the rest of the components.

The cartridge filter cannot be recycled, so it always requires its replacement. **It is recommended to proceed to the same every 1500 or 2000 hours of operation of the machine.**

MAINTENANCE AND CLEANING

To achieve optimal performance and a long machine life, we recommend:

- Cleaning the outer casing of the appliance should always be done with a cloth dampened in soapy water. Do not use solvents or any other flammable product.
- Check from time to time the filling level of the filter-bag. **Do not wait for the Obstruction / Filter Full detection device to activate.**
- Periodically proceed to replace the filters. The cartridge filter has an estimated life of about 1500/2000 operating hours. Always use original MESTRA spare parts.
- Carry out the filter change operations very carefully, trying to dirty the inside of the suction chamber as little as possible. Make sure that the filter-bag is well seated in the suction nozzle (A) and that it rests on the lower tray of the suction chamber.
- In the event of an overload in a Schuko socket (4) the thermal switch (5) will actuate. If this circumstance occurs, wait a few seconds and reset it by pressing its button.
- Check that the suction tubes in your laboratory installation do not have leaks due to cracks, breaks, or poor connections.
- It is not convenient for the suction motor for the machine to work continuously for long periods of time (over 2 hours). In any case, the machine has a protection device against this eventuality.

PRECAUTIONS

- Before connecting the appliance, make sure that it is a 230 V, 50/60 Hz socket with earth.
- Do not connect micromotors or other devices with a consumption greater than 400 W to a Schuko socket.
- Do not allow children or unqualified personnel to handle the appliance.
- Choose a rigid and well-leveled surface for the location of the suction unit.
- Do not open the top cover of the suction unit when the machine is running.
- During the tasks of cleaning the machine and replacing filters, brushes or motor, disconnect it from the network.

TECHNICAL SPECIFICATIONS

Voltage:	230 V, 50/60 Hz
Power:	750 W
Height:	585 mm
Width:	200 mm
Depth:	485 mm
Weight:	17 kg

WARRANTY CONDITIONS

Requirements for the validity of the warranty

- All **MESTRA** machines are guaranteed against any manufacturing defect for one year from delivery.
- Delivery is considered to have been made from the moment the product is made available to the customer by the carrier.
- To service any warranty, it will be an essential requirement to provide the purchase invoice.

Warranty exclusions

The device has been designed and built for use by professionals in a dental clinic, dental laboratory or jewelry workshop. No other application is foreseen by the manufacturer, who consequently rejects all responsibility for damage of any kind caused by improper use of the device.

Maintenance work specific to each device, regardless of its frequency, is not covered by the warranty.

The following situations are excluded from the guarantee:

- Parts subject to natural wear (spray nozzles, tubes, lamps, filters, thermocouples, batteries, etc.).
- Breakdowns caused by physical breakage such as casings, plastics, paints, enamels, glass and similar.
- Breakdowns caused by fire, flood, weather agents, impact, crushing, inadequate ventilation or application of inadequate voltage or energy.
- Breakdowns caused by misuse by not following the safety, installation, operation, handling, cleaning, maintenance or installation instructions that are attached with each device.
- Breakdowns caused by the use of inappropriate or non-original accessories or spare parts.
- Appliances that do not have an identified factory number or that have been altered or deleted, or that do not match the number expressed on the invoice.
- Devices that have been modified, repaired or manipulated by personnel not authorized by **MESTRA**.
- Natural wear and tear from use does not represent any defect.
- **MESTRA** is not responsible for any damage that a device may cause to external elements.
- **MESTRA** reserves the right to modify the design, construction and appearance of its products without prior notice.
- All demands that exceed the limits of this guarantee are excluded, especially demands for compensation due to loss of profit, compensation for use, as well as indirect damages, as long as they are not mandatory liability according to law.

Damage caused by transportation

- In cases where transportation is the responsibility of **MESTRA**, the customer must examine the product upon receipt. If any damage has occurred due to transport, the client must immediately note it on the carrier's delivery note and notify **MESTRA**. Otherwise, **MESTRA** will not be able to claim from the carrier.

Repair conditions

MESTRA will repair or replace the machine as it deems appropriate, under the conditions stipulated below.

- Due to this guarantee, our company is only responsible for repairing or replacing the product or part recognized as defective and returned following the agreed procedure.
- During the period of repairs or changes no equivalent equipment will be supplied.
- Defective parts or components removed will become the property of **MESTRA**.

Repairs outside the mestra facilities

- **MESTRA** will not carry out any repairs outside its workshop and is exempt from the transportation costs this may cause.
- Repairs made at the customer's home, whether under warranty or not, are entrusted to the technical service of the distributor that made the sale or to its contracted technicians.

Send conditions

- Any product return invoking the guarantee must be previously notified and accepted. To this end, the customer will contact our after-sales service.
- The defective product will preferably be sent in its original packaging or, failing that, perfectly packaged.
- The costs and risks derived from transporting the product are the responsibility of the client and these will always be carried out through a courier company so that the merchandise is insured against possible damage during the journey.
- The equipment to be repaired must be accompanied by a note indicating the defect observed, name, address and telephone number of the user. To speed up the process, it is necessary to add an exact description of the defect.

If the client requests the replacement of a component under warranty, the amount will initially be charged. The customer will send the defective component to the factory for review and payment of the amount.

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