Ref. 100291

ES.....	PULIDOR ELECTROLÍTICO	5
EN.....	ELECTROLYTIC POLISHER.....	7
FR.....	POLISSEUR ELECTROLYTIQUE	9
DE.....	ELEKTROLYTISCHES GLÄNZGERÄT	11

MESTRA®


TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60


Tel. +34 944530388 - Fax +34 944711725

info@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA


A	Tapa / Lid / <i>Couvercle</i> / Deckel
B	Conector cátodo / Cathode fitting / <i>Connecteur cathode</i> / Kathodenstecker
C	Conector ánodo / Anode fitting / <i>Connecteur anode</i> / Anodenstecker
D	Panel de control / Cathode fitting / <i>Panneau de contrôle</i> / Bedienfeld
E	Cuba / Tank / <i>Cuve</i> / Behälter
F	Placas cátodo / Cathode plates / <i>Plaques cathode</i> / Kathodenplatten
G	Gancho / Hook / <i>Crochet</i> / Haken
H	Pinza / Clip / <i>Pince</i> / Klemme
I	ON-OFF
J	Marcha-paro / Run-stop / <i>Marche-Arrêt</i> / Start-Stopp
K	Cambio display / Display change / <i>Changement display</i> / Anzeigewechsel
L	Display
M	Memoria / Memory / <i>Mémoire</i> / Speicher
N	Subir-bajar / Up-down / <i>Monter-Descendre</i> / Auf-Ab

A: **Ánodo (+)**, superficie del metal a pulir.
Ánodo (+), surface of the metal.
Anode (+), surface du métal à polir.
 Anode (+), Fläche des zu polierenden Metalls.

m: **Átomo de metal / Metal atom / Atome de metal / Metallatom**

i: **Iones portadores / Carrier ions / Ions porteurs / Trägerionen**

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
Current displacement.

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Stromfluss.

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The R-100291 is a practical electrolytic polisher especially designed for use in jewellery work and dental mechanics. So that it may work properly and perform polishing of high quality. We recommend that you carefully read the following instructions.

OPERATION

Any metal studied under the microscope exhibits a more or less wrinkled surface, with projections of every sort. If we put the same metal into an electrolytic bath (fig. 1) and connect it to a component with positive electric voltage (anode), there is a curious effect: the electrolyte withdraws, giving rise to charged molecules, called ions, that tear atoms off the surface of the metal to combine with them. The resulting particles have an electric charge, for which reason, when they are submitted to the difference in potential between cathode and anode (fig. 5), there is a smooth displacement toward the negative pole (cathode). The metal particles are deposited on the cathode (fig. 3), and there is again a charged ion that starts the cycle off anew. The process comes to an end when the electric flow is cut.

The parts of the metal surface on which there appear the most pronounced projections are those that have a large surface contact with the electrolyte, hence those that lose most metal atoms during polishing. In consequence, when polishing is over the surface of the metal has much less of an escarpment appearance (fig. 4), with the rough edges smoothed down, so that light can be reflected more homogeneously. The metal is therefore brighter. The excess of metal removed from the anode (+) ends up deposited on the surface of the cathode.

INSTALLATION

- Choose a flat and firm surface far from heat sources or vibrations.
- Connect the machine to 230 V, 50/60 Hz earthed power source.

USE OF THE UNIT

1. Fill the tank (E) (80 % approx.) with electrolyte. Be careful to ensure that the electrolyte is of high quality, for in great

measure the efficiency of the polishing will depend on this factor.

2. BLACK corresponds to the cathode and should be connected to the annular electrode at the bottom of the tank (E). RED corresponds to the anode and should be connected to the clip (H) on top of the tub (4).
3. Connect the metal parts that you wish to polish.
4. Take care to ensure that the electrical contact with the electrode is good. The piece must be totally submerged in the bath.
5. Close the lid (A).
6. Press ON-OFF (I) to start the unit. Display will light.
7. Two left digits will show the intensity (amperes).
8. Two right digits will show the working time (minutes).
9. Select required intensity and time. Press A to alternate between both values and ▲▼ to modify.
10. It is probable that during polishing the electrolyte in the tub will heat up. This is perfectly normal. We recommend not using continuously for more than 30 minutes and rest for 10 minutes before new use.
11. Also, it is possible that, by accident, pieces for polishing that are positioned on the anode hooks touch the annular cathode of the tank. When this happens, there is a short circuit and display will show the message "ERR". In this case, remove the pieces causing the short circuit and wait three minutes before resetting the circuit breaker.

MAINTENANCE AND CLEANING

To ensure that the unit functions at its best and has a long life, we recommend the following:

- Always use electrolytes of recognized quality. This is especially important if you want optimum polishing.
- Regularly renew the electrolyte according to the maker's instruction.

- Periodically empty the tank and remove its electrodes in order to clean thoroughly. Use a scourer and soap.
- Do not let the electrolyte come in contact with the control panel.
- Do not overload the unit.
- Avoid short circuits between the electrodes of the unit.

TECHNICAL DATA

Power supply:	230 V, 50/60 Hz
Power:	100 W
Timer:	0-30 minutes
Intensity:	0.5 - 3 A
Tank capacity:	1 litro
Tank dimensions:	100 x 100 x 100 mm
Height x width x depth:	145 x 210 x 250 mm
Weight:	2.9 kg

MEMORY

You can save a program in the memory.

In order to use a program, just press the key (M2).

In order to save a program, first select the required intensity and time. Then keep pressing the key (M) for a time until a beep sounds.

PRECAUTIONS

- ◆ Make sure that voltage is correct (230 V, 50/60 Hz, and earthed).
- ◆ Do not use the unit in excessively humid places.
- ◆ Make sure that the unit is not used by children or by not specialized people.
- ◆ Antes de limpiar el aparato asegúrese que se encuentra desenchufado.
- ◆ Before cleaning it, make sure the unit is unplugged.
- ◆ Do not submerge the unit in any liquid for cleaning purposes.
- ◆ Do not put any liquid into the tub other than an electrolyte of high quality.
- ◆ While the unit is operating, the electrolyte will probably heat up. Certain fumes will result. Do not inhale them.
- ◆ Carefully read the instructions of the electrolyte manufacturer. Replace the electrolyte as often as specified. Handle it with care to avoid spillage. In addition, avoid prolonged skin contact.
- ◆ To dispose of the electrolyte when it has become degenerate, follow the instructions of the manufacturer in relation to respect for the environment.
- ◆ Do not overload the unit.
- ◆ Keep the top of the tub closed as long as possible, for in this way you will prevent uncontrolled spillage of electrolyte.

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