

Ref. 100295

ES.....	3	Manual original
EN.....	5	
FR.....	7	

A	Conector cátodo / Cathode fitting / <i>Connecteur cathode</i>
B	Conector ánodo / Anode fitting / <i>Connecteur anode</i>
C	Panel de control / Control panel / <i>Panneau de contrôle</i>
D	Pinza / Clip / <i>Pince</i>
E	Gancho / Hook / <i>Crochet</i>
F	Cuba / Tank / <i>Cuve</i>
G	Placas cátodo / Cathode plates / <i>Plaques cathode</i>
H	Tapa / Lid / <i>Couvercle</i>
I	ON-OFF
J	Marcha-paro / Run-stop / <i>Marche-Arrêt</i>
K	Piloto / Pilot / <i>Pilote</i>
L	Pantalla / Display / <i>Ecran</i>
M	Subir-bajar / Up-down / <i>Monter-descendre</i>

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The compact titanium anodizing unit is a simple device that allows you to process titanium dental elements quickly, easily and economically.

Anodizing is an electrochemical process that creates a protective oxide layer on the surface of the metal, in our case titanium. This layer – known as the anodic layer – is thicker than the natural oxide layer that forms on titanium when exposed to air, which offers a number of advantages:

Corrosion resistance:

The anodic layer protects the titanium from corrosion, even under adverse conditions (use in the mouth).

Wear resistance:

The anodic layer is harder than the base titanium, making it more resistant to wear by friction.

Better adhesion and surface finish:

Any layer of another material that is superimposed (e.g. zirconium) will be easier to adhere to the anodized base element (titanium).

Reduced release of metal ions:

The titanium oxide layer minimizes the release of metal ions into the physiological medium.

Appearance:

The anodic layer can be colored in different shades, which gives the titanium a more attractive appearance than the original brownish gray. This advantage is especially interesting in the case where the anodized element serves as a support for a translucent material, for example, zirconium.

To operate the machine, you only need to fill the tank with a solution of sodium bicarbonate (commonly sold in supermarkets) in water and then select the working voltage.

Depending on the voltage, different shades will be achieved. For example, at 60 V, an attractive golden shade is obtained, which is highly appreciated in the dental sector. The colouring process takes just 5 seconds.

INSTALLATION

- Choose a rigid, flat surface for the appliance, away from sources of vibration or heat.
- Connect the electrical plug to a 220-240 V, 50/60 Hz socket with an earth connection.


PREPARING THE TANK

- Pour into the tank (F) a solution of:
 - 600 ml of water.
 - At least 20 g of sodium bicarbonate (NaHCO₃).

Other electrolytes can be used. In this case, follow the manufacturer's instructions and take appropriate safety measures.

- Make sure to mix well.
- Place the tank (F) in its housing in the machine.
- Connect the metal plates (G) to the black terminal (A).
- Connect the clamp (D) to the red terminal (B).

INSTRUCTIONS FOR USE

1. Place the objects to be anodized on the holder (E).
2. Connect the holder (E) to the clamp (D), so that the objects are completely submerged.
3. Replace the cover (H).
4. Press the button  (I).
5. Use the   buttons (M) to select the voltage on the display (L).
(60 V for gold color)
6. Press the  button (J).
7. The indicator light (K) will light up red for 5 seconds and then turn green.

MAINTENANCE AND CLEANING

To ensure optimum operation of the appliance and a long life, we recommend:

- Always use electrolytes (liquids) of recognised quality. This factor is of vital importance to achieve optimum polishing.
- Periodically renew the electrolyte, following the manufacturer's instructions.
- Periodically empty the tank and remove its electrodes to carry out a thorough cleaning. Use a scouring pad and soap for this.
- Avoid the electrolyte coming into contact with the control panel.
- Do not overload the appliance.
- Avoid short circuits between the appliance's electrodes.

PRECAUTIONS

- Make sure that the mains connection voltage is adequate (230 V, 50/60 Hz with earth connection).
- Do not use the appliance in excessively humid rooms or locations.
- Prevent children or non-specialized personnel from handling the appliance.
- Before cleaning the appliance, make sure it is unplugged.
- Do not immerse the appliance in any liquid to clean it.
- Do not put any liquids other than quality electrolytes into the tank.
- During operation of the appliance, the electrolyte may heat up, producing certain vapours. Avoid inhaling them.
- Carefully read the electrolyte manufacturer's instructions and replace it as frequently as indicated. Handle it carefully to avoid spills. Also avoid prolonged contact with the skin.
- To dispose of the electrolyte once degenerated, follow the manufacturer's instructions (care of the environment).
- Do not overload the appliance.
- Keep the tank covered to avoid electrolyte splashes.
- Do not touch the electrodes with the machine on.
- Protect yourself adequately with gloves, glasses and a mask.

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