
Ref. 030422

ES.....	4	Manual original
EN.....	6	
FR.....	8	
DE.....	10	
PT.....	12	
IT.....	14	

N.º	Ref.	CONCEPTO	DESCRIPTION
1	030420-01	Manguera	Hose
3	-	Tapa	Lid
4	030420-04	Válvula de aire	Air valve
5	030420-05	Junta	Gasket
6	030420-53	Mango inferior	Lower handle
7	-	Cuerpo de la olla	Kettle body
8	030420-54	Asa inferior	Handle
9	030420-22	Tuerca de válvula - 1	Valve nut - 1
10	030420-23	Tuerca manómetro	Manometer nut
11	030420-11	Manómetro	Manometer
12	030420-56	Válvula apertura seguridad - 1	Open safe valve - 1
13	030420-13	Junta cobre manómetro	Manometer cooper gasket
14	030420-14	Junta cobre válvula	Valve cooper gasket
17	030420-52	Mango superior	Upper handle
22	030420-41	Válvula seguridad - 2	Safe valve - 2
23	030420-38	Junta cobre válvula	Valve cooper gasket
24	030420-39	Tuerca válvula - 2	Valve nut - 2

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The Ref. 030422 is a practical stainless steel pneumatic polymerising kettle, used mainly to polymerise dentures and associated items in dentistry. To get the best performance from the unit, use it as follows:

GENERAL CONDITIONS

Limitation of Liability

The manufacturer assumes no liability for secondary damages due to:

Failure to read these instructions.

Improper use.

Use of the product by unqualified personnel.

Technical modifications to the appliance.


Use of unauthorized spare parts.

Spare Parts

Counterfeit or defective spare parts may cause appliance deterioration and safety malfunctions. Use only original manufacturer parts and spare parts.

Unpacking and Testing

It is advisable to test the appliance and its accessories immediately after delivery, checking for any damage that may have occurred during transport.

- A Safety valve 1
- B Air inlet
- C Pressure gauge
- D Sliding bolt
- E Conector
- F Safety valve 2

USING THE UNIT

1. Pour 3 or 4 litres of hot water (50-60 °C) into the kettle.

Do not heat the polymerising water in the kettle, especially with pressurised air in it.

2. Put in the dentures to be polymerised and close the lid of the kettle, matching the thread on the top of the kettle body to that on the lid and turning clockwise until the handles meet. A slight click indicates that the kettle is hermetically sealed.

Do not force the lid open or close. If necessary, press the lid down with one hand while turning it with the other.

3. Connect the connector (E) to a compressor outlet.
4. Check that the top end of the valve on the handle is closed (turned clockwise).
5. Fit the blower nozzle of the connector (E) to the air inlet valve (B). During this operation, watch the needle of the pressure gauge (C). The optimum pressure for polymerising is 1.5-2 atmospheres. A green strip on the pressure gauge scale marks this area. Pressures of more than 2 atmospheres are indicated in red. Above 2 atmospheres the safety valve begins to lose air, and it opens completely at between 2.3 and 2.6 atmospheres.
6. Polymerisation takes between 5 and 20 minutes, depending on the type of resin used. After this time, turn the upper end of the safety valve located on the handle (A) to release the air inside the kettle.
7. Once all the air is removed, open the lid by sliding the red bolt (D) on the handle towards the centre of the kettle, and then turning the lid anticlockwise.
8. Remove the work pieces from the inside of the kettle.

PRECAUTIONS

- Read the instructions in full.
- Never exceed a pressure of 2.5 atmospheres.
- Ensure that the air outlet used is properly protected by at least one discharge valve. Consult your installer.
- Do not use the kettle for anything except polymerising (e.g. cooking).
- While cleaning the unit, never immerse the lid in liquid.
- Do not heat the polymerising water in the kettle, especially with pressurised air in it.
- Do not pour water into the kettle at more than 90 °C, as this could damage the pressure gauge and other components.

- Do not attempt to open the lid of the unit while there is still pressure inside. Release pressurised air by pulling the safety valve ring.
- If air leaks through the edges of the lid, replace the rubber seal as soon as possible.
- Do not attempt to override the kettle's safety measures.
- Do not allow children near the pressure cookers when in use.
- Do not place the pressure cooker in proximity to a source of heat.
- Extreme care should be exercised moving the unit when it is pressurised. The handles and grips must be used. Protect yourself where necessary.
- Never force the pressure cooker open. Do not open until you are sure the pressure has been completely released.


Do not force the lid open or close. If necessary, press the lid down with one hand while turning it with the other.

- Never use the pressure cooker without water.
- Keep these instructions.

MAINTENANCE AND CLEANING

- After each use, clean any residues of wax or resin off the inside of the kettle with an ordinary detergent for dishes. Clean the lid with a moist cloth.
- Replace the rubber seal in the lid regularly, using original MESTRA spares.
- While cleaning the lid, take care not to knock it, especially around the pressure gauge and safety valve.
- If the lid is very reluctant to turn or bad closing, the rubber seal can be cleaned and slightly lubricated with talcum powder.
- If the air inlet valve is leaking, it may just have come loose. Use a tire valve wrench or similar to tighten it.

TECHNICAL CHARACTERISTICS

Height:285 mm
 Width:420 mm
 Diameter:245 mm
 Temperature: 50 ~ 60 °C
 Capacity: 8 L
 Use pressure: 1.5 ~ 2 bar
 Safe valve:2.3 bar
 Material:Stainless Steel

WARRANTY

- The machine is guaranteed against any manufacturing defects for one year from the date of purchase.
- To comply with any warranty, MESTRA will require the purchase invoice.
- Any breakdowns caused by improper installation or misuse of the machine are excluded from the warranty.
- Furthermore, MESTRA will not perform any repairs outside its workshop and is exempt from any shipping costs incurred.
- Repairs performed at the customer's premises, whether or not under warranty, are entrusted to the technical service of the distributor who made the sale or to its contracted services.

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TALLERES MESTRAITUA S.L. is a company manufacturing tools and appliances for the dental prothesis sector, from the business address shown at the bottom of the page.

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I HEREBY STATE: That the **PNEUMATIC POLYMERISING KETTLE Ref. 030422**, with the following characteristics:

Alto / Height:	285 mm
Ancho / Width:	420 mm
Diámetro / Diameter:	245 mm
Dimensiones interiores / Inner dimensions:	↑210 mm x Ø220 mm
Temperatura / Temperature:	50 ~ 60 °C
Capacidad / Capacity:	8 L
Presión de uso / Use pressure:	1.5 ~ 2 bar
Material / Material:	Acero inoxidable / Stainless steel

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Meets the requirements set out in directive 2014/68/UE concerning pressure equipment.

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In witness whereof, the company issues this declaration of conformity in Sondika, on 30 July 2026.

Firma / Signed:

Ignacio Mestraitua,
Director Gerente / *Managing Director*

MESTRA®

TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tel. +34 944530388 - Fax +34 944711725

48150 SONDIKA - BILBAO - ESPAÑA

info@mestra.es - www.mestra.es