
Ref. 080229

ECLIPSE T

ES.....	5	Manual original
EN.....	7	
FR.....	9	

Nº	Ref.	CONCEPTO	DESCRIPTION
1	080233-02	Interruptor general	Main switch
2	080214-01B	Tapa cristal	Glass lid
3	080214-01C	Soporte cristal	Glass support
4	080214-01D	Soporte interruptor	Switch support
5	080214-04	Barra de LED	LED bar
6	080218-18	Difusor LED	LED diffuser
7	080218-19	Junta para difusor LED	Gasket for LED diffuser
8	080214-07	Rejilla ventilación pequeña	Small ventilation grill
9	080253-02	Bisagra	Hinge
10	080214-02	Cristal	Glass
11	080093-10	Burlete de goma	Rubber strip
12	080214-01E	Marco cristal	Glass frame
13	080220-10	Rejilla ventilación grande	Large ventilation grill
14	-	Cuerpo superior	Upper body
15	080229-02	Sufridera de goma	Rubber
16	080010-05	Pasa cables	Cable bushing
17	080614-20	Fuente de alimentación 24V	Power source 24 V
18	080216-02	Cubierta fuente de alimentación	Power source cover
19	080490-06	Base de conexión	Connection base
20	000312	Racor codo plástico Ø6 1/8	1/8 elbow Ø6 plastic connector
21	000337	Racord recto metálico H Ø6 1/8	Ø6 1/8 F straight metal connector
22	080270-30	Brida	Flange
23	080242-28	Manguito	Sleeve
24	080228-02	Cuerpo inferior	Lower body
25	080235-13	Codo plástico	Plastic elbow
26	080228-10	Tubo Ø40	Ø40 pipe
27	080270-31	Porta boquillas	Nozzle holder
28	080270-26	Boquilla 3 mm	3 mm nozzle
29	080270-32/C	Pieza intermedia	Intermediate part
30	080270-47	Inyector pistola	Gun injector
31	080270-48	Casquillo pistola	Gun bushing
32	030425-22	Racor codo metálico HH	FF Metal elbow connector
33	080262-27	Alargador	Extension
34	000336	Racor recto metálico Ø6 1/8	Metal straight fitting Ø6 1/8
35	080235-13	Tubo succión	Suction tube
36	080259-07	Válvula selectora	Selector valve
37	080259-15	Protector válvula selectora	Selector valve protector
38	080259-28	Arandela válvula selectora	Washer for selector valve
39	080259-17	Espiga válvula selectora	Dowel for selector valve
40	080614-04	Tapa tolva	Hopper lid
41	080228-06	Travesaño corto	Crossbar short
42	080228-07	Larguero derecho	Crossbar right
43	080228-08	Larguero izquierdo	Crossbar left
44	080228-04	Bandeja trasera	Rear tray
45	080228-05	Bandeja delantera	Front tray
46	100048-07	Contera pequeña	Small tip
47	080228-03	Chapa inferior	Lower sheet
48	080508-18	Pata	Leg
49	000101	Tubo Ø4x6 azul	Ø4x6 blue tube
50	000250	Cable 230 V	Cable 230 V
	000252	Cable 115 V	Cable 115 V
51	080220-12	Protector cristal (10 U.)	Glass protector (10 units.)

A	Interruptor luz / Light switch / <i>Interrupteur de la lumière</i>
B	Manguitos de tela / Sleeves / <i>Manchons</i>
C	Boquilla estática / Static nozzle / <i>Buse fixe</i>
D	Llave boquilla estática / Static nozzle switch / <i>Commutateur de buse fixe</i>
E	Tubo de entrada aspiración exterior / External suction inlet tube / <i>Tube d'entrée de l'aspiration</i>
F	Entrada aire a presión / Pressure air inlet / <i>Entrée d'air sous pression</i>

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The Sandblaster Eclipse is specially designed for precision-blasting small parts used in the dental prosthesis industry.

It has a fix nozzle that recycles the abrasive of the big tank under the cabin (particles up to 300 microns).

For optimum output and performance of this unit, we advise you to read the following instructions with great care.

Installing the unit

- Carefully unpack the unit.
- A rigid flat surface must be chosen for installation of the unit, at a safe distance from any sources of vibrations or heat, and the surrounding atmosphere must not contain excess moisture, since this could cause compacting of the emery particles inside the chamber, and lead to obstructions in the microblasting nozzles. The unit must also be positioned on a level which allows objects inside the chamber to be handled with ease.
- Hook up the air connection tube supplied with the unit to the connection at the rear of the unit. The other end must be fitted to a compressed air facility, with at least 6 bars of pressure and minimum flow rate 100 litres/minute.
- Hook up the electric plug to 115-230 V, 50/60 Hz, which must be earthed.
- Fill the tank under the cabin with aluminium oxide of a maximum diameter of 300 micron. To obtain a maximum autonomy, load the hopper up to a maximum height of 2 cm below the tray.

NOTE: the unit is optimized for 200-micron particles and getting an stable jet. With smaller particle diameters (100 microns) the jet may be more irregular.

Regulating the pressure

The blowing pressure of the blasting guns can be regulated depending on the nature and characteristics of the operation to be carried out. Use the pressure-adjustment knob (B) on the top of the unit – pull the knob upwards, and rotate it to the required pressure. The unit has a pressure gauge (C) for easy selection of operating pressure. When the pressure has been adjusted, lock out the regulator knob by pushing it down to the initial position.

Connection to the suction unit

The blasting chamber must be permanently connected to forced suction via an outgoing tube on the top left of the unit (G). Suction must be able to produce a slight depression inside the chamber for continuous gentle

circulation of the powder generated towards the outgoing tube.

Depending on the requirements and geometry of the external suction unit, it is possible to use the PVC elbow fitting supplied for specific needs.


Using the unit


1. Open the lid of the chamber, and put in the part you wish to blast. Parts may also be placed inside the chamber via sleeve.
2. Hold with your hand the pieces to sandblast inside the chamber.
3. Close the lid of the chamber, and ensure it is tightly shut (the glass should fit snug against the housing).
4. Suitable protection gear must be used: goggles, gauntlets and a mask.
5. Start up the forced suction facility.
6. The LED light in the chamber may be activated for a better view of the part, if required, using switch (G).
7. To use the fixed gun, move the slide on the valve (K). Move it upward to make the gun blow. Down position cuts the air flow off.
8. Move the workpiece so that it is pointed properly towards the blasting jet.

Maintenance and cleaning

To ensure a lengthy service life, we recommend the following:

- Always use approved quality emery - **Talleres Mestraitua, S.L.** recommends **MESTRA** emery. The mobile blasting guns with nozzle diameters 1.25 mm and 1 mm can respectively blow emery of maximum diameter 250 and 150 microns. The mobile blasting gun with nozzle diameter 0.75 mm can blow a maximum diameter of 75 microns (see the table).
- The gun (Ø 3 mm nozzle) can blast aluminium oxide up to particles of Ø 300 micron.
- In case of obstruction of the gun, dismount the nozzle holder. Clear the nozzle by blowing in the opposite direction to the air flow, or, if needed, by using a pin or a thin wire. Then set the gun again, making sure it is correctly centred.

If the nozzle came loose, glue it to its holder with cyanoacrylate. The nozzle should be level with the inner surface of the nozzle holder, to make sure that the distance from the injector to the nozzle is 4.5 mm and the Venturi effect takes place.

- It is preferable to use dry-piston air compressor units, since oil in the air could cause the sand to compact and obstruct the gun nozzles.
- Clean out the emery inside the chamber on a periodic basis, using the practical emptying plug inside.
- Use a soft cloth for general cleaning of the unit.

Precautionary measures

- Never exceed 6 bars of pressure. Optimum working pressure is between 2 and 5 bars.
- Operations should always be carried out with the lid of the chamber closed, with suitable protection. Masks and goggles should be worn.
- Connect the chamber's outgoing tube to a suction facility.
- The protection sleeves (H) must be replaced immediately following any tearing or damage.
- Do not allow children or unqualified staff to operate the unit.
- Use only guaranteed-quality materials as emery (we recommend **MESTRA** emery). Poor quality emery may be toxic when inhaled, and may also cause obstructions in the blow nozzles.
- Connect the power supply cord to 115-230 V, 50/60 Hz, earthed.
- Consult your fitter immediately if any leakages are observed in the compressed air circuit.
- Before carrying out any operation, first ensure that the micro blasting module covers are closed tightly.

Technical characteristics

Height:485 mm
 Width:360 mm
 Length:310 mm
 Weight:6,3 kg
 Working pressure: 2-6 bar
 Voltage: 115-230 V, 50/60 Hz
 Power consumption: 7 W
 Air connection:Ø6 mm

Warranty

- The machine is guaranteed against any manufacturing defects for **one year** from the date of purchase.
- To comply with any warranty, MESTRA will require the purchase invoice.
- Any breakdowns caused by improper installation or misuse of the machine are excluded from the warranty.
- Furthermore, MESTRA will not perform any repairs outside its workshop and is exempt from any shipping costs incurred.
- Repairs performed at the customer's premises, whether or not under warranty, are entrusted to the technical service of the distributor who made the sale or to its contracted services.

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