



**MESTRA®**  
 RECORTADORA DE MODELOS  
**MODEL TRIMMER**  
 TAILLE-PLATRE  
 R-080080-1-2-3-4

# MESTRA®

RECORTADORA DE MODELOS

MODEL TRIMMER

TAILLE-PLATRE

R-080080-1-2-3-4

Nº	Ref.	DESCRIPCION	DESCRIPTIPON
1*	080080-01	Manguera conexión agua	Water connection hose
2	000250	Cable 230V	Cable 230V
	000252	Cable 115V	Cable 115V
3	080080-03	Disco de carburo de silicio	Silicon carbide disc
	080081-01	Disco de diamante	Diamond disc
4	080251-09	Tornillo patas	Legs screws
5	080490-06	Base de conexión	Connection base
6	080080-40	Interruptor de seguridad	Safe switch
8	000236	Fusible 7A	Fuse 7A
	000237	Fusible 15A	Fuse 15A
	000238	Fusible 20A	Fuse 20A
9	080080-09	Condensador 230V	Capacitor 230V
	080080-09B	Condensador 115V	Capacitor 115V
10	080080-41	Motor con condensador 1500 rpm 230V	Motor with capacitor 1500 rpm 230V
	080080-41B	Motor con condensador 1500 rpm 115V	Motor with capacitor 1500 rpm 115V
	080082-01	Motor con condensador 3000 rpm 230V	Motor with capacitor 3000 rpm 230V
	080082-01B	Motor con condensador 3000 rpm 115V	Motor with capacitor 3000 rpm 115V
11	080080-42	Tornillos sujeción motor	Motor camping screws
12	030230-04	Tuercas del motor	Motor nuts
13	080080-44	Tapa trasera	Rear cover
14	080080-14	Junta tórica tapa trasera	Rear cover O-ring
15	080080-15	Porta retén	Seal holder
16	040250-03	Tornillos de porta retén	Seal holder screws
17	080080-45	Tapa cubre motor	Motor hood
18	080080-18	Espárragos tapa delantera	Front cover allen screws
19	080080-19	Tapa delantera	Front cover
20	080080-20	Tornillos laterales tapa delantera	Front cover lateral screws
21	080080-21	Pieza basculante	Tilting device
22	080080-22	Tuercas ciegas	Tilting device blind nuts
23*	080080-23	Electroválvula 230V	Electrovalve 230V
	080080-23B	Electroválvula 115V	Electrovalve 115V
24*	080080-24	Manguera grifo	Electro-tap hose
25	080080-46	Tapa de registro	Inspection cover
26	080080-26	Patatas de goma	Rubber legs
30*	080080-30	Racor de conexión grifo	Tap conector
31	080080-47	Tornillos unión de tapa trasera	Rear cover joining screws
32*	080080-32	Cuerpo de grifo (1500 rpm)	Tap body (1500 rpm)
	080082-02	Cuerpo de grifo (3000 rpm)	Tap body (3000 rpm)
33*	080080-33	Arandela de cobre	Copper washer
34*	080080-34	Mando de grifo	Tap control knob
35	080080-35	Tornillo baquelita	Bakelite screw
36	080080-56	Acoplamiento para motor	Motor coupling
37	080080-04	Reten para porta reten	Seal to seal holder
38	080080-54	Manguera desagüe	Sink hose
39	080080-49	Retén acoplamiento motor	Motor coupling seal
40	080080-53	Llave fija	Spanner
41	080080-55	Llave de tubo	Socket
42**	080084-01	Base empotrable	Base

\* No existe en 080084 / Does not exist in 080084

\*\* Sólo en 080084 / Only in 080084