


**MESTRA®**  
CAJA DE TRABAJO  
**WORKING BOX**  
BOX DE TRAVAIL  
Ref. 080214



# MESTRA®

CAJA DE TRABAJO "ORION"

WORKING BOX "ORION"

BOX DE TRAVAIL "ORION"

Ref. 080214

Nº	Ref.	Descripción	Description
1	080214-01B	Tapa superior	Upper cover
2	080214-02	Cristal	Glass
3	080214-01C	Soporte cristal	Glass support
4	080214-07	Rejilla ventilación	Ventilation grill
5	080253-02	Bisagra	Hinge
6	080214-01E	Marco cristal	Glass frame
7	080214-08	Tapón grande	Large plug
8	080214-01D	Soporte interruptor	Switch support
9	080233-02	Interruptor	Switch
10	080214-09	Tapón pequeño	Small plug
11	080614-20	Fuente de alimentación	Power source
	000230	Fusible 2A	Fuse 2A
12	080490-06	Base de conexión	Connection base
13	080214-01A	Cuerpo principal	Main body
14	080614-16	Tubo aspiración	Suction tube
15	080508-07	Codo aspiración	Suction elbow
16	080214-03	Chapa inferior	Lower sheet
17	080214-01G	Cajón derecho	Right drawer
18	080214-01H	Cajón izquierdo	Left drawer
19	080093-10	Burlete de goma	Rubber strip
20	080214-04	Barra de LED	LED bar
21	080214-05	Difusor	Diffuser
22	080214-11	Tapón astillera	Micromotor support plug
23	080214-01F	Astillera	Micromotor support
24	080214-10	Tubo cables	Cable tube
25	080214-12	Rejilla aspiración	Suction grill
26	000250	Cable 230 V	Cable 230 V
	000252	Cable 115 V	Cable 115 V
27	080508-18	Pata	Leg