


MESTRA[®]

PULIDOR ELECTROLITICO
ELECTROLYTIC POLISHER
POLISSEUR ELECTROLYTIQUE
R-080490



MESTRA®
PULIDOR ELECTROLITICO
ELECTROLYTIC POLISHER
POLISSEUR ELECTROLYTIQUE
R-080490

Nº	Ref.	DESCRIPCIÓN	DESCRIPTION
1	080490-01	Carcasa	Housing
2	080490-02	Cuba	Tank
3	080490-03	Cátodo - (negro)	Cathode - (black)
4	000231	Fusible 4A (230 V)	Fuse 4A (230 V)
	000236	Fusible 7A (115 V)	Fuse 7A (115 V)
5	080490-06	Base con porta fusible	Base with fuse holder
6	000250	Cable conexión 230 V	Connection cable 230 V
	000252	Cable conexión 115 V	Connection cable 115 V
7	080490-10	Lámpara neón 230 V	Neon lamp 230 V
	080090-19B	Lámpara neón 115 V	Neon lamp 115 V
8	080490-11	Potenciómetro 230 V	Potentiometer 230 V
	080490-11B	Potenciómetro 115 V	Potentiometer 115 V
9	080490-12	Disyuntor	Circuit breaker
10	080030-18	Pomo de potenciómetro	Potentiometer knob
11	080490-14	Clavija hembra + (rojo)	Female socket + (red)
12	080490-15	Clavija macho + (rojo) cable con 2 terminales	Male pin + (red) cable with 2 terminals
13	080490-16	Clavija macho - (negro) cable con 2 terminales	Male pin + (black) cable with 2 terminals
14	080490-17	Clavija hembra - (negro)	Female socket + (black)
15	080490-21	Ánodo + (rojo)	Anode + (red)
16	080490-22	Temporizador	Timer
17	080490-23	Amperímetro	Amperimeter
18	080490-25	Disparador	Trigger
19	080490-27	Regleta	Terminal board
20	080490-29	Rectificador	Rectifier
21	080490-30	Transformador 230 V	Transformer 230 V
	080490-30B	Transformador 115 V	Transformer 115 V
22	080490-31	Pata	Leg