


# MESTRA®

CALENTADOR DE CERA ALTA CAPACIDAD (1L)

HIGH CAPACITY WAX HEATER (1L)

RECHAUFFEUR DE CIRE GRANDE CAPACITE (1 L)

R-080076



# MESTRA®

CALENTADOR DE CERA ALTA CAPACIDAD

**HIGH CAPACITY WAX HEATER**

RÉCHAUFFEUR DE PLAQUES BASE

R-080076

Nº	Ref.	Descripción	Description
1	080078-01	Carcasa	Housing
2	080078-02	Tapa inferior	Lower lid
3	080076-02	Control 230V	Control 230V
4	080078-04	Fibra aislante	Insulation fiber
5	080078-05	Sonda	Termal probe
6	080076-01	Cuba	Container
7	040210-07	Junta tórica	O-ring
8	080077-04	Resistencia 230V	Resistor 230V
9	080077-09	Tapa superior	Upper lid
10	080520-16	Base de enchufe	Connection base
11	000250	Cable 230V	Cable 230V
12	050364-04	Pata	Leg
13	080077-05	Junta de silicona	Silicone gasket
14	000231	Fusible 4A	Fuse 4A