


MESTRA®

ARTICULADOR PARA PRÓTESIS BALANCEADAS
ARTICULATOR FOR BALANCED PROSTHESES
 ARTICULATEUR LEGER POUR PROTHESES BALANCEES
 R-010070



Nº	Referencia	Descripción	Description
1	010090-01	Prisionero registro vertical	Vertical register bolt
2	010070-02	Tornillo condilar	Condilar bolt
3	010070-03	Tornillo	Bolt
4	010070-04	Plano oclusal	Occlusal plane
5	010070-05	Pivote muelle	Spring pin
6	010070-06	Muelle	Spring
8	010090-06	Aguja referencia	Reference needle
10	010070-10	Pasador cubetas	Tank pin
11	010070-11	Cubeta de plástico	Plastic tank
12	010070-12	Base de metal	Metal base
13	010070-13	Parte inferior	Lower part
14	010070-14	Parte superior	Upper part
15	010070-15	Columna	Column
16	010070-16	Guía incisal	Incisal guide