


# MESTRA®

BATIDORA PEQUEÑA Y GRANDE

SMALL & LARGE MIXEL CUP

BATTEUR PETIT - GRAND

R-060220 - R-060230



# MESTRA®


BATIDORA PEQUEÑA  
**SMALL MIXEL CUP**  
*BATTEUR PETIT*  
R-060220

Nº	Ref.	CONCEPTO	Description
1	060220 - 01	Manilla	Handle
2	060220 - 02	Racor	Connector
3	060220 - 03	Cuerpo	Body
4	060220 - 04	Junta tórica cuerpo	Body O-ring
5	060220 - 05	Tapa metacrilato Ø102 mm	Methacrylate lid Ø102 mm
6	060220 - 06	Junta tórica tapa Ø86 mm	Lid O-ring Ø86 mm
7	060220 - 07	Junta tórica tornillo	O-ring screw
8	060220 - 08	Tornillo tapa	Lid screw
9	060220 - 09	Junta tórica pala	Blade O-ring
10	060220 - 10	Pala	Blade
11	060220 - 11	Bol Ø102 mm	Bowl Ø102 mm

BATIDORA GRANDE  
**LARGE MIXEL CUP**  
*BATTEUR GRAND*  
R-060230

Nº	Ref.	Descripción	Description
1	060220 - 01	Manilla	Handle
2	060220 - 02	Racor	Connector
3	060220 - 03	Cuerpo	Body
4	060220 - 04	Junta tórica cuerpo	Body O-ring
5	060230 - 05	Tapa metacrilato Ø147 mm	Methacrylate lid Ø147 mm
6	060230 - 06	Junta tórica tapa Ø134 mm	Lid O-ring Ø134 mm
7	060220 - 07	Junta tórica tornillo	O-ring screw
8	060220 - 08	Tornillo tapa	Lid screw
9	060220 - 09	Junta tórica pala	Blade O-ring
10	060230 - 10	Pala	Blade
11	060230 - 11	Bol Ø148 mm	Bowl Ø148 mm

**MESTRA®**  
MICROBOMBA DE VACIO  
**MICRO-VACUUM PUMP**  
MICRO-POMPE A VIDE  
R - 080001



# MESTRA®

MICROBOMBA DE VACIO

**MICRO-VACUUM PUMP**

MICRO-POMPE A VIDE

R-080001

Nº	Ref.	DESCRIPCIÓN	DESCRIPTION
1	080001-01	Carcasa	Housing
2	080001-02	Bomba completa 230V	Complete pump 230V
	080001-02	Bomba completa 115V	Complete pump 115V
3	080001-06	Membrana	Membrane
4	080001-16	Cabezal completo	Head complete
5	080001-17	Goma válvula	Rubber valve
6	080001-07	Junta tórica	O-ring
7	080070-12	Pasa cables	Cable bushing
8	080001-08	Chasis	Lower housing
9	000200	Cable 230V	Cable 230V
	000202	Cable 115V	Cable 115V
10	080490-31	Pata	Legs
11	080490-32	Tornillo pata	Legs screw
12	080010-20	Interruptor	Switch
13	080001-19	Filtro	Filter
14	080001-20	Manguera	Hose