



**MESTRA®**  
 RECORTADORA DE MODELOS  
**MODEL TRIMMER**  
 TAILLE-PLATRE  
 R-080080-1-2-3-4

**MESTRA®**

RECORTADORA DE MODELOS

**MODEL TRIMMER**

TAILLE-PLATRE

R-080080-1-2-3-4

| Nº   | Ref.       | DESCRIPCION                          | DESCRIPTIPON                       |
|------|------------|--------------------------------------|------------------------------------|
| 1*   | 080080-01  | Manguera conexión agua               | Water connection hose              |
| 2    | 000250     | Cable 230V                           | Cable 230V                         |
|      | 000252     | Cable 115V                           | Cable 115V                         |
| 3    | 080080-03  | Disco de carburo de silicio          | Silicon carbide disc               |
|      | 080081-01  | Disco de diamante                    | Diamond disc                       |
| 4    | 080251-09  | Tornillo patas                       | Legs screws                        |
| 5    | 080490-06  | Base de conexión                     | Connection base                    |
| 6    | 080080-40  | Interruptor de seguridad             | Safe switch                        |
| 8    | 000236     | Fusible 7A                           | Fuse 7A                            |
|      | 000237     | Fusible 15A                          | Fuse 15A                           |
|      | 000238     | Fusible 20A                          | Fuse 20A                           |
| 9    | 080080-09  | Condensador 230V                     | Capacitor 230V                     |
|      | 080080-09  | Condensador 115V                     | Capacitor 115V                     |
| 10   | 080080-41  | Motor con condensador 1500 rpm 230V  | Motor with capacitor 1500 rpm 230V |
|      | 080080-41B | Motor con condensador 1500 rpm 115V  | Motor with capacitor 1500 rpm 115V |
|      | 080082-01  | Motor con condensador 3000 rpm 230V  | Motor with capacitor 3000 rpm 230V |
|      | 080082-01B | Motor con condensador 3000 rpm 115V  | Motor with capacitor 3000 rpm 115V |
| 11   | 080080-42  | Tornillos sujeción motor             | Motor camping screws               |
| 12   | 030230-04  | Tuercas del motor                    | Motor nuts                         |
| 13   | 080080-44  | Tapa trasera                         | Rear cover                         |
| 14   | 080080-14  | Junta tórica tapa trasera            | Rear cover O-ring                  |
| 15   | 080080-15  | Porta retén                          | Seal holder                        |
| 16   | 040250-03  | Tornillos de porta retén             | Seal holder screws                 |
| 17   | 080080-45  | Tapa cubre motor                     | Motor hood                         |
| 18   | 080080-18  | Espárragos tapa delantera            | Front cover allen screws           |
| 19   | 080080-19  | Tapa delantera                       | Front cover                        |
| 20   | 080080-20  | Tornillos laterales tapa delantera   | Front cover lateral screws         |
| 21   | 080080-21  | Pieza basculante                     | Tilting device                     |
| 22   | 080080-22  | Tuercas ciegas                       | Tilting device blind nuts          |
| 23*  | 080080-23  | Electroválvula 230V                  | Electrovalve 230V                  |
|      | 080080-23B | Electroválvula 115V                  | Electrovalve 115V                  |
| 24*  | 080080-24  | Manguera grifo                       | Tap hose                           |
| 25   | 080080-46  | Tapa de registro                     | Inspection cover                   |
| 26   | 080080-26  | Patatas de goma                      | Rubber legs                        |
| 30*  | 080080-30  | Racor de conexión grifo              | Tap conector                       |
| 31   | 080080-47  | Tornillos unión de tapa trasera      | Rear cover joining screws          |
| 32*  | 080080-32  | Cuerpo de grifo - Forma L (1500 rpm) | Tap body - L shape (1500 rpm)      |
|      | 080082-02  | Cuerpo de grifo - Forma T (3000 rpm) | Tap body - T shape (3000 rpm)      |
| 33*  | 080080-33  | Arandela de cobre                    | Copper washer                      |
| 34*  | 080080-34  | Mando de grifo                       | Tap control knob                   |
| 35   | 080080-35  | Tornillo baquelita                   | Bakelite screw                     |
| 36   | 080080-56  | Acoplamiento para motor              | Motor coupling                     |
| 37   | 080080-04  | Reten para porta reten               | Seal for seal holder               |
| 38   | 080080-54  | Manguera desagüe                     | Sink hose                          |
| 39   | 080080-49  | Retén acoplamiento motor             | Motor coupling seal                |
| 40   | 080080-53  | Llave fija                           | Spanner                            |
| 41   | 080080-55  | Llave de tubo                        | Socket                             |
| 42** | 080084-01  | Base empotrable                      | Base                               |

\* No existe en 080084 / Does not exist in 080084

\*\* Sólo en 080084 / Only in 080084