








MESTRA[®]

CHORRO DE ARENA "GALAXY"
 "GALAXY" SANDBLASTER
 SABLEUSE "GALAXY"
 R-080235



	
R-080235-50 (10 u.)	R-080235-51



MESTRA®

CHORRO DE ARENA "GALAXY"
"GALAXY" SANDBLASTER
SABLEAUSE "GALAXY"

Ref. 080235

Nº	Ref.	Concepto	Description
1A	080235-38	Bisagra cristal	Glass hinge
1B	080235-39	Tapa trasera carcasa	Housing endloader lid
1C	080235-40	Carcasa	Housing
1D	080235-67	Bandeja perforada con pomo	Slotted tray with knob
1E	080235-42	Base carcasa	Housing base
2	080235-02	Visera cristal	Glass visor
3	080235-03	Cristal	Glass
4	080262-03	Burlete cristal	Glass strip
5	080235-60	Barra de LED	LED bar
6	000250	Cable	Cable
7	080270-31	Porta boquillas	Nozzle holder
8	080235-25	Boquilla carburo de boro	Boron carbide nozzle
9	080235-07	Tuerca pistola	Gun nut
10	030350-03	Espárrago	Stud
11	080235-08	Flexo pistola	Gun arm
12	080235-09	Cuerpo pistola	Gun body
13	080235-10	Inyector	Injector
14	080235-11	Racor prolongación inyector	Prolongation injector connector
15	000312	Racor codo plástico 1/8 Ø6	1/8 elbow Ø6 plastic connector
16	080235-12	Casquillo pistola	Guns sleeve
17	080235-13	Manguera aspiración	Aspiration hose
19	080235-15	Suifideras cesta	Basket rubbers
20	080422-17	Interruptor luz	Switch
21	080251-07	Pasamuros recto doble Ø6	Double connector Ø6
22	080490-06	Base con porta fusibles	Fuse base
23	000231	Fusible 4 A	Fuse 4 A
24	080235-16	Pulsador	Pusher

Nº	Ref.	Concepto	Description
25	000310	Racor codo plástico Ø6 1/4	1/4 elbow Ø6 plastic connector
26	080235-17	Manorregulador	Pressure control
27	080246-07	Acople electroválvula	Couple electrovalve
28	080242-23	Electroválvula	Electrovalve
29	080242-24	Filtro electroválvula	Electrovalve filter
30	080270-33	Tapón vaciado tolva	Drain plug
31	080235-18	Temporizador	Timer
32	080550-11	Pomo temporizador	Time pome
33	080235-19	Tornillo cesta	Basket screw
34	080235-20	Cesta	Basket
35	080235-21	Arandela estriada cesta	Basket washer
36	080235-22	Bocallave	Socket
37	080235-52	Conjunto motor	Motor assembly
38	080010-05	Pasa cables	Cable bushing
39	080235-55	Protector manorregulador	Pressure control protector
40	080075-07	Pomo bandeja	Salver pome
42	080490-31	Pata	Leg
43	080242-28	Manguito	Wool sleeves
44	080270-30	Brida	Flange
45	080233-06	Fuente alimentación	Power source box
46	080242-25	Codo salida aspiración	Exit aspiration elbow
47	080235-62	Protector LED	LED protector
48	080235-63	Tapa protector LED	LED protector lid
49	080235-50	Protector cristal (10 u)	Glass protector (10 u)
50	080235-51	Protector mandos	Command protector
51	080235-66	Rejilla fondo	Bottom grating
52	080235-61	Varilla rejilla	Grating column