



MESTRA®

BOX PARA APLICACIÓN DE ESPRAY
BOX FOR SPRAY APPLICATION
 BOX POUR L'APPLICATION DE SPRAY
 Ref. 080602

MESTRA®

BOX PARA APLICACIÓN DE ESPRAY
BOX FOR SPRAY APPLICATION
BOX POUR L'APPLICATION DE SPRAY
Ref. 080602

Nº	Ref.	CONCEPTO	DESCRIPTION
1	080600-28	Pomo pantalla	Shield knob
2	080600-02	Pantalla lámpara	Lamp shield
3	080600-07	Junta pantalla	Shiels gasket
4	080600-16	Soporte portalámparas	Lamp holder support
5	080600-13	Portalámparas	Lamp holder
6	080600-05	Lámpara 230V	Lamp 230V
	080600-05B	Lámpara 115V	Lamp 115V
7	080602-01	Chapa trasera	Back sheet
8	080084-01	Base Schuko	Schuko plug base
9	080600-08	Carcasa	Housing
10	080600-17	Soporte base porta fusibles	Fuse holder base support
11	080490-06	Base porta fusibles	Fuse holder base
12	080600-01	Visor	Visor
13	080490-31	Pata	Leg
14	080600-22	Filtro estropajo	Cloth filter
15	080600-12	Conjunto filtros metálicos	Metal filter kit
16	080600-40	Rejilla inferior cámara	Chamber lower grid
17	080600-06	Conjunto complete sandwich	Sandwich assembly
18	080400-07	Interruptor	Switch
19	080600-09	Marco frontal	Front frame
20	080086-07	Manguito inferior de manguera	Lower hose sleeve
21	080600-04	Manguera de desagüe	Hose
22	000250	Cable 230V	Cable 230V
	000252	Cable 115V	Cable 115V