



MESTRA®

DEPÓSITO DE AGUA DESTILADA

TANK FOR DISTILLED WATER

RESERVOIR D'EAU DISTILLÉE

Ref. 080620

MESTRA®
 DEPÓSITO DE AGUA DESTILADA
TANK FOR DISTILLED WATER
RESERVOIR D'EAU DISTILLEE
 Ref. 080620

Nº	Ref.	Concepto	Description
1	080620-10	Junta tapa	Lid O-ring
2	000334	Racor recto	Straight fitting
3	080620-04	Válvula entrada	Inlet valve
4	080620-03	Regulador presión	Pressure regulator
5	080620-13	Válvula seguridad	Safe valve
6	080620-12	Acoplador	Connector
7	080250-14	Manómetro	Pressure gauge
8	080620-05	Válvula salida	Outlet valve
9	080620-15	Tapa con válvula de apertura (sin junta)	Lid with open valve (without O-ring)