



**MESTRA®**  
**MÁQUINA DE VAPOR 3L**  
**STEAM MACHINE 3L**  
**GENERATEUR DE VAPEUR 3L**  
R-100193

<b>N.</b>	<b>DESCRIPCIÓN</b>	<b>DESCRIPTION</b>	<b>Ref.</b>
411	Cubierta	Cover	
414	Pulsador	Push button switch	100193-04
415	Cubierta de goma para pulsador	Rubber cap for push button	100193-06
416	Mango de baquelita	Bakelite handgrip	
418	Pistola de vapor completa	Complete steam gun	
419	Pomo regulación electroválvula	Electrovalve regulation knob	100193-18
420	Tanque	Tank	
422	Transformador	Transformer	
423	Porta fusibles	Fuse holder	100193-30
425	Fusible 1A 24V	1A 24V Fuse	100193-25
426	Fusible 10A 220V	10A 220V Fuse	100193-26
429	Junta de teflón	Plug teflon gasket	100193-07
430	Junta de silicona	Plug silikon gasket	100193-08
431	Boquilla pistola de vapor	Steam gun nozzle	100193-05
432	Varilla porta pistola	Gun holder rod	
433	Manguera de vapor	Steam hose	100193-21
434	Resistencias 1300W	1300W Heating elements	100193-14
437	Termostato	Thermostat	100193-09
438	Conector termostato	Thermostat connector	
439	Junta resistencia	Heating element seal	100193-20
440	Conexión tubo nivel	Level pipe fitting	100193-16
441	Clip nivel	Level clamp	100193-17
442	Presostato	Pressure switch	100193-03
444	Electroválvula completa	Complete electrovalve	100193-10
445	Luz piloto	Pilot light	100193-11
446	Manómetro	Manometer	100193-24
448	Asa	Unit handle	
449	Portalámparas indicador nivel	Level indicator socket	100193-32
450	Lámpara indicador nivel 230V	Level indicator bulb 230V	100193-12
	Lámpara indicador nivel 24V	Level indicator bulb 24V	100193-27
	Lámpara indicador nivel 380V	Level indicator bulb 380V	100193-28
	Lámpara indicador nivel 30V	Level indicator bulb 30V	100193-29
452	Interruptor	Boiler switch	100193-13
454	Tapón	Boiler cap	100193-01
456	Tubo nivel	Level pipe	100193-15
458	Conexión manómetro	Manometer fitting	
460-487	Conexión presostato	Pressure switch fitting	
464	Cubierta flexible manguera	Flex White sheathing flex	100193-22
465	Cubierta electroválvula	Solenoid valve cover	
466	Cestillo manguera	Flex Hose flex basket	
469	Retén de llenado	Filler seal	100193-23
484	Soporte barra	Rod holder	
475	Tornillo	Screw	
482	Pata de goma	Rubber foot	
485	Rejilla	Boiler grill	
486	Regleta conexiones	Terminal board	
487	Zumbador	Buzzer	
488	Cable	Electric cable	
-	Pistola completa con tubo	Complete gun with hose	100193-02