


**MESTRA®**  
MECHERO DE GAS  
**GAS BURNER**  
*BEC BUNSEN*

R-050051, R-050052, R-050053, R-050054



# MESTRA®

MECHERO DE GAS

**GAS BURNER**

*BEC BUNSEN*

R-050051, R-050052, R-050053, R-050054

Nº.	Ref.	DESCRIPCIÓN	DESCRIPTION
1	050050-01	Recipiente de cera (opcional)	Wax container (optional)
2	050000-08	Boquilla con regulador	Nozzle with regulador
3	050000-01	Chiclé Ø0.3 (butano/propano)	Injector Ø0.3 (butane/propane)
	050020-01	Chiclé Ø0.5 (gas natural)	Injector Ø0.5 (natural gas)
4	050051-02	Racor delantero	Front fitting
5	050051-03	Junta de goma	Rubber gasket
6	050000-07	Tuerca ciega	Cap nut
7	050051-04	Espiga con junta (opcional)	Dowel with gasket (optional)
8	050051-01	Llave de paso	Gas tap
9	050051-05	Base	Base
10	050000-06	Tornillo	Bolt