



MESTRA



HYBRID Kolinsky-Toray Brush range

It is known that the excellence of a brush is given by four basic attributes: its ability to recover its initial shape after use (recovery), its ability to absorb liquids (capillarity), the ability to release that liquid progressively (dosing), and finally, its ability to slide with very little resistance (slip).

Due to its condition as a natural fibre, **Kolinsky** sable hair has a scaly structure. On the other hand, **Toray** has a smooth structure with holes, so both fibres have advantages and disadvantages. However, if they are combined in the right proportion, a perfect balance between recovery, capillarity, dosing and sliding can be achieved. This is what our new range of **HYBRID** brushes offers.

- ✓ **Advantages from Kolinsky hair:**
excellent liquid retention and subsequent dosing of the same.
- ✓ **Advantages from Toray fiber:**
Excellent ability to recover its shape after use, excellent abrasion resistance, extraordinary glide with hardly any resistance.
- ✓ **Other advantages:**
Ergonomic handle, variety of sizes, very competitive price.



It allows the application of ceramics with precision thanks to its fine tip.



NUMBER	REFERENCE	MATERIAL	SPECIFICATIONS
0	050370	Wood	Hybrid (Kolinsky-Toray)
1	050371	Wood	Hybrid (Kolinsky-Toray)
2	050372	Wood	Hybrid (Kolinsky-Toray)
4	050373	Wood	Hybrid (Kolinsky-Toray)
6	050374	Wood	Hybrid (Kolinsky-Toray)
8	050375	Wood	Hybrid (Kolinsky-Toray)
10	050376	Wood	Hybrid (Kolinsky-Toray)



Nº 0



Nº 1



Nº 2



Nº 4



Nº 6



Nº 8



Nº 10



MESTRA

Talleres Mestraitua SL
 (+34) 944 530 388
 info@mestra.es - www.mestra.es
 48150 SONDIKA (Spain)



Management System
 ISO 9001:2015

www.tuv.com
 ID 1100018003