

Electrovibrators

Ref. 080010 - 080020 - 080025 - 080030



ELECTROVIBRATORS



Mestra electrovibrators are reliable, long-lasting, safe and maintenance-free. They are solidly and sturdily built.

- ✓ Working area made of synthetic rubber, easy to extract and washable.
- ✓ On / Off switch and pilot lamp.
- ✓ Potentiometer for stepless vibration adjustment.
- ✓ Standing on non-slip suction pads.

LARGE ELECTROVIBRATOR

Ref. 080010

For prosthetics laboratories which require high performance in coating and plaster vibration. 270 x 380 mm working area built on 6 dampers, enabling large loads to be vibrated.

Characteristics:

Height 130 mm
Width 390 mm
Depth 300 mm
Weight 13,2 kg
Voltage 230 V - 50-60 Hz
Rater power 275 W
Working area 370 x 275 mm



SQUARE ELECTROVIBRATOR

Ref. 080020

The most versatile and the most widely sold.

- ✓ 220 x 220 mm working area.
- ✓ Allows up to 9 3x type casting rings to be vibrated simultaneously.

Characteristics:

Height 105 mm

Width 235 mm

Depth 235 mm

Weight 6,3 kg

Voltage 230 V - 50-60 Hz

Rater power 250 W

Working area 220 x 220 mm



RECTANGULAR ELECTROVIBRATOR

Ref. 080025

Its characteristics are similar to those of the Square Electro Vibrator, except for its dimensions: 220 x 160 mm.

Characteristics:

Height 105 mm

Width 235 mm

Depth 170 mm

Weight 6,2 kg

Voltage 230 V - 50-60 Hz

Rater power 250 W

Working area 220 x 160 mm



ROUND ELECTROVIBRATOR

Ref. 080030

Its high vibrating power makes this unit particularly suitable for tank emptying. 120 mm diameter working area, allowing 4-1x type casting rings to be vibrated simultaneously, or 1-9x type.

Characteristics:

Height 100 mm

Width 152 mm

Depth 152 mm

Weight 2,5 kg

Voltage 230 V - 50-60 Hz

Rater power 140 W

Working area 120 mm





At **MESTRA** we are specialized in the manufacture and distribution of dental machinery and appliances. We are a Spanish company based in Bizkaia. We currently export to clients in more than fifty countries on five continents and we have more than a thousand references in our catalogue.

We have maintained a line of constant expansion since our foundation in 1945, thanks to our best qualities: the high quality of our products that have earned the ISO 9001:2015 certification by TÜV Rheinland, a close and attentive relationship with our customers, and a constant technical innovation in the design of our products.

Our design philosophy is based on three fundamental pillars: the needs of our clients, the observation of the techniques used in dental laboratories and our more than 77 years of experience providing creative and contemporary solutions to the daily problems of professionals in the sector.

Do you want to take a virtual visit to our training and exhibition center?

Contact us and we'll show you!



Talleres Mestraitua SL
(+34) 944 530 388
info@mestra.es - www.mestra.es
48150 SONDIKA (Spain)

