



# Finishing box ECLIPSE

Ref. 080217



# Finishing box ECLIPSE

Ref. 080217



Finishing box designed for small reworking tasks. Its small dimensions make it ideal for places with little space, where it is important to keep the environment clean without the need for bulky vacuum equipment.

**Eclipse** work box fits anywhere!

It also has a practical carrying handle, which -considering the low weight of the unit- is very useful for transport and storage.

Similar in performance to the "Finishing box Eclipse" (Ref. 080216), but with the special feature of integrating its own suction unit. If required, it is possible to connect the machine to an external extraction source.

## Specifications:

Height 275 mm • Width 360 mm • Depth 400 mm

Weight 2,8 kg • Interior light LED 7W

# Suction and filtering

The machine is equipped with a suction and filter unit with enough power to prevent dust and dirt from escaping during simple jobs.

If necessary, it is possible to connect the machine to an external suction source, thus ensuring efficient particle absorption even for larger jobs.



# Features

---

As with all products in the **Eclipse** family, ergonomics have been taken care of to the utmost: the working position (shoulders/arms/spine) is optimal. The same applies to the view through the large window, which can also be folded down to make it easier to clean the inside of the unit. An LED lamp - in a neutral colour shade- provides shadow-free illumination of all corners of the working chamber and prevents visual fatigue.

The suction motor is silent, avoiding the annoying buzzing sound so characteristic of this kind of equipment. Finally, the filter change and/or cleaning operation is extremely simple: it can be carried out in just one minute without the need for tools of any kind.

In short: the **Eclipse** finishing box is the best ally for professionals who carry out simple reworking work in environments with simple reworking in a confined space.







At **MESTRA** we are specialized in the manufacture and distribution of dental machinery and appliances. We are a Spanish company based in Bizkaia. We currently export to clients in more than fifty countries on five continents and we have more than a thousand references in our catalogue.

We have maintained a line of constant expansion since our foundation in 1945, thanks to our best qualities: the high quality of our products that have earned the ISO 9001:2015 certification by TÜV Rheinland, a close and attentive relationship with our customers, and a constant technical innovation in the design of our products.

Our design philosophy is based on three fundamental pillars: the needs of our clients, the observation of the techniques used in dental laboratories and our more than 77 years of experience providing creative and contemporary solutions to the daily problems of professionals in the sector.

**Do you want to take a virtual visit to our training and exhibition center?**

**Contact us and we'll show you!**



Talleres Mestraitua SL  
(+34) 944 530 388  
info@mestra.es - www.mestra.es  
48150 SONDIKA (Spain)

