



# Sandblaster ECLIPSE AUTOMATIC

Ref. 080218



# Sandblaster ECLIPSE AUTOMATIC

Ref. 080218

Made entirely of technical plastic, a highly resistant material, easy to clean and without oxidation. It allows to place it in any corner of the laboratory.



## Characteristics:

Dimensions 430 mm x 320 mm x 310 mm

Weight 5,5 kg • Max. air consumption 90l/min.

Operating pressure 2 a 5 atm. • Inside light LED



The blowing nozzle is adjustable in the three dimensions of space.

The skeletons are deposited in a basket that rotates by the action of a motor.



Hopper capacity up to 8 kg of sand.

The basket can be easily uncoupled from the drive shaft, which simplifies the preparation of the works.



**MESTRA**

Talleres Mestraitua SL  
(+34) 944 530 388  
info@mestra.es - www.mestra.es  
48150 SONDIKA (Spain)



Management System  
ISO 9001:2015  
www.tuv.com  
ID 1100018003