



ULTRASONIC CLEANER

Ref. 100162





Device specially designed for cleaning retainers, aligners, toothbrushes, mouth guards, whitening trays, dentures, etc. It can also be used as a jewelry cleaner.

Its 28 W power ensures efficient cleaning in just a few minutes. The UV-LED lamp located on the lid is used to disinfect the items placed in the tank while they are being cleaned.

Especificaciones:

Height 122 mm • Width 122 mm • Depth 104 mm

Weight 395 g • Bowl Ø 87 x 38 mm • Capacity 180 ml

Frequency 42 kHz • Timer 5/10 min

Voltage 220 V, 50/60 Hz • Power 28 W



As an added advantage, the small size of the device allows it to be used and stored anywhere. Even in a small suitcase when travelling.

As an additional element, the cleaner has a practical digital countdown timer in seconds.



MESTRA

Talleres Mestraitua SL
(+34) 944 530 388
info@mestra.es - www.mestra.es
48150 SONDIKA (Spain)

