



ALGINATE AUTOMATIC MIXER

Ref. 100191



Alginate automatic mixer

Ref. 100191

Just fill the cup with alginate and water, put it into the machine, press the button, and in a few seconds you will achieve a completely homogeneous mixture, ready for use.

Silent, small, comfortable and above all economical and reliable. Factory supplied with a cup set, a plastic spatula and a practical water dispenser.

It can also be used as plaster mixer with astonishing results.

MAXIMUM MIXTURE: Alginate + water = 80 grs



Compact and small, it can be placed anywhere taking up minimal space. Its design allows an excellent integration within the clinic.

EXTREMELY SIMPLE OPERATION:

Just select the number of measures of alginate you want to mix and press the key. The display shows the mixing time. As simple as this.

Especificaciones

Height 305 mm	Weight 17 kg
Width 205 mm	Voltage 230 V, 50/60 Hz
Depth 265 mm	Consumption 350 W



The machine is factory supplied with four cups, a practical water dispenser, a mixing stick, a spatula and a holder.

The operation consists of four steps:



- 1 Fill the cup with the proportion of water and alginate. Premix and close.



- 2 Put the cup into the machine.



- 3 Close the lid and press the start-up key.



- 4 The mixture will be ready in just 10 seconds.



- * The mixture is easily collected from the cup.
The cleaning is fast and simple.



At **MESTRA** we are specialized in the manufacture and distribution of dental machinery and appliances. We are a Spanish company based in Bizkaia. We currently export to clients in more than fifty countries on five continents and we have more than a thousand references in our catalogue.

We have maintained a line of constant expansion since our foundation in 1945, thanks to our best qualities: the high quality of our products that have earned the ISO 9001:2015 certification by TÜV Rheinland, a close and attentive relationship with our customers, and a constant technical innovation in the design of our products.

Our design philosophy is based on three fundamental pillars: the needs of our clients, the observation of the techniques used in dental laboratories and our more than 77 years of experience providing creative and contemporary solutions to the daily problems of professionals in the sector.

Do you want to take a virtual visit to our training and exhibition center?

Contact us and we'll show you!



Talleres Mestraitua SL
(+34) 944 530 388
info@mestra.es - www.mestra.es
48150 SONDIKA (Spain)

