

ÉLITE induction micromotor

Ref.100285



- ✓ Motor Type: Brushless
- ✓ Vertical/horizontal position
- ✓ Right and left turn
- ✓ Hand/foot speed control
- ✓ Self-diagnosis
- ✓ Error showing



Specifications

Dimensions:	182 x 91 x 205 mm
Weight:	2,5 kg
Max. speed:	50000 rpm
Torque:	7.8 N·cm
Power:	230 W
Voltage:	230 V - 50/60 hz
Handpiece:	Ø 24 x 149 mm
Weight handpiece:	202 g



At **MESTRA**, we are specialized in the manufacture and distribution of dental machinery and appliances. We are a Spanish company based in Bizkaia. We currently export to clients in more than fifty countries on five continents and we have more than a thousand references in our catalogue.

We have maintained a line of constant expansion since our foundation in 1945, thanks to our best qualities: the high quality of our products that have earned the ISO 9001:2015 certification by TÜV Rheinland, a close and attentive relationship with our customers, and a constant technical innovation in the design of our products.

Our design philosophy is based on three fundamental pillars: the needs of our clients, the observation of the techniques used in dental laboratories and our more than 77 years of experience providing creative and contemporary solutions to the daily problems of professionals in the sector.

Do you want to take a virtual visit to our training and exhibition center?

Contact us and we'll show you!



Talleres Mestraitua SL
(+34) 944 530 388
info@mestra.es - www.mestra.es
48150 SONDIKA (Spain)

