




# LED lamp

Ref. 100289V




# LED lamp


Ref. 100289V

LED lamp with articulated arms.


- ✓ 180 LEDs.
- ✓ Temperature color 6500 - 200 K
- ✓ Energy consumed 30 W
- ✓ 3 types of light: cold, neutral and warm.
- ✓ Tactile adjustable light intensity.
- ✓ Flicker-free.
- ✓ With table clamp.



Tactile adjustable light intensity.






Articulated arms.



Adjustable up to 180°.

# Three types of light: cold, neutral and warm

---







At **MESTRA** we are specialized in the manufacture and distribution of dental machinery and appliances. We are a Spanish company based in Bizkaia. We currently export to clients in more than fifty countries on five continents and we have more than a thousand references in our catalogue.

We have maintained a line of constant expansion since our foundation in 1945, thanks to our best qualities: the high quality of our products that have earned the ISO 9001:2015 certification by TÜV Rheinland, a close and attentive relationship with our customers, and a constant technical innovation in the design of our products.

Our design philosophy is based on three fundamental pillars: the needs of our clients, the observation of the techniques used in dental laboratories and our more than 77 years of experience providing creative and contemporary solutions to the daily problems of professionals in the sector.

**Do you want to take a  
virtual visit to our training  
and exhibition center?**

**Contact us and  
we'll show you!**



**MESTRA**

**Talleres Mestraitua SL**  
(+34) 944 530 388  
info@mestra.es - www.mestra.es  
48150 SONDIKA (Spain)



Management  
System  
ISO 9001:2015

www.tuv.com  
ID 1100018003