

SPECIAL OFFER | Clinic

From June 1 to September 30, 2025



icanCLAVE 3F Autoclave

NEW


icanCLAVE 3F is a super-fast autoclave, especially for the sterilisation of unwrapped instruments, which greatly improves the rotation rate of surgical and the convenience of small operations such as dentistry.

- ✓ Printer port.
- ✓ USB port
- ✓ Internal memory.
- ✓ Barcode label printer port.
- ✓ Carrying case.



Specifications

Voltage, Frequency: 220/110V, 50/60Hz • Chamber volume: Aprox. 3L
Power: 2300 W • Number of trays: 1 • Tray material: Aluminium
Tray dimensions: 178 x 285 mm
Chamber dimension: 192(W) x 45(H) x 285(D)



REGULATIONS

- ✓ 97/23/EEC Pressure Equipment Directive.
- ✓ Regulation (EU) 2017/745 on medical devices (MDR)
- ✓ EN13060 Small steam sterilizers.
- ✓ EN61326, EN 61000 Electromagnetic compatibility.
- ✓ EN61010-1, UL 61010-1 Safety requirements.
- ✓ EN61010-2-040, IEC 61010-2-040 Specific requirements for steam sterilizers used for treatment of medical materials.
- ✓ EN 13485:2012/AC:2012
- ✓ ISO 13485:2003

Official technical service ICANCLAVE in Spain.



Official technical service
icanCLAVE
in Spain

REF. 100367


LAUNCH OFFER
1.995 €

18-D Desktop autoclave class B (18 liters)

- ✓ ODE valves, ULUAK vacuum pump, ULKA water pump.
- ✓ Ultra quiet (70 dB).
- ✓ Triple fractionated vacuum system with an independent steam generator, which guarantees fast and safe sterilization cycles.
- ✓ High precision control system for perfect sterilization results (EN13060).
- ✓ Auto-diagnosis with error codes.
- ✓ Fully automatic with pre-programmed sterilization cycles.
- ✓ Chamber and door in stainless steel, with electro-polished finish.
- ✓ Open type water tank, which allows easy filling and cleaning.
- ✓ With printer.

Specifications

Chamber dimensions 247×350 mm	Nominal power 1750 V
Operation temperature 5-40°C	Operation relative humidity Max. 80%
Operation atmospheric pressure 76-106 kPa	Noise < 70 dB
Circuit breaker F16A/400 V	Sterilization temperatures 121/134°C
Distilled water tank capacity 2.5 l	



MESTRA®
by icanCLAVE



Official technical service
icanCLAVE
in Spain

REF. 100365

2.852 €

Sirocco S suction unit

NEW

- ✓ Powerful induction motor.
- ✓ Easy to use and change filter.
- ✓ Exceptionally quiet.
- ✓ Interior LED lamp



Sirocco S
+ Finishing box

REF. 080534AC

699 €

Specifications

Power 350 W	Noise 55 dB
-------------	-------------



Biological espore incubators


- ✓ Simple and easy to use.
- ✓ Essential to carry out the necessary tests that confirm the proper performance of the sterilization devices.
- ✓ Both have four cameras to be able to perform as many verifications at the same time.



Elite model
Adjustable temperature
and digital display

REF. 100370


300 €



Practic model
Functional
and economical

REF. 100375

129 €



SPECIAL OFFER | Clinic

From June 1 to September 30, 2025



Pressure former PlastPress


The most compact positive pressure machine on the market.

Greater accuracy in the forming process and lower noise compared to vacuum machines.

Practical and versatile, it accepts round and square sheets up to 6 mm thickness.

It must be connected directly to a compressor (between 3 and 8 bar).

Sheets \square from 125 to 131 mm
Sheets \circ from 125 to 134 mm



Specifications

Power 450 W

REF. 080544

941 €


Vacuum forming machine

Thermoforming machine for square and round sheets.

Sheets \square from 125 to 134 mm
Sheets \circ from 125 to 134 mm
(adapter available for \circ 120, not included)

Specifications

Power 925W



REF. 080542

651 €


Pressure former

The strong pressurization allows precise molding, regardless of the cut of the model, as well as various thicknesses (0.5 to 4.0 mm). Square and round sheets (127, 125, 120 mm) can be used.

Specifications

Pressure 4-6 bar

Power 450W




REF. 080543

2.422 €

Drying oven for thermoforming sheets

A fan circulates the hot air producing an uniform drying:

- ✓ Draws out moisture and keeps dry thermoforming sheets, to avoid the formation of bubbles.
- ✓ Drying of instruments plaster models, etc.
- ✓ 6 stainless steel trays.
- ✓ Stainless steel housing.
- ✓ 24 hour timer.



Specifications

Temp. 35-75°C

Power 550-650 W



REF. 080118

268 €

Thermoforming sheets

Double packaging (inside-outside).
For all types of use.

CE Medical device



REFERENCE	HARDNESS	APPEARANCE	Dim. 128 x 128 mm thickness	QTY.	PRICE
080563	Hard	Transparent	0.5 mm	25 u.	11€
080564	Hard	Transparent	0.75 mm	25 u.	14,56€
080565	Hard	Transparent	1.0 mm	25 u.	14,56€
080566	Hard	Transparent	1.5 mm	25 u.	18,30€
080567	Hard	Transparent	2.0 mm	25 u.	27€
080568	Soft	Translucent	0.5 mm	25 u.	25€
080569	Soft	Translucent	1.0 mm	12 u.	6,98€
080570	Soft	Translucent	1.5 mm	12 u.	9,27€
080571	Soft	Translucent	2.0 mm	12 u.	13,97€
080572	Soft	Translucent	3.0 mm	12 u.	17,45€

Thermoforming sheets



CE Medical device

Standard square and round thermoforming sheets.
Made of high quality plastic with CE marking certification as medical device.

There is no soot, nor fumes during the heating.
As they do not contain recycled materials, they do not tear while subjected to vacuum.

Available in a wide variety of thicknesses and types.
Supplied in metalized package and vacuum closed.

From
4 €

REF.	TYPE	HARDNESS	APPEARANCE	Dim. 126x126 mm thickness	QTY.	PRICE
080574	Crystal	Hard	Transparent	0.3 mm	10	4€
080575	Crystal	Hard	Transparent	0.75 mm	5	4€
080576	Crystal	Hard	Transparent	1.0 mm	5	4€
080577	Crystal	Hard	Transparent	1.5 mm	5	6€
080578	Crystal	Hard	Transparent	2.0 mm	5	8€
080579	Milky	Hard	Translucent	0.6 mm	10	4€
REF.	TYPE	HARDNESS	APPEARANCE	Dim. 131x131 mm thickness	QTY.	PRICE
080580	Soft	Soft	Translucent	1.0 mm	10	7€
080581	Soft	Soft	Translucent	1.5 mm	10	10€
080582	Soft	Soft	Translucent	2.0 mm	10	12€
080583	Soft	Soft	Translucent	3.0 mm	5	9€
080584	Soft	Soft	Translucent	4.0 mm	5	13€
REF.	TYPE	HARDNESS	APPEARANCE	Dim. Ø125 mm thickness	QTY.	PRICE
080574R	Crystal	Hard	Transparent	0.3 mm	10	4€
080575R	Crystal	Hard	Transparent	0.75 mm	5	4€
080576R	Crystal	Hard	Transparent	1.0 mm	5	5€
080577R	Crystal	Hard	Transparent	1.5 mm	5	7€
080579R	Milky	Hard	Translucent	0.6 mm	10	5€
080580R	Soft	Soft	Translucent	1.0 mm	10	5€
080581R	Soft	Soft	Translucent	1.5 mm	10	9€
080582R	Soft	Soft	Translucent	2.0 mm	10	9€
080583R	Soft	Soft	Translucent	3.0 mm	5	6€
080584R	Soft	Soft	Translucent	4.0 mm	5	11€

SPECIAL OFFER | Clinic

From June 1 to September 30, 2025



Automatic bag sealer with printer



Optional roller
REF. 080436
200 x 240 x 635 mm
17 kg


REF. 080435

2.170 €

Bag sealer

- ✓ Sealing length 305 mm.
- ✓ Sealing width 12 mm.

*Coils not included



REF. 080433

**SUPER OFFER €
265 €**

Self-sealing sterilisation pouch

New bags 5x9 cm,
the smallest from the market!



Specifications

Ref. 080448: 5 x 24 cm	Ref. 080449: 7 x 26 cm
Ref. 080450: 9 x 26 cm	Ref. 080451: 5 x 9 cm
Ref. 080452: 14 x 28 cm	Ref. 080454: 19 x 36 cm
Ref. 080456: 25 x 40 cm	Ref. 080457: 30 x 40 cm
Ref. 080458: 30 x 45 cm	Ref. 080459: 5.7 x 13 cm

200 units
per package

**CHECK
PRICES**

Class I medical device (MDR 2017/745)
Conforming to standards: ISO 11607-1:2017 and EN 868-5:2009

Coils for sealers

Specifications

Ref. 080441: 5 cm x 200 m
Ref. 080442: 7.5 cm x 200 m
Ref. 080443: 10 cm x 200 m
Ref. 080444: 15 cm x 200 m
Ref. 080445: 20 cm x 200 m
Ref. 080446: 25 cm x 200 m
Ref. 080447: 30 cm x 200 m

Class I medical device (MDR 2017/745)
Conforming to standards: ISO 11607-1:2017 and EN 868-5:2009



**CHECK
PRICES**

Barrier film

This adhesive film offers a
practical and effective solution for
dental environments, providing a
layer of protection for instruments
and materials.

Each roll contains 1200 pre-cut
10 x 15 cm adhesives.

Blue colour.



REF. 080439

**LAUNCH OFFER
10 €**




Composite Heater

Device to keep up to 8 composite
syringes warm.

Specifications

Temperature 40-50-70°C
Holes Ø 16 mm



REF. 080077

**SUPER OFFER
225 €**

Tank for distilled water

Made of stainless steel, it features
a pressure regulator and a double
safety valve.

Capacity 18 l



REF. 080620

307 €



Distiller decalcifier

4 liters



REF. 100210

230 €



SPECIAL OFFER | Clinic

From June 1 to September 30, 2025



Ultrasonics bath 2.8, 3.4 and 5.8 liters

- ✓ Temperature.
- ✓ Drain.
- ✓ Metal basket.
- ✓ Made of stainless steel.



5.8 L

REF. 100152

452 €



Inner dimensions:
Height 155 mm
Width 298 mm
Depth 152 mm



3.4 L

REF. 100151

421 €



Inner dimensions:
Height 100 mm
Width 298 mm
Depth 152 mm



2.8 L

REF. 100150

328 €



Inner dimensions:
Height 100 mm
Width 238 mm
Depth 140 mm

Ultrasonic bath Calypso 4 liters




- ✓ With heating
- ✓ With drain valve
- ✓ Includes crystal glass
- ✓ Includes cutter cleaning pliers (ref. 070520)
- ✓ Includes stainless steel basket

Specifications

Cap. tank 4 l
Frequency 40 kHz

REF. 080422

1.124 €




Ultrasonic bath 180 ml with UV-light

NEW



Specifications

Cap. tank 180 ml	Frequency 42 kHz
Cuve Ø 87*38 mm	Timer 5/10 min
Tension 220 50/60 Hz	Power 28 W



REF. 100162

**LAUNCH OFFER
47 €**


Ultrasonic bath for denture cleaning



*The battery is not included.

Specifications

Power supply 2 batteries AA



REF. 100160

CHECK PRICES

Ultrasonic bath 3, 2.5 and 1.4 liters


NEW



2.5 L

REF. 100140

262 €



3 L

REF. 100142

170 €



1.4 L

REF. 100141

170 €



SPECIAL OFFER | Clinic

From June 1 to September 30, 2025



Alginate mixer

- ✓ Transportable.
- ✓ Ergonomic.
- ✓ Two speeds.
- ✓ Spatule included.

CE Medical device
MDR 2017/745



230 mm
175 mm
250 mm
3 kg

REF. 080515

403 €



Gift - One cup
valued in 27,69 €
Ref. 080515-02

Alginate automatic mixer

- ✓ Comfortable automatic mixing.
- ✓ Quick mixing obtention in a few seconds.
- ✓ Easy cleaning of all components.
- ✓ Max mix 80 grs.



Gift - Sealed alginate
container valued
at 93,20 €
Ref. 100189



305 mm
265 mm
205 mm
17 kg

REF. 100191

819 €

Apex locator

- ✓ Automatic calibration.
- ✓ LCD display clear and accurate image.
- ✓ Different colours.
- ✓ Sterilizable accessories.
- ✓ Rechargeable battery.



REF. 100330

SUPER OFFER
291 €

110 mm
80 mm
138 mm
0,4 kg

Motor for endodontics

Specifications

Speed 100-600 rpm
Torque 1.4 N. cm
Power supply batteries



215 mm
85 mm
120 mm
0,3 kg

CE
0120

REF. 100350

SUPER OFFER
612 €

Intraoral camera with autofocus system

- ✓ 21.5-inch LCD touch screen.
- ✓ High definition camera (4K).
- ✓ Autofocus.
- ✓ Image and video recording on USB memory.
- ✓ QR code for downloading images via WIFI.
- ✓ Support bar.
- ✓ Android system.
- ✓ Display of 1 or 4 images on screen.



280 mm
480 mm
400 mm
Screen: 2,52 kg
Camera: 0,1 kg

REF. 100336

749 €

280 mm
480 mm
400 mm
Monitor: 2,52 kg
Cámara: 0,1 kg

Tabletop ultrasonic scaler

CE
0197

80 mm
134 mm
200 mm
0,66 kg



D3

2 Functions
Scaling - Periodontics
Satelec type handpiece
Detachable and autoclavable
with 5 tips.

REF. 100313

250 €

81 mm
159 mm
188 mm
0,65 kg



D5

3 Functions
Scaling - Periodontics
- Endodontics
Satelec type handpiece
Detachable and autoclavable
with 6 tips.

REF. 100316

SUPER OFFER
300 €

230 mm
190 mm
240 mm
1,23 kg



D7

3 Functions
Scaling - Periodontics
- Endodontics
Satelec type handpiece. Detachable and
autoclavable at 135° with 8 tips. With
liquid reservoir or drug treatment.

REF. 100318

SUPER OFFER
622 €

SPECIAL OFFER | Clinic

From June 1 to September 30, 2025



MESTRA

Whitening Lamp

- ✓ Wide spectrum of illumination (from 420 to 472 nm)
- ✓ High power (up to 2000 mW/cm²)
- ✓ Easily adjustable in height and orientation.
- ✓ Its digital display shows the power, times, etc. of the whitening process.
- ✓ Includes remote control 20 disposable mouth splints, two protective glasses and a lip retractor.



REF. 100340
SUPER OFFER
490 €




Whitening Lamp

- ✓ Wide spectrum of illumination.
- ✓ 3 different types of LED light sources: blue, red, (red + blue).
- ✓ 5 blue lights from 430 to 490 nm.
- ✓ 2 high power red lights 610 to 650 nm.
- ✓ Light output power up to 2327 mW/cm².
- ✓ 8 purple lights with sterilization effect.
- ✓ Made of metallic material, in appearance elegant and easy to assemble.
- ✓ Adjustable bleaching time.



REF. 100341
691 €



Gutta Percha cutter

The tip is flexible design to let you cut smoothly.

Press the power button and the tip can heat up to 500°C immediately. Release the button and the tip will cool down quickly.



REF. 100349
LAUNCH OFFER
90 €




Silicone mixer

- ✓ Suitable for 99% of silicones.



REF. 100360
930 €



Specifications

Power 150 W Push force 450 kg ± 10
Environmental temp. 5-40°C Power supply 240V-50Hz

Dry piston compressors with dryer

Family of compressors consisting of two sizes of 80 and 160 litres/minute to which a dryer with automatic drain can be coupled. Extremely quiet operation (60 dB), easy maintenance (dry piston) and excellent reliability.

They are supplied with two pressure gauges (boiler and service pressure), outlet pressure regulator and complete set of valves for draining. The service pressure can reach 8 bars. **Not suitable for milling machines.**




1. Ref. 110340 is also available with HEPA filter (H14) Ref. 110340H
2. The activated carbon filter Ref. 110300-53 is available as an accessory for all compressors.




80 l/min
(with dryer)
REF. 110310

1.334 €



160 l/min
(with dryer)
REF. 110340

1.974 €



CERTIFICATE by TÜV SÜD
Type certificate N° 12202200491258001
Directive 2014/29/UE (pressure vessels)

SPECIAL OFFER | Clinic


From June 1 to September 30, 2025



Eclipse finishing box with built-in suction

This box has similar features to the "Eclipse revision booth" (Ref. 080216), but with the particularity of integrating its own suction unit.

Its small dimensions make it ideal in places with little space, where it is important to keep the environment clean without the need for bulky suction equipment. If necessary, the machine can be connected to an external extraction source.



Front



Back

NEW




REF. 080217

464 €

Eclipse work box

- ✓ Docking coupling with Jeb type suction.
- ✓ Ergonomic work position.
- ✓ LED lamps
- ✓ Easy cleaning.




REF. 080216

378 €



Mini sandblaster for extraoral use

- ✓ Particle size: 50 to 100 microns.
- ✓ Essential for brackets cleaning.
- ✓ Adequacy of micro retentions and crowns.
- ✓ Adjustment of interproximal contacts.
- ✓ Occlusal sandBlasting.
- ✓ Inlay-Onlay.



Specifications


Work pressure 4 to 5.5 bares
Particle size Up to 100 microns

REF. 100122

SUPER OFFER
130 €

Eclipse I Sandblaster - 1 tank

- ✓ One tank and one mobile gun for microblasting.
- ✓ Blower.
- ✓ Operating pressure 2 to 6 bares.
- ✓ Internal LED light.
- ✓ Air consumption 80 l/min



REF. 080221

679 €



NEW

Plutón Sandblaster - 1 tank

- ✓ One high capacity cup.
- ✓ Blowing pressure up to 6 bar.
- ✓ Stainless steel body.
- ✓ High power LED light.
- ✓ Wide view inside.
- ✓ Air consumption 80 l/min
- ✓ Working pressure 2-5 atm



REF. 080254

684 €



NEW

Portable micromotor brushless intramatic

Autonomy 10 hours



CONTROL BOX

Height 45 mm	Width 105 mm
Depth 145 mm	Weight 0,4 kg

REF. 100282

Height 20 mm	Width 20 mm	Depth 155 mm
Weight 0,2 kg	Power 30W	Torque 250 gf-cm
Speed 30000 rpm		

REF. 100282

583 €

SENSE A3 (white)

Brushes

- ✓ Motor with brushes.
- ✓ Manual or pedal speed control.
- ✓ ON/OFF pedal.
- ✓ Displayed speed.
- ✓ Right & Left turning ability.
- ✓ Overload protection function.
- ✓ Fault warning function.



Specifications

Max. speed 35000 rpm	Voltage 0-33V
Torque 3.0 N-cm	Handpiece 150 x Ø 15 mm
Power 70W	Handpiece weight 215 g

REF. 100238

210 €

NEW

SPECIAL OFFER | Clinic

From June 1 to September 30, 2025



MESTRA

Magnifying glasses



LAUNCHING OFFER

NEW

Lightweight, handy and with excellent panoramic and clear vision.

Manufactured with CR-39 resin, higher resistance against abrasion and excellent transparency.

- ✓ Available in 2.5x or 3.5x magnification.
- ✓ Packing dimensions: (WxDxH) 25.5 x 16.5 x 3.5 cm.
- ✓ Weight: 65g.



Augmentation 2,5 x

REF. 100076

LAUNCH OFFER

36 €

Augmentation 3,5 x

REF. 100077

LAUNCH OFFER

36 €

Magnifying glasses

- ✓ Equipped with practical lighting to avoid shadows.

- ✓ Four different magnifications: 1,7X, 2X, 2,5X, 3,5X



REF. 100079

99 €

Safety glasses

Vent valve

REF. 100014

10 €

Anti-fog

REF. 100011

15 €



UV rays

REF. 100013

15 €


Great Vision

REF. 100012

15 €

Face protection

Folding + 10 units



REF. 100022

SUPER OFFER

12 €

Face protection

Frame + 10 units



REF. 100021

SUPER OFFER

13 €

Flask

for duplicating with alginate

- ✓ The inner measurements can be reduced by using rubber pieces



REF. 020206


97 €

Round electro vibrator

Ø120 mm



- ✓ Made of white colour anodized aluminium, and rests on suction pads.
- ✓ Recess of 30 mm for inserting the tray.
- ✓ It is equipped with a potentiometer, a selector switch (+, -) for stepped adjustment of vibrations, and an "ON" light.



REF. 080030

277 €

Polymerizing pneumatic Pot




Specifications

Pressure max. 2 bar

Temp. max. 60°C

Capacity 6 l



Inner dimensions:
Ø 220x170 mm Height



Hose with automatic closing

Double safety valve

REF. 030420

300 €

Water bath for wax plates




For intraoral records

- ✓ Stainless tank and lid.

Specifications

Temperature 30-75°C

Capacity 1l



REF. 080078

402 €

SPECIAL OFFER | Clinic

From June 1 to September 30, 2025



Trolley

- ✓ Made of stainless steel.
- ✓ Pull-out tray (495 x 345 mm)



REF. 100451
SUPER OFFER
145 €

885-1340 mm
550 mm
495 mm
4,5 kg



NEW

Trolley

- ✓ Made of stainless steel.
- ✓ Pull-out tray (495 x 345 mm)




REF. 100452
SUPER OFFER
145 €

785-1220 mm
550 mm
495 mm
4 kg



NEW

Scanner trolley



875 mm
1055 mm (tray)
550 mm (wheels)
600 mm
19 kg



REF. 100450
380 €

NEW

Electric mobile trolley

- ✓ Made of steel sheet and enameled white.
- ✓ Two sockets for equipment connection.
- ✓ Three trays.
- ✓ Four wheels (smooth and balanced movement).

760 mm
290 mm
400 mm
12,5 kg



REF. 100305
297 €

Electric mobile trolley

- ✓ Made of steel sheet and enameled in white.
- ✓ Three trays.
- ✓ Three sockets to power the equipment.
- ✓ Four wheels with smooth and balanced movement.
- ✓ Weight supported per tray: 5 kg

760 mm
300 mm
460 mm
11 kg



REF. 100307
SUPER OFFER
216 €

Electric mobile trolley

- ✓ Made of steel sheet and enameled in white.
- ✓ Three trays.
- ✓ Three sockets, one drawer.
- ✓ Four wheels with smooth and balanced movement.
- ✓ Weight supported per tray: 5 kg.

765 mm
340 mm
400 mm
15 kg



REF. 100308
SUPER OFFER
366 €

Electric mobile trolley

- ✓ Made of steel sheet and enameled white.
- ✓ Two sockets for equipment connection.
- ✓ Four wheels (smooth and balanced movement).
- ✓ Three glass trays.

740 mm
360 mm
460 mm
20 kg



REF. 100306
SUPER OFFER
420 €

Working table

- ✓ Metallic table, wood countertop.
- ✓ Working chair
- ✓ Polishing lathe and two trays.
- ✓ LED lamp.
- ✓ Micromotor.

Gift:

- ✓ Magnifying glass
- ✓ Magnetic bur holder

Work chair

Allows four different regulations.
Wheels with safety brakes
(EN 1335-2:2009)



REF. 100301K2
1.867 €

REF. 100300
266 €



SPECIAL OFFER | Clinic

From June 1 to September 30, 2025



MESTRA

Semiadjustable arcon type articulator Arquimedes Professional A7 Plus

REF. 010099

200 €

Robust structure made of aluminum, ensuring greater stability. Design that provides greater visibility.

- ✓ Semiadjustable / Arcon.
- ✓ Fixed intercondylar distance (110 mm).
- ✓ Curve condylar guide.
- ✓ Regulable condylar guide.
- ✓ Regulable Bennet angle.
- ✓ Central lock.
- ✓ Stabilization system of the condilar guide movement with silicone connection.
- ✓ Support for upper branch in open position.



MESTRA
by bio-art



Standard Professional Set

REF. 010102

SUPER OFFER
265 €

- ✓ Articulator A7 PLUS
- ✓ Standard face bow
- ✓ Case



Elite Professional Set

REF. 010101

379 €

- ✓ Articulator A7 PLUS
- ✓ Elite face bow
- ✓ Case



Semiadjustable arcon type articulator Arquimedes Protrusal A7 Plus E

REF. 010103

SUPER OFFER
285 €

It has exclusive aluminum condylar guides with millimeter adjustment.

- ✓ Stabilization system of the condilar guide movement with spring.
- ✓ Support for upper branch in open position.
- ✓ Fixed intercondylar distance (110 mm).
- ✓ Regulable condylar guide.
- ✓ Regulable Bennet angle.
- ✓ Flat condylar guide.
- ✓ Bilateral locking system.



MESTRA
by bio-art



Standard Protrusal Set

REF. 010104

SUPER OFFER
330 €

- ✓ Articulator A7 PLUS-E
- ✓ Standard face bow
- ✓ Case



Elite Protrusal Set

REF. 010106

SUPER OFFER
450 €

- ✓ Articulator A7 PLUS-E
- ✓ Elite face bow
- ✓ Case




Semiadjustable arcon type articulator + Case 4000 S

REF. 010107

240 €

- ✓ Adjustable intercondylar distance (expanding spindle).
- ✓ Curve condylar guide.
- ✓ Regulable condylar guide.
- ✓ Regulable Bennet angle.
- ✓ Stabilization system of the condilar guide movement with O-ring.



MESTRA
by bio-art



Standard 4000S Set

REF. 010107S

331 €

- ✓ Articulator 4000S
- ✓ Standard face bow
- ✓ Case



Elite 4000S Set

REF. 010107E

418 €

- ✓ Articulator 4000S
- ✓ Elite face bow
- ✓ Case



SPECIAL OFFER | Clinic

From June 1 to September 30, 2025



MESTRA

Arcon type articulator Arquimedes Junior A7 Fix

- ✓ Semiadjustable / Arcon.
- ✓ Fixed intercondylar distance (110 mm).
- ✓ Condylar guide set at average 30°.
- ✓ Bennet angle set at average 15°.
- ✓ Curve condylar guide.
- ✓ Magnetic stabilization system of the condylar guide movement (built-in).
- ✓ Support for upper branch in open position.



REF. 010098

112 €



Standard face bow

An apparatus used to record the spatial position of the patient's upper dental arch in relation to the skull and transfer the record to the articulator. Transfer of the patient record is performed by setting the complete face bow in the upper frame of the articulator.



REF. 010102-50

95 €



CE Producto sanitario

Elite face bow

Articulated and fixation system for the fork that allows a faster registration process. Eliminates the needs for face bow in the assembling process of the casts in the articulator.



REF. 010100-55


198 €



CE Producto sanitario

Orthodontic and denture boxes

Made of plastic.
25 units, assorted colours.




REF. 070600

**CHECK PRICES
PER QUANTITY**

Boxes for orthodontic plates

Made of plastic.
25 units, assorted colours.




REF. 070590

**CHECK PRICES
PER QUANTITY**

Denture boxes

Made of plastic.
100 units.


Blue
Ref. 070000
Pink
Ref. 070001



**CHECK PRICES
PER QUANTITY**

Boxes for orthodontic plates

Blue
Ref. 070002
Green
Ref. 070006
Light green
Ref. 070007
Orange
Ref. 070008



Made of plastic in
different colours

**CHECK PRICES
PER QUANTITY**

Denture boxes



Made of plastic
5 u. assorted

REF. 070622

**CHECK PRICES
PER QUANTITY**

Denture boxes




Made of plastic.
5 u. assorted.

REF. 070618

**CHECK PRICES
PER QUANTITY**

Denture boxes

Made of plastic.
5 u. assorted




REF. 070616

**CHECK PRICES
PER QUANTITY**

Box for prosthesis (with drain)

Made of plastic.
4 u.




REF. 100267

**CHECK PRICES
PER QUANTITY**

Transport and cleaning box

Made of plastic.



REF. 100265

**CHECK PRICES
PER QUANTITY**

Boxes with membrane (10 unit bags) safety lock.

- ✓ Specially designed for the transport and delivery of high quality elements.
- ✓ The inside is cushioned by a high flexibility membrane.
- ✓ Made of highly resistant transparent material with safety lock.
- ✓ Supplied in bags of 10 units.



Ref. 100262



Ref. 100263




Ref. 100264


Made of plastic.

CHECK PRICES
PER QUANTITY

Boxes with membrane (10 unit bags) safety lock.



Ref. 100259



Ref. 100260




Ref. 100261

CHECK PRICES
PER QUANTITY

Transport box




Ref. 100255



Ref. 100256

CHECK PRICES
PER QUANTITY

LED lamp




- ✓ 84 large LEDs.
- ✓ Color temperature 6000 K
- ✓ Tactile dimmable light intensity.
- ✓ Specially rigid arm.
- ✓ With clamp for fastening to the table.

REF. 100288

SUPER OFFER
90 €



LED lamp



- ✓ 90 large LEDs.
- ✓ Color temperature 6000-7000 K
- ✓ Tactile dimmable light intensity.
- ✓ With clamp for fastening to the table.

REF. 100289

SUPER OFFER
90 €

NEW

Accessories for lamp and magnifying glass

REF. 100287 - REF. 100288



Wheeled stand
REF. 100288-10

57 €



REF. 100288-02
Clamp support



REF. 100288-03
Wall bracket




REF. 100288-04
Screwable
table support

25 cm



REF. 100288-07
Desktop stand

LED lamp



- ✓ 117 LEDs.
- ✓ Color temperature 5600-6000 K
- ✓ With clamp for fastening to the table.

REF. 100290

SUPER OFFER
90 €

Magnifying glass



- ✓ 90 white LEDs.
- ✓ Ø 190 mm lens.
- ✓ Clamp to fix it to the table.

3 diopter lens


REF. 100287

5 diopter lens

REF. 100287-5x

SUPER OFFER
76 €

LED lamp



- ✓ 140 LEDs.
- ✓ Colour temperature 3000/6500 K
- ✓ 3 types of light: cold, neutral and warm.
- ✓ Tactile dimmable light intensity.
- ✓ With clamp for fastening to the table.

REF. 100289V

108 €

Swivels up to 180°



Cold light
1 row



Neutral light
2 rows



Warm light
1 row

NEW



MESTRA

Talleres Mestraitua SL
(+34) 944 530 388
info@mestra.es - www.mestra.es
48150 SONDIKA (Spain)

Offer valid unless typographical error or out of stock. The manufacturer reserves the right to make changes without prior notice.



Sistema de
Gestión
ISO 9001:2015
www.tuv.com
ID: 1100018003